X-RAY GENERATION

LIST OF DOCUMENTATION IN THIS BINDER:

- ⊗ SUBSYSTEM MANUAL OPTIMUS R/F
- ⊗ UNIT MANUAL Converter R/F
- O UNIT MANUAL Extension set for an additional tube assembly WG/GWB
- O UNIT MANUAL 26 V DC / 230 V AC Adapter
- UNIT MANUAL Handswitch for OPTIMUS
- O UNIT MANUAL Extension Photomultiplier (SEV)
- O UNIT MANUAL Mains group EWD
- O UNIT MANUAL Power Distribution Unit (PDU)

Note: \otimes indicated document present

LIST OF ALL BINDERS FOR X-RAY GENERATION:

- SUBSYSTEM MANUAL OPTIMUS R/F (this binder)

PHILIPS

INTRODUCTION AND TECHNICAL DATA

Philips Medical Systems DMC GmbH

SERVICE MANUAL 742 **SUBSYSTEM**


INSTALLATION

OPTIMUS R/F

9890 000 02182

FAULT FINDING

REPLACEMENT



PROGRAMMINGS

ADJUSTMENTS

ACCEPTANCE

SERVICE INFORMATION

CAN-controlled X-ray generator of the converter type

PARTS LIST

DMC Hamburg

Printed in Hamburg, Germany

© 2002 Philips Medical Systems ALL RIGHTS RESERVED

SCHEMATIC DRAWINGS

SERVICE MANUAL - SUBSYSTEM

OPTIMUS R/F Author: U. Wente

Type No: 9890 000 02182 Techn. No: Basis 4512 104 70625

Release: 3.5

In case there are any questions concerning this manual, please send this LOPAD via fax to 49/(0)40/5078 2481

File: OPTIMUS_RF_25672_AA

List of pages and drawings (LOPAD)

Manual Order No: 4512 984 25672

released: 09/2002

0.1 1 3.1 3.2

1-0.1 0.2	(d/02.1)		
1-1 16	(d/02.1)		
Z-1.1	(01.0)	A4	4512 982 0092.
Z-1.2	(01.0)	A4	4512 982 0092.
Z-1.3	(01.0)	A4	4512 982 0099.
Z-1.4	(01.0)	A4	4512 982 0099.
Z-6.1	(01.0)	A4	4512 982 0099.
Z-6.2	(01.0)	A4	4512 982 0099.
Z-7.1	(01.0)	АЗ	4512 982 0088.
Z-7.2	(01.0)	АЗ	4512 982 0088.
Z-7.3	(01.0)	АЗ	4512 982 0088.
Z-7.4	(01.1)	АЗ	4512 982 0088.
Z-7.5	(01.0)	АЗ	4512 982 0088.
Z-7.10	(01.0)	АЗ	4512 982 0010.

2-0.1 0.3	(c/02.1)		
2-1 86	(c/02.1)		
2Z-2.0 3x	(01.0)	A4	4512 983 05591
2Z-2.1	(a/01.0)	A4	4512 983 05601
2Z-2.2	(a/01.0)	A4	4512 983 05611
2Z-2.3	(a/01.0)	A4	4512 983 05621
2Z-2.4	(a/01.0)	A4	4512 983 05631
2Z-2.5	(a/01.0)	A4	4512 983 05641
2Z-2.6	(01.0)	A4	4512 983 05651
2Z-2.7	(01.0)	A4	4512 983 05661
2Z-2.8	(01.0)	A4	4512 983 05671
2Z-2.9	(01.0)	A4	4512 983 05681
2Z-2.10	(01.0)	A4	4512 983 05691
2Z-3	(97.0)	АЗ	4512 983 05721
2Z-4	(97.0)	A4	4512 983 05771
2Z-5	(97.0)	A4	4512 983 05731
2Z-10	(02.0)	АЗ	4512 983 05701

3-0.1	(b/02.0)	
3-1 93	(b/02.0)	
3Z-1	(a/01.0)	A4
3Z-21	(97.1)	A4
4-0.1	(a/02.1)	
4-1 25	(a/02.1)	
5Z-1	(a/01.0)	A3
5Z-2	(a/01.0)	A3
6-0.1	(b/02.0)	
6-1 24	(b/02.0)	
7-0.1	(a/01.0)	
7-1 4	(a/01.0)	
8-1	(00.0)	FCO-Checklist
8-2	(00.0)	FCO-Checklist
P-List	OPTIMUS R/F	(840061)

Z1-1.1 (98.0) A3 Z1-1.2 (97.0) A3 4512 983 05561 Z1-2.1 (b/97.0) A3 4512 983 05581 Z1-2.2 (c/97.0) A3 4512 983 05931 Z1-2.3 (99.0) A3 Z1-3.3 (a/01.0) A3 Z1-4.1 (a/01.0) A3	
Z1-4.2 (a/99.0) A3 Z1-5.1 (97.1) A3 4512 983 06551 Z1-6 (02.0) A3 4512 983 05531 Z1-11.1 (02.0) A3 4512 983 05521	
Z1-11.2 (96.0) A4 Z1-13.2 (c/00.0) A3	
Z1-14.1 (b/98.0) A3 4512 983 05541 Z1-14.2 (c/97.1) A3 4512 983 05551 Z1-15.1 (97.0) A3 4512 983 05571 Z1-15.2 (97.1) A3 4512 983 05711	
Z0-2 (02.0)	
Z2-1.0 (96.0) A4 Z2-1.1 (d/99.0) A3 4512 983 06461 Z2-1.2 (e/99.0) A3 4512 983 06471 Z2-1.3 (c/98.0) A3 4512 983 06481 Z2-2.1.1 (97.0) A3 4512 983 06491 Z2-2.1.2 (97.0) A3 4512 983 06501 Z2-2.2.1 (96.0) A3 4512 983 06511 Z2-2.2.2 (a/97.0) A3 4512 983 06521 Z2-5.1 (a/02.0) A3 Z2-5.2 (b/02.0) A3 4512 983 05511 Z2-5.3 (a/02.0) A3 Z2-5.4 (02.0) A3 Z2-13 (d/00.0) A3 Z2-14.1 (b/99.0) A3 Z2-14.2 (a/96.0) A3 Z2-14.3 (a/96.0) A3 Z2-14.4 (97.0) A3 Z2-15.1 (97.0) A3 Z2-16 (97.0) A3 Z2-17 (97.0) A3	

(nur für Fabrik-Auslieferung)

Service software No:

OMA: 4512 116 0220xOMB: 4512 116 0230x

INTRODUCTION AND TECHNICAL DATA TEXT

	Contents	1-0.1
1.	Product information	1-1
1.1.	Applications	1-1
1.2.	Options	1-1
1.2.1.	Hardware options	1-1
1.2.2.	Software options	1-2
2.	Compatibility	1-3
2.1.	Generator components	1-3
2.2.	Tubes	1-3
3.	Mechanical data	1-4
4.	Environmental data	1-5
4.1.	Electrical environment	1-5
4.2.	Climatic conditions	1-5
4.3.	Emission	1-5
5.	Electrical data	1-6
5.1.	Power data and mains conditions	1-6
5.2.	Power supply for applications	1-8
5.3.	Operating data	1-8
5.4.	Power supply	1-9
5.4.1.	Type of power supply	1-9
5.4.2.	Calculating the mains resistances	1-10
5.4.3.	Earth-leakage circuit breaker	1-11
5.4.4.	Emergency-OFF device	1-11
6.	Tools	1-12
7.	Traceable items	1-13
8.	Preparation	1-14
8.1.	Installation material	1-14
8.2.	Cables	1-14
8.3.	Manpower	1-15
9.	Planned maintenance	1-16

DRAWINGS

Generator cabinet Z-1.1
Room layout
Operating panel Z-1.3
Power Distribution Unit (PDU) Z-1.4
Connection of generator Z-6.1
Connection of generator with PDU
Connection diagram
Connection diagram Z-7.2
Connection diagram
Connection diagram Z-7.4
Earthing diagram Z-7.5
Legend for earthing and cabling diagram Z-20.1

1. Product information

The Optimus family of generators for radiography is based on computer-controlled converter technology. The converter operates in the non-audible frequency range.

Applicable options are essentially achieved by releasing software modules using customized PAL ICs. Control between the internal **F**unction **U**nits (FUs) and the external online equipment takes place by a CAN bus. Safety-relevant signals are transferred directly on the so-called "Signal bus".

Units without any CAN interface are operated by the "Adapter for 4 auxiliary units WA" option.

1.1. Applications

- Radiography
- Tomography
- Fluoroscopy

1.2. Options

Component overview according to the commercial catalogue.

Only the versions in the current commercial catalogue can be ordered.

If an existing generator is to be upgraded the commercial department must order:

- MGR0011 (upgrade of an existing configuration)
- + MGRxxxx
- + S/N

1.2.1. Hardware options	+ MGR0011 + S/N
 Dual-speed rotor control Mains transformer PDU: 400 - 480V; 50 / 60Hz, also for 400V mains supply without neutral lead N 	
with taps for 400 / 440 / 460 / 480V	
343 / 380 / 390V max. 50kW!	
- Adapter decade cable for 4 aux. units R/F WB	
- Option rack	
- Extension set for one additional tube	9890 000 0234x
- Tube extension WG	9890 000 0238x
- Operating panel	9890 000 0240x
- Operating module Optimus	9890 000 0278x
- Operating desk data cable 10m, 20m, 30m	9890 000 0241x / 2x / 3x
- Stand for operating panel	9890 000 0244x MGR1482

-	Wall mounting of operating panel	9890 000 0245x
-	26VDC / 230VAC adapter	9890 000 0246x MGR2281
-	Handswitch for Optimus	9890 000 0249x
-	Patient Data Organizer (PDO)	9890 000 0255x MGR2091
-	Decade cable set 14 x 4m top decade \rightarrow AMP decade	9803 704 2010x
-	Extension photo pick-up Optimus	9890 000 0258x

1.2.2. Software options

+ MGR0011 + S/N

Software options are provided by the function key (see also 5Z-1, EZ 139 Central Unit D38). Additional hardware components are not required.

-	Automatic Exposure Control (AEC)	9890 000 0281x ¹⁾	MGR2171
-	Anatomic Programmed Radiography / Fluorography (APR/F) $ \ldots $	9890 000 0282x ¹⁾	MGR2181
-	Automatic Tomo Time Input (TTI)	9890 000 0222x	MGR2121
-	Tomo Density Control (TDC)	9890 000 0223x	MGR2122
-	VARIOFOCUS	9890 000 0227x	MGR2101
_	Area dose calculator	9890 000 0256x	MGR2141

^{1) =} Options only for base 9890 000 0218x Options are always included in base 9890 000 0216x

2. Compatibility

2.1. Generator components

-	Base OPTIMUS	9890 000 0218x
-	H.V. transformer R/F 1 tube, 50kW	9890 000 0270x
-	H.V. transformer R/F 2 tubes, 50kW	9890 000 0271x MGR2061 (Upgrade 1> 2) tubes
-	H.V. transformer R/F 1 tube, 65 / 80kW	9890 000 0272x
-	H.V. transformer R/F 2 tubes, 65 / 80kW	9890 000 0273x MGR2062 (Upgrade 1> 2) tubes
-	50kW extension - R/F	9890 000 0274x
-	65kW extension - R/F	9890 000 0275x
-	80kW extension - R/F	9890 000 0276x
-	Firmware Rel. 3.5	9890 000 0251x

2.2. Tubes

Recommended standard tubes:

- RO 17 50
- SRO 25 50
- SRO 33 100

Further compatible tubes:

-	RO 30	-	SRO 09 51	-	SRO 20 55
-	RO 12 30	-	SRO 13 30	-	SRO 22 50
-	RO 16 48	-	SRO 20 50	-	SRO 32 100

- RO 30 50

Compatible tube housings:

- ROT 350
- ROT 351

The latest information on further tubes which are connectable is available at the service center Hamburg.

Note

When the generator is retrofitted it is important to use the screened cable $3 \times 1.31 \text{mm}^2$ (0722 215 02054) as the stator cable.

If necessary exchange the old stator cable.

3. **Mechanical data**

For installation dimensions and weights see drawings Z-1.1.

Transport data:

		Weigh	ts [kg]	Dimensions [cm]		
Case No.	Contents	net	gross	length	width	height
1	- Generator cabinet - Operating panel - Cables	178	226	210	82	84
	1-tube version:	73	100	77	07	00
2	2-tube version:	88	115	77	67	80
	Contents: H.V. generator					

4. Environmental data

The environmental data comply with to PMS standard UXW 13600.

4.1. Electrical environment

Class S0 - Dedicated mains supply, 3 phases and neutral. Thus single phase voltage is also available.

A low impedance, permanently installed connection, fed in by the step down transformer of the hospital to supply large systems like in MR, CT and X-ray departments is required.

Note

Use always a mains cable with 4 wires and concentric PE-shield, type NYCY.

4.2. Climatic conditions

Ambient temperature	10°C - 40°C
Relative humidity	15% - 90%; no condensation
Relative atmospheric pressure	70kPa - 110kPa

4.3. Emission

Heat dissipation	max. 1200W; average per hour
Noise level, normal operation	≈ 46dBA
Noise level, maximum power operation	≈ 55dBA
EMC	IEC 950

To avoid any possible annoying noise of the implemented fans it is advisable to install the generator cabinet outside the examination room.

5. **Electrical data**

5.1. Power data and mains conditions

	Voltage			
	50kW	65kW	80kW	
Mains voltage	3 x 400	V ±10% (≅ 415V ^{+6%} / 380\	/ _{-5%})	
	3 x 400	3 x 400 / 440 / 460V ±10% *		
	3 x 480	V +6% -10% *		
	* = with external transformer PDU (option)			
	The following connection c	ables are recommended:		
	Input: 3 +1 x 10mr	m ² (L1, L2, L3, PE)		
	Output: 4 x 4mm ² Generator supply 3 x 400V 4 + 1 x 4mm ² Device supply 3 x 220V (option) 3 x 1,5mm ² Switch control and temperature supervision			
	3 x 190 343V ±10% **			
	** = with external transformer; max. 50kW (option)			

	Frequency		
	50kW 65kW 80kW		
Mains frequency		49 61Hz	

	Max. mains current		
Voltage	50kW	65kW	80kW
Exposure: 400V	145A	190A	230A
440V	135A	180A	215A
460V	125A	170A	210A
480V	120A	160A	205A
190V	300A	-	-
Short-time power consumption [I x U x √3]	100kVA	132kVA	160kVA
Fluoro: 400V	-	8A	
480V	-	7A	
Fuse protection	35A	50A	
(slow blow)	100A at ≤ 240V	-	-
Connected load [I _{Fuse} x U x √3]	25kVA 35kVA		·VA
Emergency power supply:			
static (Inverter)	Short-time power consumption [I x U x $\sqrt{3}$]		
dynamic (Diesel generator with flywheel mass)	Connected load [I _{Fuse} x U x √3]		

	Mains resistance		
Voltage	50kW	65kW	80kW
400V	$\leq 300 m\Omega$	≤ 200mΩ	
440V	≤ 3 50mΩ	≤ 24 0mΩ	
460V	≤ 3 50mΩ	≤ 240mΩ	
480V	≤ 400m Ω	≤ 300 mΩ	
480V valid for DOD only	\leq 300m Ω	\leq 240m Ω \leq 180m Ω	
	Note 500m Ω is the absolute max. mains resistance.		

Power supply for applications 5.2.

	Generator power				
Supply	50kW 65kW 80kW				
Output 1	230V / 400V; gene	erator switched and protected	d> I ₁		
Output 2	230V / 400V; gene	230V / 400V; generator switched and protected> I ₂			
Output 3	230V / 400V; generator not switched and protected> I ₃				
	Sum value I _{max (1 +2 +3)} ≤ 10A				
Only with external transformer PDU					
Output 4	127V / 220V generator switched and protected> I ₄				
Output 5	double socket 127V / 15A to phase L1; generator not switched and protected> I ₅				
	Sum value I _{max (4 +5)} ≤ 16A				

5.3. Operating data

	Generator power			
Data	50kW	65kW	80kW	
Exposure: Tube current	1 650mA	1 900mA	1 1100mA	
Tube voltage	4	40 150kV in kV- or %-steps		
mAs product		0,5 850mAs		
Exposure time	1ms 6s / 16s			
Exposure frequency	≤ 12exp./s			
Interfacing option for	door conta	door contact, external radiation warning indicator		
Fluoro: Tube current		0,25 6mA		
Tube voltage	40 110kV			
kV/mA curves		3		
Setting time	≤ 1s			


5.4. Power supply

5.4.1. Type of power supply

3 phase WYE L1 X-ray installation L3 (N) \[\frac{1}{2} \]


- 400V
- 440V / 460V / 480V with external mains transformer PDU 9890 000 0260x.
- Neutral not required if the external mains transformer PDU 9890 000 0260x is ordered.
- 190V ... 343V with external mains transformer 9803 720 8100x (max. 50kW).

3 phase DELTA, balanced earth or floating



- External mains transformer PDU 9890 000 0260x is required.
- 400V / 440V / 460V / 480V
- 190V ... 343V with external mains transformer 9803 720 8100x (max. 50kW).





- External mains transformer PDU 9890 000 0260x is required (requires modification at the EMC-filter of the kV power unit).
- 400V / 440V / 460V / 480V

Caution!

Ensure the sequence of phases in the wall junction box corresponds to designations L1, L2, L3.


5.4.2. Calculating the mains resistances

Note

The cross section of lead l_3 must not exceed 25mm². (See figure below).

If possible the sum of R₀, R₁, R₂ and R₃ should be smaller than the R_{XG} requires.

With higher internal mains resistances the generator output is reduced correspondingly.



- designates the mains resistance of the distributor transformer R_0
- depends on the length of lead I₁ between distributor transformer and mains distributor R_1 and on the selected cross section as well:

==>
$$R_1 = I_1 \times R_{Cu}$$
 R_{Cu} from table below

- consists of upstream elements such as: R_2
 - Emergency-OFF switch 4.0m Ω
 - Earth-leakage circuit breaker 5.5mΩ
- depends on the length of lead I₃ between mains distributor and wall junction box R_3 and on the selected cross section as well:

==>
$$R_3 = I_3 \times R_{Cu}$$
 R_{Cu} from table below

The resistances consider the go and return lines so that the calculation can be based on simple cable lengths.

Copper cross section [mm ²]	Resistance R _{Cu} [mΩ/m]
16	2.19
25	1.4
35	1.0
50	0.7
70	0.5
95	0.38
120	0.30
150	0.24

Note

 $500m\Omega$ is the **absolute max.** mains resistance.

ALL RIGHTS RESERVED

5.4.3. Earth-leakage circuit breaker

To be provided between mains fuse and X-ray installation depending on local regulations.

Siemens earth-leakage circuit breaker:

- Order No.: 5SZ3 466 OKG00
- Rated fault current 30mA
- Rated current 63A
- Connection terminals for wire cross sections of up to 25mm²

5.4.4. Emergency-OFF device

To be provided depending on local regulations.

There are 2 possibilities:

- 1. All the Emergency-OFF buttons are connected in series and looped into the switch-ON circuit (12VDC) of the generator.
- The Emergency-OFF circuit acts on an external mains contactor which switches OFF the power before it is fed into the generator.

6. **Tools**

- Service engineer standard tool kit
- Service-PC:

IBM-compatible, 640kB RAM, 3.5" floppy disk drive, ≥ 1 serial port

- Installation and service software OMA: 4512 116 022xx / OMB: 4512 116 023xx. Supplied on two floppy disks within this generator service manual.
- PC-hardkey (DIAGGEN):

Necessary to carry out the installation and to run the service software (special programmings, fault finding).

- 0-modem cable:

Minimum length is distance between generator cabinet and operating desk.

Male 25-pole D-Sub connector at the generator side.

A 5m data cable of bucky controller can be used: 4512 130 5693x

- Mains resistance measuring instrument
- Dose measuring instrument
- mAs-meter
- Multimeter
- Digital oscilloscope with 2-beam memory
- Recommended PLCC extraction tool (AMP 822154-1): 2422 487 89772

7. Traceable items

The following items have serial numbers of the following format when delivered ex factory:

- 1. Generator cabinet 6 digit serial number
- 3. Operating desk 8 digit serial number

Preparation 8.

Connection of the generator:	see drawing Z-6.1
Connection of the generator with P	DU: see drawing Z-6.2
Operating panel:	see drawing Z-1.3
Connection diagram:	see drawing Z-7.1/.2/.3/.
Earthing diagram:	see drawing Z-7.5
Legend for earthing and cabeling:	see drawing Z-7.10

Installation material 8.1.

To be ordered from the service department of PMS Hamburg:

- inclusive connection block (25mm²) for mains supply and connection block (10mm²) for unit supply.
- Relay for radiation warning indicator 4512 100 4523x 1 interface relay with a floating contact (230V/1A) is included in the scope of delivery of the generator.

8.2. **Cables**

H.V. cables

with O3 / O3 plugs:	. 9806 402 6xx02
length:	. 6m - 30m in steps of 2m
capacity:	. 155pF/m
diameter:	. 16.5mm

The cable length is indicated by the 9th and 10th digit of the numeric code.

Thermal contact cable

-	2-wire screened for 1 excess temperature switch	4512 100 66151
_	10-wire screened for additional supervision like	
	temperature alarm switch, buzzer, selection indicator	0722 215 19005

Stator cable

$3 \times 1.31 \text{mm}^2$	screened	 0722 215 02054
J ∧ 1.51111111 ,	SCIECTICU	 0122 213 02034

Note

The above described cables are part of the pre-installed systems.

AMPLIMAT cable

with D-Sub and 3-Plus plug:

12m	 9890 000 01721
16m	 9890 000 01731
20m	 9890 000 01741
24m	 9890 000 01751

Caution!

or

AMPLIMAT cables 9803 507 0xx02 (for hybrid measuring chambers 9803 509 xxxxx) with 3-Plus plugs at both ends must be connected in the generator by the following adapter for each cable:

Adapter for AMPLIMAT cable: 4512 108 09042. The generator includes 1 adapter.

The hybrid measuring chambers 9803 509 xxxxx require connection (chassis) between contacts:

```
D-Sub end GND (13) <---> RF 0V (8) (generator input)
```

3-Plus end GND (N) <---> RF 0V (J) (generator input)

This connection is established by the adapter for the AMPLIMAT cable. See drawing Z1-6 "Basic interface".

In case a hybrid measuring chamber 9803 509 xxxxx is not operated with the required

```
AMPLIMAT cable . . . 3-Plus / 3-Plus . . . . 9803 507 0xx02 but with
```

AMPLIMAT cable ... D-Sub / 3-Plus 9890 000 017xx

make sure to establish this connection (13 <---> 8) in the D-Sub connector!

For ALC measuring chambers 9890 000 016xx connection GND <---> RF 0V is not permitted. Therefore, ALC measuring chambers AMPLIMAT cables 9890 000 017xx should always be used.

Operating desk

Note

Use the shortest cables. Noise immunity increases.

cable set	10m	9890 000 02411
	20m	9890 000 02421
	30m	9890 000 02431

8.3. Manpower

At least two persons are necessary to insert the H.V. tank in the generator cabinet. The weight of the 2-tube version is about 88kg.


9. Planned maintenance


The technical documentation for carrying out maintenance work in compliance with the applicable regulations are available at the responsible authority of Philips Medical Systems.


The importance of having maintenance implemented is pointed out to the operator in the operating instructions.

It must be guaranteed that the person carrying out maintenance work knows about the respective national regulations and that this person observes these regulations throughout all steps of maintenance work.

1-16 OPTIMUS R/F (d/02.1)OPTIMUS_RF_1_d021







- 1) Wall junction box
- 2)Lateral clearance unless there is an adjacent cabinet
- 3) Filler panel
- 4) Wall-cabinet spacing angle


weight: 210 kg

Generator Cabinet Mechanical dimensions

4512 982 0092.


(01.0)


Z - 1.1



1) With no other cabinets beside them

E= Control cabinet (E not for OPTIMUS CD) C= Operating desk





Room layout

(01.0)

Z - 1.2

(01.0)


4512 982 0099.


© Philips Medical Systems

→202 **→**

825

Z - 1.3






 1) Screw connection (screws 7 x 60, dowels 8) only when requested. It is actually not needed





Power distribution unit (PDU) 9890 000 0260x Dimensions and weight




1) Space with no other cabinets beside them.

Connection of generator

4512 982 0099. © Philips Medical Systems


(01.0)

Z - 6.1




*) Space with no other cabinets beside them.


Connection of generator with Power Distribution Unit (PDU)



OPTIMUS R/F Connection diagram

(01.0)





4512 982 0088. (01.0)


OPTIMUS R/F Connection diagram

Z - 7.3


© Philips Medical Systems

1	WT			Stand (s)
device interface 1				<u> </u>
	△ 1WAX1	dx	release 1	· · · · · · · · · · · · · · · · · · ·
	<u> </u>	dx	release 2	; ; > free attacha
	△ 1WAX3	dx	release 3	
	△ 1WAX4	dx	release 4	
	△ 1WAX11	dx	format contacts bucky 1 (bucky/tomo okay)	· ·
	△ 1WAX12	dx	format contacts bucky 2 (bucky/tomo okay)	·
	△ 1WAX14	dx	EXON	-
	△ 1WAX21	dx	tomo times ←	
	△ 1WAX22	dx	tomo trajectories →	·
	<u> </u>	dx	APR extension (2xRGDV, 6xAPR)	
	△ 1WAX24	dx	PSC external patient size compensation	
				<u> </u>
device interface 2				:
	△ 2WAX1	dx	release 1	
	△ 2WAX2	dx	release 2	
	△ 2WAX3	dx	release 3	
	△ 2WAX4	dx	release 4	
	△ 2WAX11	dx	format contacts bucky 1 (bucky/tomo okay)	<u></u>
	△ 2WAX12	dx	format contacts bucky 2 (bucky/tomo okay)	
	△ 2WAX14	dx	EXON	•
	△ 2WAX21	dx	tomo times ←	: :
	△ 2WAX22	dx	tomo trajectories —>	
	△ 2WAX23	dx	APR extension (8xAPR)	<u> </u>
	△ 2WAX24	dx	PSC external patient size compensation	: :
Bucky/tomo or examination units without CAN—interface				
	△ 1WBX1	dx	release 1	
	△ 1WBX2	dx	release 2	· · · · · · · · · · · · · · · · · · ·
	△ 1WBX3	dx	release 3	:
	△ 1WBX4	dx	release 4	······································
	△ 1WBX11	dx	format contacts bucky	
	△ 1WBX12	dx	fluoro: Dose rate limit/lock in	· ·
	△ 1WBX13	dx	Format contacts serial changer	······································
	△ 1WBX14	dx	EXON	
	△ 1WBX21	dx dx	Change over cassette/camera	
	<u>Δ</u> (WBλ2) Δ 1WBX22	dx	Tubeshield decade	······································
				 :
	△ 1WBX23	dx	External APRT	 :
	△ 1WBX24		PSC external patient size compensation	·

OPTIMUS R/F Connection diagram



OPTIMUS R/F Earthing diagram





(01.0)

© Philips Medical Systems

 Heading symbols - X-ray beam symbols image intensifier image distributor I I camera collimator film plane TV camera X—ray tube electronics cabinet angio electronics cabinet DSI electronics cobinet anglo
electronics cobinet DSI
operators console
electronics cobinet N-ray generator
high tension tank frontal
H.V. changeover switch
hard capy unit
electronics cobinet TOMO
laser hard copy unit
web operators console
web electronics cabinet TOMO
monitor
web electronics cabinet TOMO
monitor
web electronics cabinet TOMO
monitor
web electronics cabinet DMO
monitor
web electronics cabinet web electronics
connection box wall bracket
overhead connection box ceiling crane
electronics cabinet
bucky table
ceiling crane longitudinal carriage
wall bracket
ceiling crane
vertical DIAGNOST 1/2/4
bucky DIAGNOST 1/2/4
bucky DIAGNOST YE/VT
web bucky DIAGNOST
wall connection box

- Cable symbols





Legend for earthing and cabling diagram

INSTALLATION

TEXT

	Contents	2-0.1
1.	Installing the wall junction box	2-1
2.	Preparatory work	2-2
2.1.	Mounting of the H.V. generator in the cabinet	2-2
2.2.	Electrical connection of the H.V. generator	2-3
3.	Installing the operating panel	2-5
3.1.	Desk version	2-5
3.2.	Stand version	2-7
3.3.	Wall mounted version	2-8
3.4.	Supporting angle version	2-9
3.5.	Additional release switch	2-10
4.	Electrical connection	2-11
4.1.	Earthing	2-11
4.2.	Mains connection	2-11
4.2.1.	Mains connection of the generator	2-11
4.2.2.	Mains connection of the PDU	2-12
4.3.	Stator connection	2-13
4.3.1.	Shielding	2-13
4.3.2.	Connection	2-14
4.4.	Signal cables	2-17
4.4.1.	Room decade cable	2-17
4.4.2.	Tube supervision	2-19
4.4.3.	CAN interface	2-19
4.4.4.	Adapter for 4 auxiliary units	2-20
4.4.5. 4.4.6.	Dose inputs	2-23 2-23
4.4.0. 4.5.	H.V. cables generator side	
4.5. 4.6.	Emergency-OFF circuit	2-24
5.	Hardware programming	2-25
6.	Switching-ON of the generator	2-26
7.	Installation software XRGSCOPE	2-27
7.1.	PC and generator settings to avoid problems during up/downloading of CU complete files	2-27
7.1.1.	Preparation of the service PC to guarantee a safe loading process	2-27
7.1.2.	Preparation of the generator	2-27
7.2.	Installation procedure	2-29
8.	Setting-to-work overview	2-32
8.1.	Configuration	2-33
8.1.1.	Date and time	2-33
8.1.2.	Mains data	2-33
8.1.3.	Tubes	2-34
	Tube data set	2-34
	Tube speed selection	2-35
	Tube limits	2-35
	Capacitance of tube connection	2-36
8.1.3.5.	Tube operating modes	2-37

2-0.2

8.1.3.6.	Disable tube	2-37
8.2.	Registration devices	2-38
8.2.1.	Data set A B	2-38
8.2.2.	Interface assignment	2-42
8.2.3.		2-43
8.3.	Tube adjustment	2-44
8.3.1.	Tube conditioning	2-44
8.3.1.1.		2-44
		2-45
	Tube adaptation	2-49
	·	2-49
		2-50
		2-52
8.3.3.	Final tube adjustment work	2-53
8.3.4.	-	2-54
8.4.	Dose rate control	2-55
8.4.1.		2-55
8.4.2.	•	2-55
		2-56
	, , , , , , , , , , , , , , , , , , , ,	2-57
	· · · · · · · · · · · · · · · · · · ·	2-58
8.4.3.		2-59
8.4.4.		2-61
8.4.5.		2-61
	1,3	
8.5.	Application limits	2-62
8.5.1.	X-mode limits	2-62
8.5.2.		2-63
8.5.3.	•	2-63
8.6.		2-64
8.6.1.	1 0 0 0	2-65
8.6.2.	Language	2-66
8.6.3.	1 0 0 1 7	2-66
8.6.4.	1 0 0 (1)	2-69
	•	2-70
	Creating sub menus	2-70
	3 3	2-71
8.6.5.	0 0	2-72
	, i	2-76
	\	2-76
8.6.6.	0. 17 0	2-77
8.6.7.	č	2-78
8.6.8.	1 3	2-79
8.6.9.	9	2-80
8.7.	Option: Tomo Density Control TDC	2-82
8.8.	Option: VARIOFOCUS	2-83
8.9.	Option: Area Dose Calculator	2-83
8.10.	Acceptance test	2-83
8.11.		2-84
8.12.	-	2-84
9.	Labels	2-85
10.	Final installation work	2-86

DRAWINGS

RGDV programming (3x)	2Z-2.0
RGDV programming example: 1	2Z-2.1
RGDV programming example: 2	2Z-2.2
RGDV programming example: 3	2Z-2.3
RGDV programming example: 4	2Z-2.4
RGDV programming example: 5	2Z-2.5
RGDV programming example: 6	2Z-2.6
RGDV programming example: 7	2Z-2.7
RGDV programming example: 8	2Z-2.8
RGDV programming example: 9	2Z-2.9
RGDV programming example: 10	2Z-2.10
List of characters	2Z-3
Data sets of chambers	2Z-4
Programming of device interfaces	2Z-5
Labelling	2Z-10

1. Installing the wall junction box

Mount the wall junction box at the place where the generator is intended to be installed.
 (See drawing "Connection of generator" in section 1 and manual UNIT 4512 103 75380 for wall junction boxes).

- If necessary, install the optional surge arrester WN inside the wall junction box. (See surge arrester documentation.)
- If applicable, mount the filler panels of the generator to the wall junction box.
- Have the mains cable present at the clinic connected to mains terminal MEX by a person who is authorized for this job.
- Check the phase sequence of L1, L2 and L3.



Warning!

Switch OFF the mains supply present at the clinic and make sure that it cannot be switched ON again accidentally.

2. Preparatory work

2.1. Mounting of the H.V. generator in the cabinet

Caution!

Do not tilt the H.V. generator while transporting it.

In case of a tilting angle larger than 45°, the setting-to-work of the generator can be started not before a waiting time of about 8 hours has passed. Otherwise the H.V. generator may be destroyed by electrical sparkover.


Unpack generator cabinet E.

 In case the packing material is strongly soiled with oil check the oil level. Repair it if necessary.

Tolerance: $\pm 2mm$

Oil: Shell Diala G in 2.5l container

4512 148 43172



· Remove the deaerating screw completely from the cover of the H.V. generator. Only this way the precision of the high voltage measuring divider corresponds to the specification. In case of return shipment of the H.V. generator this screw must be fixed again. Therefore, keep the screw laying on top of the cover.

Caution!

Make sure that no foreign matter falls into the oil. Otherwise the transformer must be exchanged.

- Take the two transport bars from the rear side of the cabinet.
- Lift the H.V. generator into the generator cabinet with the transport bars. The 4 connecting bolts GX1001 to 1004 must point at the front of the generator cabinet.

ALL RIGHTS RESERVED


2.2. Electrical connection of the H.V. generator

• Connect the H.V. generator electrically:

Generator Connection		on		
version	from	<>	to	Remark
	E1 (GND)	<>	GX1100 (GND)	Ground
50/65/80kW	ZX12	<>	G100X15	Route the cables along the front and left-
	ZX35	<>	G100X14	hand edge of the H.V. generator. Fix them.
	QC13:1	<>	GX1003	Twist the cables! Note The sequence of the connecting bolts is not in numerical order. See drawing page 2-4.
50kW	QC03:1	<>	GX1002	Push the screening cap forward over the connecting bolts and tighten it. Attach the converter cables including the screening to the screening cap with cable ties.
	GX1001	<>	GX1003	The 50kW version might have direct links on each side or a link on one side and a choke of 1 6 loops on the other side for the reason of kV symmetry.
	GX1004	<>	GX1002	Note Do not change these links or chokes.
	QC13:1	<>	GX1001	Twist the cables!
ec/oolan	QC03:1	<>	GX1002	The sequence of the connecting bolts is not in numerical order. See drawing page 2-4.
65/80kW	2QC13:1	<>	GX1003	Push the screening cap forward over the
	2QC03:1	<>	GX1004	connecting bolts and tighten it. Attach the converter cables including the screening to the screening cap with cable ties.
	WGX61	<>	GK1:1	
50/65/80kW	WGX67	<>	GK1:2	
2 nd tube	WGX62	<>	GK2:1	
	WGX68	<>	GK2:2	

• Turn the two earthing angles of the H.V. generator outward and screw them on to the members of the cabinet.





3. Installing the operating panel

3.1. Desk version

See "Operating panel" in section 1.

Accessories:


- 2 feet for the unit
- 2 elastic buffers, black
- 5 insert strips for the RGDV buttons
- sheet with RGDV symbols
- release switch
- · Unpack the desk carefully.
- Mount the release switch on the left-hand or right-hand side of the desk:


Mount the holding bracket to the edge of the desk (1) with the two M4x10 countersunk screws.

For visual reasons the release button should be in line with the +/- buttons on the control desk. Use the appropriate holes in the bracket.

Slide the release switch onto the bracket. Fasten it parallel to the desk edge with the two M4x10 cheese-head screws, securing rings and washers (2).

- Screw in the 2 feet for the unit at the bottom of the desk.
- Glue the 2 black elastic buffers to the front edges of the bottom of the desk such that they are acting as the front feet.







• Define the assignment of the RGDV buttons 1 ... 8. Glue the respective symbols to the insert strips which are provided with subsidiary lines (1).

- Raise the keyboard about 3mm above the desk. Use an allen key (2).
- Push the insert strips under the keyboard foil. Press the angulated, protruding end of each insert strip into the housing of the desk (3).
- Lower the keyboard to its initial position (4).
- · Remove the cable cover at the rear side of the desk.
- · Connect the cables:

Release switch

Supply cable EZX20 C200X1 EZX6 earth Data cable EZX46 C300X1





- Check the function programming plug for X44 as shown in Z2-5.2 and put it into socket EZX44.
- Provide drag relief for the supply and data cables with the clamp present on the desk.

C300X3

· Screw on the cable cover. Make sure that the cable drag relief device of the release switch (1 cable tie) remains under the cover.

ALL RIGHTS RESERVED

3.2. Stand version

See "Operating panel" in section 1.

Additional accessories:

- 4 dowels S10
- 4 hexagon cap screws 8 x 60mm
- 4 washers
- · Position the desk stand according to the respective room layout.
- · Mark the fixing holes on the floor.
- Set the 4 dowels supplied into the floor (drill bit: 10mm).
- · Screw on the desk stand with 4 screws and washers.
- Route the supply and data cables from the bottom to the top within the desk stand.
 Provide the cables with drag relief.
 Cable ends including plugs should protrude beyond the edge of the desk by about 500mm.
- Mount the release switch as described in chapter 3.1.
- Assign the RGDV buttons 1 ... 8 to the desired symbols as described in chapter 3.1.
- Connect the cables to the desk as described in previous chapter 3.1.
- Screw on the cable cover.
 Ensure the cable drag relief device of the release switch (1 cable tie) remains under the cover.
- · Mount the operating panel on the stand.

3.3. Wall mounted version


See "Operating panel" in section 1.

Additional accessories:

- 4 ball-head bolts
- 4 dowels S8
- 4 hexagon cap screws 5 x 30mm
- 4 washers
- 2 screws 4 x 10mm
- 2 angle plates
- 4 nuts
- · Screw on the angle plates into the wall frame. The short ends of the angles must be pointing upwards.
- Screw the 4 ball-head bolts into the wall support.
- · Mark the 4 fixing holes of the wall frame at the respective place on the wall.
- · Set the dowels supplied into the wall (drill bit: 8mm).
- · Screw on the wall frame with the 4 hexagon cap screws and washers.
- · Provide drag relief for the supply and data cables in the wall frame.
 - Cable ends including plugs should protrude beyond the edge of the desk by about 500mm.
- Mount the release switch as described in chapter 3.1.
- Assign the symbols to the desired RGDV buttons 1 ... 8 as described in chapter 3.1.
- · Connect the cables to the desk as described in previous chapter 3.1.
- · Mount the operating panel on the wall frame and fix it with the left two screws.
- · Screw on the cable cover. Ensure the cable drag relief device of the release switch (1 cable tie) remains under the cover.

The wall frame is designed symmetrically. In case connection cables come from above the frame can be mounted upside down.

Only the ball-head bolts and the angle plates still keep their position.




Supporting angle version 3.4.

• Screw the ball-head bolts (2) into the supporting angles (1):

> Left angle ---> on the left at the bottom Right angle ---> on the right at the bottom

- Press the ball-head bolts (2) into the respective snap bushing of the desk (4).
- Fix the supporting angles (2) to the desk (4) with the 2 screws M4 (3).



3.5. Additional release switch

An optional second release switch is supplied with a longer spiral cable: 9890 000 0249x

The scope of delivery includes various wall hooks and an adapter cable. Electrical connection is made in parallel with the existing release switch which is mounted on the desk itself.


• To do this, plug the pins of the adapter cable into the D-Sub connector of the existing release switch.

Sequence:

Adapter: connector pin	<>	D-Sub: connector pin
1	<>	6
2	<>	9
3	<>	7



Reference: Figure below and drawing Z1-11.1 "Operating panel C" in section Z1 "Schematic drawings".



ALL RIGHTS RESERVED

4. Electrical connection

4.1. Earthing

See "Earthing diagram" in section 1.

4.2. Mains connection

4.2.1. Mains connection of the generator



Warning!

Switch OFF the mains supply present at the clinic and make sure that it cannot be switched ON again accidentally.

See "Connection diagram" in section 1.

• Measure the internal mains resistance at the terminal MEX with a suitable measuring instrument.

Required max. mains resistance at generator input:

Meine veltege	Mains resistance		
Mains voltage	30kW	50kW	65/80kW
190V *	-	40m $Ω$	-
220V *	130mΩ	60mΩ	-
240V *	160mΩ	80mΩ	-
380V	500mΩ	300m $Ω$	200mΩ
400V	500mΩ	300m $Ω$	200mΩ
440V	500mΩ	350mΩ	240mΩ
460V	500mΩ	350mΩ	240mΩ
480V	500mΩ	400mΩ	300mΩ

^{*} with external mains transformer (max 50kW)

Maximum permissible internal mains resistance: $500m\Omega$

Internal resistance of Power Distribution Unit PDU: $20m\Omega$ at 50Hz

 $23m\Omega$ at 60Hz

Caution!

Connect phase wires in correct phase sequence.

• Connect the mains cable of the generator to terminal MEX: L1 / L2 / L3 within the wall connection box. If the optional Power Distribution Unit PDU WT is fitted, connect the cables at that point to terminal WTX2.

• Connect the examination unit supply (max. 10A) to (230V / 400V):

Output 1 Geometry power ENX2101/...5

Output 2 I.I. / TV; GECO; UI; SYSCO ENX2201/...9

Output 3 Imaging ENX2301/...9

Output 3 not switched by generator

4.2.2. Mains connection of the PDU

See delivered Unit Manual "Power Distribution Unit".

OPTIMUS_RF_2_c021

4.3. Stator connection

4.3.1. Shielding

Caution!


To suppress interferences of the high-speed rotor control, the stator connections must be provided with a 360° screen at the tube and generator end.

General remarks:

- Always use screened cables: 0722 215 02054.
- Shorten the stator cable to the required length.
 Do not accommodate excess lengths at the generator.
- Keep stator cable separate from all the other signal cables to avoid interference.
- · Earth the screen at both cable ends.

Screening procedure:

- Remove any enamel or dirt from the clamp providing drag relief in the tube housing to make sure the clamp is conductive.
- Remove the plastic covering around the clamp, about 1cm (0.5").
- Wrap copper foil around the visible screen of the cable until the original diameter of the cable is obtained.
- Remove the present red wire going from the screen end to the earthing point of the tube housing.
- Fix the screen of the stator cable with the clamp.
 Ensure that the clamp is secured and the ground contact works!




4.3.2. Connection

Caution!

Do not mix up the phases, otherwise components of the rotor control may be destroyed.

At the tube end

• Place the jumpers across terminals 100 and 200 according to the figure.



• Connect the stator cable:

```
wire 1
        ---> phase U
wire 2
        ---> phase V
wire 3
        ---> phase W
```

• Earth the screening of the stator cable at the tube housing with the metallic clamp.

At the generator end: One-tube version

See "Connection diagram" in section 1.

- Connect the stator cable to the terminal EX1100 (U-V-W).
- Check the stator connection by measuring the resistances:

```
U - V = wire 1 - 2 \approx 11\Omega

U - W = wire 1 - 3 \approx 20\Omega


V - W = wire 2 - 3 \approx 9\Omega
```

 If an inductance meter is available, measure the following inductance values:

$$U - V = wire 1 - 2 = 57mH \pm 10\%$$

 $V - W = wire 2 - 3 = 34mH \pm 10\%$

- Fix the screen below the screening clamp.
- Relieve the tension of the stator cable by a cable tie.



At the generator end: Two-tube version

See "Connection diagram" in section 1.


- Connect the stator cables to the terminals EWG:K11/K12.
- Check the stator connections of both tubes by measuring the resistances:

 If an inductance meter is available, measure the following inductance values:

$$U - V = wire 1 - 2 = 57mH \pm 10\%$$

 $V - W = wire 2 - 3 = 34mH \pm 10\%$

- · Fix the screen below the screening clamp.
- Relieve the tension of the stator cables by cable ties.



2-15

At the generator end: Three-tube version

See "Connection diagram" in section 1.


· Connect the stator cables to the terminals E1WG:K12 and E2WG:K11/K12.

· Check the stator connections of all 3 tubes by measuring the resistances:

• If an inductance meter is available, measure the following inductance values:

```
U - V
        = wire 1 - 2 = 57mH \pm10%
V – W
        = wire 2 - 3 = 34mH \pm 10\%
```

- Fix the screens below the screening clamps.
- · Relieve the tension of the stator cables by cable ties.



ALL RIGHTS RESERVED

4.4. Signal cables

See: - "Connection diagram" in section 1.

- Z1-6 "Basic interface" in section "Schematic drawings"

4.4.1. Room decade cable

· Connect the door switches at the generator:

or


In case no switch is present link: pin 8 <---> pin 10

EZ150 K1:

max. switching and loading current = 1A

Caution!

Make sure the polarity of the relay is correct.



• In case of tube 2 or 3 is located in room 2 or 3 or while room decades

> WG X3 or 1WG X3 or 2WG X3

are intended to be used for room supervision the soldering link WG W1 must be placed.

Only when this link is placed it is guaranteed that relay WG K3 pulls up and room decade WG X3 is activated when switch-over of the tube takes place.

See Z1-14.xx "Tube extension".


• If needed connect an external relay for each examination room to control external radiation warning devices.

One relay inclusive cable is part of delivery, Additional ones can be ordered: 4512 100 4523.

A mounting place is reserved on the mains connection terminal MEX of the wall junction box.

Caution!

Make sure the polarity of the relay is correct.




4.4.2. Tube supervision

- · Connect the thermal switch or the thermal sensor of the tube housing assembly.
- · Check polarity of contacts 9 and 10 to CAN adapter to remote systems like D96.

For U.S.A. and U.K. only:

· Connect the HHS-lamp to indicate the selected tube housing assembly.



4.4.3. CAN interface

Only for examination units which are provided with a CAN system interface.

· Connect the following plugs:

System	Connector			
	EZX23 signal bus	EZX42 or EZX42-1 system CAN	EZX43 or EZX43-1 system CAN	
BukyDiagnost TH / TH2	Χ		X	
DigitalDiagnost	Х	Х	Х	
Thoravision	Х	Х	Х	
EasyDiagnost with Bucky unit	Х	Х	Х	

4.4.4. Adapter for 4 auxiliary units

Adapter for 4 auxiliary units RAD (Bucky, Tomo) WA / 1WA / 2WA

Used for Bucky examination units which provide their control signals individually via decade cables.

Each of the release circuits and the Bucky decades can be assigned to one or several of the RGDV buttons 1...8 via software programming.

Survey: "Connection diagram" in section 1

Z1-1.2 "Block diagram Expansions" in "Schematic drawings" section

Detail: Z1-15.1 "Adapter 4 aux. units WA / 1WA / 2WA" in "Schematic drawings" section

It provides:

- 1. 4 release decades to be used for
 - grid / syncmaster auxiliaries
 - HHS cassette present interlock
- 2. 2 Bucky decades to be used for
 - format size contacts (collimator, side fields ON/OFF)
 - Bucky-tomo remote switch-over
 - tomo ready condition
 - Bucky ready condition
 - HHS cassette present interlock to be activated by SW programming (see chapter 8.2)

After activation via SW: Signals not provided by the system must be simulated by jumpers.

- 3. 1 decade EXON signal for the system via opto coupler
- 4. 1 decade (output) for 8 tomo trajectories
- 5. 1 decade for 8 tomo time inputs
- 6. 1 decade for 8 external APRT
 - WA backpanel programmed as 1WA = 2 AUX + 6 APRT
 - WA backpanel programmed as 2WA = 8 APRT
- 7. 1 decade for PSC (Patient Size Correction)

ALL RIGHTS RESERVED

Adapter for 4 auxiliary units R/F (Fluoro, Bucky) WB

Used for fluoroscopy and Bucky examination units which provide their control signals individually via decade cables.

Each of the release circuits and the Bucky decades can be assigned to one or several of the RGDV buttons 1...8 via software programming.

Survey: "Connection diagram" in section 1

Z1-1.2 "Block diagram Expansions" in "Schematic drawings" section

Detail: Z1-15.2 "Adapter 4 aux. units WB" in "Schematic drawings" section

It provides:

- 1. 4 release decades to be used for
 - grid / syncmaster auxiliaries
 - HSS cassette present interlock
- 2. 1 Bucky decade to be used for
 - format size contacts (collimator, side fields ON/OFF)
 - Bucky system ready condition
- 3. 1 fluoroscopy decade supporting
 - dose rate limiter (short SID)
 - lock-in
- 4. 1 decade EXON signal for the system via opto coupler
- 5. 1 decade cassette / camera remote switch-over
- 6. 1 decade for image intensifier and Scopomat measuring chamber formats

- 7. 1 system decade providing external signals
 - generator ready display
 - prep or fluoro display
 - 5 minute fluoro warning display / buzzer
 - under exposure display
 - system reset for under-exposure / 5min buzzer
 - EXON signal for the system (potential free relay contact)
 - door contact
- 8. 1 decade for 8 external APRT
- 9. 1 decade for PSC (Patient Size Correction)

4.4.5. Dose inputs

Connect the measuring chambers to the D-Sub connectors EZX21 / 22 / 31 / 32 / 41.

There are no assignement restrictions because the measuring chambers are allocated to the auxiliaries in SW programming.

• Withdraw pins 101-102-103 or A-D-H for measuring field selection at the junior / extremity measuring chamber.

These measuring chambers have only one measuring field. The terminal for the left-hand field is used in other configurations for switching over intensification and must not be connected here.

Connect the photosensor dose output of DSI from Central Part E (BHX 12) to EZX41.
 Only EZX41 sends the necessary signal REL (or REL_CH = RELease CHamber) for automatic and non-automatic techniques.

Note

or

AMPLIMAT cables 9803 507 0xx02 (for hybrid measuring chambers 9803 509 xxxxx) with 3-Plus plugs at both ends **must** be connected in the generator by the following adapter for each cable:

Adapter for AMPLIMAT cable: 4512 108 09042. The generator includes 1 adapter.

The hybrid measuring chambers 9803 509 xxxxx require connection (chassis) between contacts

```
D-Sub end GND (13) <---> RF 0V (8) (generator input)
```

3-Plus end GND (N) <---> RF 0V (J) (generator input)

This connection is established by the adapter for the AMPLIMAT cable. See Z1-6 "Basic interface" in section "Schematic drawings" of the generator manual.

In case of a hybrid measuring chamber 9803 509 xxxxx is not operated with the required

```
AMPLIMAT cables . . . . 3-Plus / 3-Plus . . . . 9803 507 0xx02
```

but with

```
AMPLIMAT cables . . . . D-Sub / 3-Plus . . . . 9890 000 017xx
```

make sure to establish this connection (13 <---> 8) in the D-Sub connector!

For ALC measuring chambers 9890 000 016xx connection GND <---> RF 0V is not permitted. Therefore, ALC measuring chambers AMPLIMAT cables 9890 000 017xx should always be used.

4.4.6. Patient Data Organizer PDO (option)

See operator's manual Patient Data Organizer.

4.5. H.V. cables generator side

See "Connection diagram" in section 1.

- · Mark the H.V. cables at the generator and the tube end with the correct polarity.
- Fix the H.V. cables on the left-hand side of the wall junction box on the middle rail to provide drag relief for the cables. The short ends of the H.V. cables which are going to the H.V. generator must be routed in downward direction in this area.

The free cable lengths including plugs should be about 1.5m.

 Twist the H.V. cables counter-clockwise by one turn and connect them to the H.V. generator. The twisting of the cables allows that the H.V. cables can be put into a loop when the cabinet is placed against the wall.

The H.V. sockets should always be filled with some oil. At least the lower half of the plugs must be wet with oil.

Caution!

Do not use a silicone washer.

Do not grease the plugs with silicone.

The union nuts of the high-voltage connectors must be tightened up to ensure good electrical contact for screening.

Only high-voltage connectors which have threaded flange halves may be used.

Older high-voltage cables still have connectors where the flange halves are kept together with a spring washer.

In such cases the modification kit 4512 103 8085x is required.

Emergency-OFF circuit 4.6.

• Connect the emergency-OFF buttons to EZX4:1/2. If not necessary, link pins 1 - 2.

See Z1-2.1 "Power supply" in "Schematic drawings" section and Z2-5.2 "Backpanel Basic rack-2Z" in the "Wiring diagrams" section.

2-24 OPTIMUS R/F (c/02.1)OPTIMUS_RF_2_c021

5. Hardware programming

Programmings on PCB EZ150 basic interface:

Note

Never change jumper W1.

 Voltage supply for the amplifiers of connected measuring chambers:

Voltage\Soldering link	EZ 150 W2	EZ 150 W3
15V default	OFF	ON
40V	ON	OFF

Working voltage range for ALC measuring chambers: 15 ... 45V

Working voltage range for hybrid measuring chambers: 40 ... 45V


ALC measuring chambers can be recognized by the code No. 4512 104 xxxxx.

Hybrid measuring chambers are based on code No. 4512 102/103 xxxxx.

- Set the gain factor for AEC techniques with jumper EZ150:W4:
 - Factor 1 ==> W4 in position 3 = default
 For screen/film combination with at least one system speed ≤ 200.
 - Factor 4 ==> W4 in position 1
 For screen/film combinations with all system speeds > 200.
- The software programming has to be set accordingly.

The rest of the generator hardware has been properly programmed at the factory.

If required, refer to section 5: PROGRAMMINGS.



Switch-ON of the generator 6.

- Switch ON the fuses present at the clinic.
- Switch ON automatic circuit-breakers ENF1, ENF2 and ENF3.

The yellow LED on EN100 power ON circuit must be illuminated.

7. Installation software XRGSCOPE

7.1. PC and generator settings to avoid problems during up/downloading of CU complete files

Optimus R/F release 2.x and 3.x CMOS data are up/downloaded in one string without handshake.

Any kind of interruption can cause the loading process to fail.

Problems occur mainly during the download to the PC.

A download file which is not complete cannot be used as a safety backup file.

7.1.1. Preparation of the service PC to guarantee a safe loading process

Start XRGSCOPE always from DOS if possible.

When using any WINDOWS version:

- · Switch OFF all screensavers.
- · Do not run other programs.
- · Do not insert any CD in the drive.

Any kind of power management of the PC hardware (BIOS) as well as the windows power management should be switched OFF.

If connected to mains power some of these might be automatically OFF.

7.1.2. Preparation of the generator

Preparation of generators without a CAN interface:

• Switch ON the generator.

The loading process can be started once relay ENK1 has been energized.

Preparation of generators which are connected via a CAN interface:

- BuckyDiagnost TH and TH2
- DigitalDiagnost
- Thoravision
- EasyDiagnost with Bucky unit
- Switch OFF the generator.
- · Disconnect the following plugs:

System	Connector			
	EZX23 signal bus	EZX42 or EZX42-1 system CAN	EZX43 or EZX43-1 system CAN	
BukyDiagnost TH / TH2	Χ		X	
DigitalDiagnost	Х	Х	Х	
Thoravision	Х	Х	Х	
EasyDiagnost with Bucky unit	X	X	X	

· Switch ON the generator.

Note

The download procedure must not be started before relay ENK1 has been energized at least 2 minutes after the generator has been switched ON.


Note

In case of mailing CU complete download files or any other xxx.TDL files use a zipped file format. These files are ASCII files which might be destroyed while being mailed.

7.2. Installation procedure

Provide the service PC with the hardware key and switch it ON.
 The hardware key provides access to special program settings and to menu "Faultfind".
 Standard programming is possible without a hardware key.

 Connect the PC to X5 on EZ139 CENTRAL UNIT CU via a serial data cable: (A 5m long data cable can be ordered via 12NC: 4512 130 56931)



• Insert the floppy disks containing the self-unpacking exe-files of the firmware in the disk drive of the PC:

OMA: 4512 116 022xxOMB: 4512 116 023xx

For unpacking on the harddisk of the PC about 5MB are needed.

- Generate a directory e.g. [C:\OPT_R_35] on the PC by typing <md C:\OPT_R_35> or use WIN commands.
- Copy the firmware from both floppy disks to the PC into the same directory C:\OPT_R_35 by typing
 <copy A:*.* C:\OPT_R_35> or use WIN commands.
- Start unpacking the programs by typing <OMAxxxxx.exe> and <OMBxxxxx.exe> or doubleclick on <*.exe> -files.

The programs unpack all files needed for the update of the firmware and the newest service tools.

- After unpacking [OMAxxxxx.exe] and [OMBxxxxx.exe] can be deleted on the harddisk by typing
 del OMAxxxxx.exe> and <del OMBxxxxx.exe>.
- For the current contents of [OMA/OMBxxxxx.exe] read [OMA/OMBxxxxx.txt] on the floppy disks.
- For a new installation of the generator firmware see section 4 REPLACEMENT, chapter 3 "Exchange of firmware ...".
- Call the installation program by typing <**xrgscope**> or <**xrgscope** lcd> for PCs with LCD screen.

 Start with a XRGSCOPE-screen CUSTOMER.tdl with data which are actually stored in CUSTOMER.tdl. Type "XRGSCOPE customer"

Whenever data screens like 'error log index' are saved to a xxx.TDL file (function save <F3> appears in the bottom line), 'customer data' being saved in the CUSTOMER.tdl file are attached to the saved data screen. It helps to separate saved screen files of different sites, customers or rooms in the same hospital.

Site data must be stored in CUSTOMER.tdl file, only data of this file are attached to the saved screens. One can save site specific customer data in self-made files using <F3>. To recall site data use the load <F4> function.

Procedure:

Either the CUSTOMER.tdl screen is open or open the 'Customer Data' screen. Push <F4> and select a site data file. The old data screen comes up. Now save this screen with <F3> typing in CUSTOMER as file name.

- Customer Name:
- City / State:
- Country:
- Generator Location:
- Generator Serial Number:
- Generator 12 NC:
- Memo1:
- Memo2:
- Memo3:

After <ESC> the following menu line appears:

Select "OPTIMUS".

The following menu line appears:

Pro	gram	Adjust	Accep	ot Faultfind	Quit

ALL RIGHTS RESERVED

General information:

- Button <F1> <help> Call help / cancel help.

- Button **<F2> <transmit>** Store screen contents / data set in the generator ==> transmit to generator.

- Button <F3> <save> Store data screen on disk.

For an open data screen the path desired can be selected.

Button <F4> <load> Load data set from disk. The desired path can be selected.

Button <ESC>
 Commands one step back. Can be used repeatedly.


Fields with ↓ Select the possible range of values by pushing <RETURN>.

The data are specified by the generator as fixed values.

- Fields with [...] Input of data via the keyboard.

Error numbers which appear at the beginning of the programming procedure must be erased from the screen by pushing the **<RETURN>** key.

Current data files, for instance, for online help, tube types, APR programming are available in the PHILIPS-Intranet. Use path: *http://technet.best.ms.philips.com*/ and pull down menu as shown below.



If the installation program is called with **<xrgscope** ?> the possible starting parameters for the service program are listed.

Setting-to-work overview 8.

Note

The programming of a generator must take place in the sequence specified below.

As long as no tube or RGDV has been assigned there is no display at all on the desk except "PHILIPS OPTIMUS".

- 8.1. Configuration
- Registration devices 8.2.
 - RESET the generator
- Tube adjustment
 - · RESET the generator
- 8.4. Dose rate control
 - · RESET the generator
- 8.5. Application limits
 - · RESET the generator
- Human interface 8.6.
 - · RESET the generator
- Option: Tomo Density Control TDC 8.7.
- 8.8. Option: VARIOFOCUS
- Option: Area Dose Calculator
- 8.10. Acceptance test
- 8.11. Interlock facility for APR modification
- 8.12. Backup of all configuration data

8.1. Configuration

· Switch the generator ON.

8.1.1. Date and time

· Select menu:

PROGRAM/ DATE AND TIME

• Enter the respective local data.

8.1.2. Mains data

· Select menu:

PROGRAM/ MAINS DATA

• Select the nominal value of the mains voltage U.

Range: 380V, 400V, 440V, 480V

Default: 400V

If 460V is present program 480V. If 415V is present program 400V.

• Enter the maximum internal mains resistance Ri.

Range: $0 \dots 500 \text{m}\Omega$

Depending on the internal mains resistance and the mains voltage the generator calculates the maximum possible output.

8.1.3. Tubes

Note

Generators which are connected by a CAN interface have to be prepared as described below.

Preparation for:

- BuckyDiagnost TH and TH2
- DigitalDiagnost
- Thoravision
- EasyDiagnost with Bucky unit
- · Disconnect the following plugs:

System	Connector			
	EZX23 signal bus	EZX42 or EZX42-1 system CAN	EZX43 or EZX43-1 system CAN	
BukyDiagnost TH / TH2	X		X	
DigitalDiagnost	Х	Х	Х	
Thoravision	Х	Х	Х	
EasyDiagnost with Bucky unit	Х	Х	Х	

· Switch ON the generator.

Note

The download procedure must not be started before relay ENK1 has been energized at least 2 minutes after the generator has been switched ON.

Do not care if the generator is not programmed yet.

8.1.3.1. Tube data set

· Select menu:

PROGRAM/ TUBES/TUBE 1 ... 3/ TUBE 1 ... 3 DATA SET

- Start the displayed file TUBExxx.tdl with <RETURN>.
 - All the permitted combinations of tube type and housing type are listed in a window.
- Select the respective combination of tube type and housing type from the list and push <RETURN>.
- · RESET the generator.

Then the data which have been configured up to now are read by the processor when the system is started.

ALL RIGHTS RESERVED

8.1.3.2. Tube speed selection

Depending on the type of tube loaded the anode speed is automatically programmed.

Caution!

Wrong programming can cause tube problems.

· Select menu:

PROGRAM/ TUBES/TUBE 1 ... 3/ TUBE 1 ... 3 SPEED SELECTION

RPM \ tube type	RO	SRO
Exposure rotation [RPM]	3000	9000
Fast exposure rotation [RPM]	n/a	n/a
Fluoroscopy rotation [RPM]	3000	3000

8.1.3.3. Tube limits

· Select menu:

PROGRAM/ TUBES/ TUBE LIMITS

· Program the maximum working voltage which is indicated on the data label:

Max. tube voltage limit Range: 40 ... 150kV Default: 150kV

Adaptation of the tube takes place only up to this limit.

Note

If older tubes are to be operated on this generator, it is urgently recommended that the max. kV used in practical operation so far be specified instead of the theoretically possible value.

The max. kV value should be determined during the conditioning procedure as described in chapter 8.3.1.

After adaptation of a tube the upper kV limit is displayed for each focus of each tube under:

Adapted to [kV]: e.g. 125kV

All the other limit programmings are performed by the generator automatically and do not usually have to be observed.

8.1.3.4. Capacitance of tube connection

· Select menu:

PROGRAM/ TUBES/ CAPACITANCE TUBE CONNECTOR

Range: 2.000 ... 10.000nF

The total capacitance for each tube connected is indicated:

$$C = \frac{1}{2} \left(C_{H.V. generator} + C_{H.V. cable} \right)$$

= 4.550nFDefault for H.V. generator + 20m H.V. cable (155pF/m)

$$C_c \times L \qquad \qquad C_c = \text{ specific cable capacitance in [pF/m]} \\ C [nF] = 3 + ----- \\ L = \text{ single cable length in [m]} \\$$

Cinale length [rel	Capacitance tube connection [nF]			
Single length [m]	For 155pF/m cable	For 200pF/m cable		
14	4.085	4.400		
16	4.240	4.600		
18	4.395	4.800		
20	4.550	5.000		
22	4.705	5.200		
24	4.860	-		
26	5.015	-		
28	5.170	-		
30	5.325	-		

The high-voltage cables type 9806 402 6xx02 currently being supplied have a capacitance of 155pF/m.

8.1.3.5. Tube operating modes

· Select menu:

PROGRAM/ TUBES/ TUBE OPERATING MODE

- Intermediate boost:

Select ... Disable = During preparation the rated filament current is applied (default).

Enable = During preparation a reduced filament current is applied.

After the release of exposure boosting takes place for a short time before the exposure is released. Effective with tube currents > 80% of max. value.

- Rotation prolongation after PREP:

Select ... Disable = The tube is braked as soon as preparation has been cancelled.

Enable = After cancellation of preparation the tube is only braked after 30s. Within this

time preparation can be repeated as often as necessary. Recommended for paediatrics and casualty rooms.

The enable mode works with high speed rotor control unit only.

8.1.3.6. Disable tube

For correction of the configuration.

· Select menu:

PROGRAM/ TUBES/ DISABLE TUBE

When the tube is disabled the above stored data set of the tube is erased. To enable the tube the data set has to be loaded again.

Registration devices 8.2.

8.2.1. Data set A ... B

Select menu:

PROGRAM/ REGISTRATION DEVICES/ RGDV 1 ... 8/ DATA SET A ... B

Program the data set A and B of RGDV 1 ... 8 for all exam. / aux. units desired.

Data set A Room: Room number of the exam. / aux. unit for room decade (radiation warning display and door contact). Tube: Tube assignment for the exam. / aux. unit. Release circuit number: Number of the release decade of the release circuit adaptation unit WA, WB programmed (e.g.: 1 for X1 etc., see Z1-1.2). Enable handswitch at generator **Enable release switch Enable external Programming** desk: at generator desk release switch 1) NO Χ YES Χ Х X Syncmaster present: NO Free cassette (without cassette present interlock) - Bucky or tomo synchronous contact 2) YES - WA/WB 4:1-2 (20/21) Exposure request instantly with preparation 1) Exposure switch type: Single step Individual preparation request and exposure request Double step = Density correction in steps of 6%. Range: -8 ... +8 Bucky format density correction: Correction during collimation, input at WAX11/12 pin 1-2 or WBX11 pin 1-2, side fields active when contact closed (<24x24cm), center field only when opened. With GALILEO or NICOL collimator via Bucky controller with CAN. Cone density correction: Density correction in steps of 6%. Range: -8 ... +8 Input at WBX13 pin 4-5. Valid only when the assigned cone contact is closed. Dose measurement input: Measuring chamber respectively at input EZX 21, 22, 31, 32, 41 No measuring chamber assigned. none For free cassette or tomography without TDC.

1) For this function the WA/WB option is required. Ignore for RGDV with CAN: BuckyDiagnost TH/TH2, Thoravision, DigitalDiagnost, EasyDiagnost with Bucky unit.

2) All RGDV with CAN: BuckyDiagnost TH/TH2, Thoravision, DigitalDiagnost, EasyDiagnost with Bucky unit.

Data set A continued:

Dose measurement sensor type:

Bucky AMPLIMAT: Input via EZX21/22/31/32/41, measuring field selection on control desk possible.

Scopo AMPLIMAT: Input via EZX21/22/31/32/41, only center field displayed on the control desk.

It has to be programmed in case of a "Junior Diagnost" or "Extremities"

measuring chamber to prevent from side field selection.

Photo sensor: Input via EZX61, if "none" is programmed at "dose measurement input", in this

case the exposure ON signal has to be connected to the central part via

signal bus.

Photo sensor / AMPLIMAT input: Photo sensor input via EZX21/22/31/32/41.

Only EZX41 supports DSI exposure series with non AEC techniques.

Exposure series /

Tomo movement: NO = Instant brake after exposure end.

NO = To be programmed for tomo systems via system CAN.

BuckyDiagnost TH/TH2, DigitalDiagnost,

EasyDiagnost with Bucky unit

YES = More than one exposure possible with one PREP.

- For tomo units released by 1WA/2WA, PREP must be kept active at the release decade to get the tomo stand back

to the start position.

- DSI or multiple exposures for SCOPO.

Release delay (automatic techniques): enable = Must be enabled for all AEC techniques.

Automatically disabled if non AEC techniques

selected.

disable = Not to be programmed.

Keep release delay always on enable.

Data set A continued:

Mounted radiographical controller:

Must be programmed if any release circuit adaptation none =

> unit 1WA, 2WA, 1WB is assigned to this RGDV. No CAN controlled system is assigned to this RGDV.

Bucky controller 1 / DigitalDiagnost CAN controlled system

- BuckyDiagnost TH/TH2, DigitalDiagnost

RGDV1 ... 4 only

- EasyDiagnost with Bucky controller

RGDV5 ... 8 only

Bucky controller 2 No function yet.

Thoravision Can controlled system

> - Thoravision RGDV1 only

Assignment of the release unit 1WA, 2WA, 1WB Release circuit adaptation unit:

> none free cassette or in case of a CAN driven examination unit

Mounted tomo extension: no tomo unit installed none

> 1WA (1)WAX21 valid as tomography time input 2WAX21 valid as tomography time input 2WA

Medium II format kV corr. (dose equiv. steps): Range = 0...8 dose equiv. + kV correction steps

Medium II format density corr. (-6% steps): Range = 0...8 -6% density steps Medium II format mAs corr. (-6% steps): Range = 0...8 -6% mAs steps

Small II format kV corr. (dose equiv. steps): Range = 0...8 dose equiv. + kV correction steps

Small II format density corr. (-6% steps): Range = 0...8 -6% density steps Small II format mAs corr. (-6% steps): Range = 0...8 -6% mAs steps

An RGDV must not be assigned to a "mounted radiographical controller" and a "release circuit adaptation unit" together.

	Data Set B (underlined items = default)				
Used for tomo:	NO YES	=		of the tomography time e.g. via WAX21 or from a	•
Used for fluoroscopy:	NO YES	=	Ok Fluoro unit		
CT add on:	<u>NO</u>	=	Ok (No function avail	able yet)	
Disable time override:	NO YES	=	Disables time overrid Automatically disable	e function at desk. d with "Used for tomo =	YES".
Tube power factor:			1 <u>100%</u>		
kV steps:			•	kV-grading in steps of 6% kV steps ≙ 25% d	
mA steps:			step width in <u>25</u> , 12 o	or 6%	
Time steps:			exposure time step width in 25, 12 or 6%		
Density steps:			step width in 25, <u>12</u> or 6%		
Density correction (6% steps)	:		$-8 \dots \underline{0} \dots +8$ correction steps For correction see chapter 13.		
Underexposure display (non-automatic techniques):	YES YES NO	= =	Underexposure is also indicated with techniques without AMPLIMAT. To be programmed for tomo systems via system CAN: BuckyDiagnost TH/TH2, DigitalDiagnost EasyDiagnost with Bucky unit. Must be set for all non-CAN tomo systems.		
Tube overload protection:	ON OFF	= =	Overload protection active (default): red = no PREP possible Exposures are possible irrespective of load status. Must not be programmed.		
			desk display	tube load	
			green	100%	
			green - yellow	100%	
			yellow	80%	
			yellow - red	64%	

ALL RIGHTS RESERVED

red

0%

8.2.2. Interface assignment

· Select menu:

PROGRAM/ REGISTRATION DEVICES/ RGDV INTERFACE ASSIGNMENT/ BUCKY/TOMO 1WA...2WA

Note

There must be no programming here if the diagnostic unit is connected via the CAN interface: BuckyDiagnost TH/TH2, Thoravision, DigitalDiagnost, EasyDiagnost with Bucky unit.

 Assign the format and ready contacts of the decade connector WAX11 or WAX12 to a Bucky or tomography RGDV. Refer to Z1-15.1.

- Decade Bucky 1 (X11) See following table.

- Decade Bucky 2 (X12) Program the functions as for the first Bucky decade.

- Tomo Time 0.1 ... 6000ms for each trajectory.

One tomography unit can be programmed for each device interface.

Decade Bucky 1 2				
Tomo mode switch:	disable =	Input "tomo mode" is not activated. Remote changeover Bucky / tomography not possible via the examination unit.		
	enable =	Input "tomo mode" is activated. Remote changeover Bucky / tomography possible. Bucky and tomo RGDV must be defined.		
Bucky RGDV - switch related	: none/ RGD	V 1 8		
	The inputs "format contacts" and "Bucky ready" are activated.			
	When the tomo mode switch is enabled, this RGDV is activated when the remote tomo mode switch is open.			
Bucky RGDV:	none/RGD\	/ 1 8		
	The inputs '	format contacts" and "Bucky ready" can be assigned to any		

RGDV button.

The inputs "format contacts" and "tomo ready" are activated.

Tomo RGDV - switch related: none/ RGDV 1 ... 8

When the tomo mode switch is enabled, this RGDV is activated when the

remote tomo mode switch is closed.

· Select menu:

PROGRAM/ REGISTRATION DEVICES/ RGDV INTERFACE ASSIGNMENT/ BUCKY/ SCOPOMAT 1WB

Bucky RGDV: none or

RGDV 1 ... 8 activates inputs at WBX 11 for RGDV:

- X11:9-10 R/F controller ready-closed.

- X11: 1 - 2 Format size correction contact (Bucky chambers only).

- · RESET the generator.
- Fill in the program settings in table "RGDV programming" 2Z-2.0 at the end of this section.

8.2.3. Examples for RGDV programming

Example No.	System	Refer drawing
1	Diagnost 76: Exposure Scopo / BV-DSIBucky wall standFree cassette	2Z-2.1
2	BuckyDiagnost TH with Bucky controllerDiagnost 76: Exposure Scopo / BV-DSI	2Z-2.2
3	Diagnost 96: Cassette + DSI with / without tomoBucky wall standFree cassette	2Z-2.3
4	HDH with / without tomo time inputBucky wall standFree cassette	2Z-2.4
5	 BuckyDiagnost TH/TH2, Digital Diagnost with Bucky controller and CAN interface Bucky wall stand Free cassette 	2Z-2.5
6	BuckyDiagnost TSBucky wall standFree cassette	2Z-2.6
7	 EasyDiagnost: Exposure Scopo / BV-DSI CAN controlled CS Bucky unit Bucky wall stand Free cassette 	2Z-2.7
8	 Bucky TH any version with Bucky controller Generator equipped with or without decade adaptation unit WA Auxiliary for MCS (only) = RGDV4 in combination with free cassette 	2Z-2.8
9	 Bucky TH any version Generator equipped with or without decade adaptation unit WA Auxiliary for MCS (only) = RGDV5 8 	2Z-2.9
10	 Bucky TH any version Generator equipped with or without decade adaptation unit WA Auxiliary for Trauma Diagnost (only) = RGDV5 8 Auxiliaries RGDV1 4 must not be used with a Bucky TH system via CAN 	2Z-2.10

8.3. **Tube adjustment**

8.3.1. Tube conditioning



Warning!

Radiation is released during the conditioning procedure!

8.3.1.1. Preconditions / Programmings

· Switch OFF the generator.

Preparation of generators which are connected via a CAN interface:

- BuckyDiagnost TH and TH2
- DigitalDiagnost
- Thoravision
- EasyDiagnost with Bucky unit
- · Switch OFF the generator.
- · Disconnect the following plugs:

System	Connector		
	EZX23 signal bus	EZX42 or EZX42-1 system CAN	EZX43 or EZX43-1 system CAN
BukyDiagnost TH / TH2	X		Х
DigitalDiagnost	Х	Х	Х
Thoravision	Х	Х	Х
EasyDiagnost with Bucky unit	Х	Х	Х

• Switch ON the generator.

Note

The programming procedure must not be started before relay ENK1 has been energized at least 2 minutes after the generator has been switched ON.

· Perform the following programmings temporarily for each tube connected to one of the assigned RGDV = Free cassette

Select menu XRGSCOPE:

PROGRAM/ REGISTRATION DEVICES/ RGDV#/ DATA SET A

Programming	Temporarily	Original Tube
Enable handswitch	YES	
Syncmaster present	NO	
Exposure switch type	Double Step	Verify the customized entries in 2Z-2.x
Exposure series / Tomo	YES	
Mounted radiographic	NONE	

- · RESET the generator.
- Select appropriate programmed RGDV = "Free cassette" for the tube to be conditioned.

8.3.1.2. Procedure

· Select large focus only.

Note

The generator must be in the READY state.

• Run reconditioning procedure for an adapted tube, refer to following table, left column TUBE ADAPTED.

or

- Run conditioning procedure for a new or non-adapted tube, refer to following table right column TUBE NOT ADAPTED.
- · It is recommended that the high tension be monitored during conditioning.

Connect the scope:

Channel1: kV AV HT at EZ130 X3 (1V/div), scale: 20kV/V Trigger external: CTRL_X_C/ at backpanel EZ X74, negative slope

Time base: 2ms/div

• In case of problems like tube arcing see the following flowchart EXPOSURE SEQUENCE as an example. The flowchart applies for applicable kV range only, e. g.:

109kV is the max. kV value for normal application, perform just up to next higher kV step = 117kV.

Note


Refer to flowchart EXPOSURE SEQUENCE.

If the tube arcs at a certain kV value, switch another 3 exposures with same parameters and 10s pause between subsequent exposures. In case of success (no arcing anymore) continue with next kV step of the following table.

If the last exposure still arcs go one kV step back and follow normal procedure. If this routine has been performed three times without improvement: ==> Replace the tube!

	Exposure parameters for conditioning				
	Tube adapted		# exposures	Tube no	t adapted
kV	mA	ms		kV	mAs
80	10	50	<1>	80	0.5
80	10	500	<1>	80	5
80	200	250	<1>	80	50
	10 seconds pause			10 seco	nds pause
80	max. mA	100	<1>	80	100
	1 minute pause			1 minu	te pause
90	10	50	<1>	90	0.5
90	10	500	<1>	90	5
90	200	250	<1>	90	50
	10 seconds pause			10 seco	nds pause
90	max. mA	100	<1>	90	100
	1 minute pause			1 minu	te pause
100	10	50	<1>	100	0.5
100	10	500	<1>	100	5
100	200	250	<1>	100	50
	10 seconds pause			10 seconds pause	
100	max. mA	100	<1>	100	100
	1 minute pause			1 minute pause	
110	10	50	<1>	110	0.5
110	10	500	<1>	110	5
110	200	250	<1>	110	50
	10 seconds pause			10 seco	nds pause
110	max. mA	100	< 1 >	110	100
	1 minute pause			1 minu	te pause
120	10	50	<1>	120	0.5
120	10	500	< 1 >	120	5
120	200	250	< 1 >	120	50
	10 seconds pause			10 seco	nds pause
120	max. mA	100	<1>	120	100
1 minute pause			1 minu	te pause	
130	10	50	<1>	130	0.5
130	10	500	<1>	130	5
130	200	250	<1>	130	50
	10 seconds pause			10 seco	nds pause
130	max. mA	100	<1>	130	100
	1 minute pause			1 minu	te pause

	Expos	ure parame	eters for condition	ning	
Tube adapted		# exposures	Tube no	t adapted	
kV	mA	ms	-	kV	mAs
140	10	50	< 1 >	140	0.5
140	10	500	< 1 >	140	5
140	200	250	< 1 >	140	50
	10 seconds pause			10 secor	nds pause
140	max. mA	100	< 1 >	140	100
	1 minute pause			1 minut	te pause
145	10	50	< 1 >	145	0.5
145	10	500	< 1 >	145	5
145	200	250	< 1 >	145	50
	10 seconds pause			10 seconds pause	
145	max. mA	100	< 1 >	145	100
1 minute pause			1 minute pause		
148	10	50	<1>	148	0.5
148	10	500	< 1 >	148	5
148	200	250	< 1 >	148	50
	10 seconds pause			10 secor	nds pause
148	max. mA	100	< 1 >	148	100
	1 minute pause			1 minut	te pause
150	10	50	< 1 >	150	0.5
150	10	500	< 1 >	150	5
150	200	250	< 1 >	150	50
	10 seconds pause			10 secor	nds pause
150	max. mA	100	< 1 >	150	100
	1 minute pause		1	1 minut	te pause



Note

If a tube arcs at any kV value which is not required for application the max. kV (e.g.117kV) program this new limit value by XRGSCOPE:

PROGRAM/ TUBES/ TUBE LIMITS/ MAX. TUBE VOLTAGE LIMIT [kV]/ [117]

As the limit value decreases for this reason, a following re-adaptation procedure sets the field ADAPTED TO [kV] to this value as well.

- Set RGDV programming to original status if no adaptation procedure has to be executed.
- · RESET the generator.

8.3.2. Tube adaptation



Warning!

Radiation is released during the adaptation procedure!

8.3.2.1. General information

Tube adaptation is an automatic process which includes:

- 1. The measurement of the mA offset value that is caused by:
 - the kV measuring circuit
 - the emission current feedback circuit (VCO)
- 2. The measurement of the individual standby filament current (based on 100µA).
- 3. The emission current characteristic as f (kV, filament current).
- 4. The dynamic behavior (positive and negative boost adaptation) where the inertia of the filament with respect to heating up and cooling down is registered.

For more information refer to section 3: FAULT FINDING.

Note

In case of problems check the symptom / solution list at the end of this adjustment chapter. Repeat the adaptation for this particular focus.

8.3.2.2. Preconditions / Programmings

· Switch OFF the generator.

Preparation of generators which are connected via a CAN interface:

- BuckyDiagnost TH and TH2
- DigitalDiagnost
- Thoravision
- EasyDiagnost with Bucky unit
- · Switch OFF the generator.
- · Disconnect the following plugs:

System	Connector		
	EZX23 signal bus	EZX42 or EZX42-1 system CAN	EZX43 or EZX43-1 system CAN
BukyDiagnost TH / TH2	Χ		X
DigitalDiagnost	Х	Х	X
Thoravision	Х	Х	X
EasyDiagnost with Bucky unit	Х	Х	X

· Switch ON the generator.

The adaption procedure must not be started before relay ENK1 has been energized at least 2 minutes after the generator has been switched ON.

- The tube must be conditioned as described in chapter 8.3.1 "Tube conditioning".
- · Check the upper kV limit

Select menu XRGSCOPE:

PROGRAM/ TUBES/ TUBE LIMITS/ MAX. TUBE VOLTAGE LIMIT [kV]

The programmed value should match the nominal value of the tube connected or in case of older tubes the upper kV limit should be set to the max. application kV.

Once an adaptation is completed the new limit value indicates as ADAPTED TO [kV].

 Perform the following programmings temporarily for each tube connected to one of the assigned RGDV = Free cassette

Select menu XRGSCOPE:

PROGRAM/ REGISTRATION DEVICES/ RGDV#/ DATA SET A

Programming	Temporarily	Original Tube
Enable handswitch	YES	
Syncmaster present	NO	
Exposure switch type	Double Step	Verify the customized entries in 2Z-2.x
Exposure series / Tomo	YES	
Mounted radiographic	NONE	

8.3.2.3. Procedure

· RESET the generator.

• It is recommended that the high tension be monitored during adaptation.

Connect the scope:

Channel1: kV AV HT at EZ130 X3 (1V/div), scale: 20kV/V CTRL X C/ at backpanel EZ X74, negative slope Trigger external:

Time base: 2ms/div

• Select the RGDV = Free cassette for the tube to be adapted.

 Select menu XRGSCOPE: OPTIMUS XRG/ ADJUST/ TUBE ADAPTATION

• Select the tube and focus to be adapted, start with small focus!

Note

To avoid any malfunction make sure that READY is displayed on the desk before transmitting data by pushing <F2>

READY state disappears, ADAP is displayed on the desk, WAITING is displayed on the PC screen. Wait until the generator turns back to READY state.

Start the adaptation process by pushing the handswitch in PREP and EXP position continuously.

The generator switches about 125 exposures for each focus. The radiation sign at the desk indicates exposures but there is no beep at the end of each exposure.

The actual kV parameters are displayed during adaptation.

The generator carries out the adaptation automatically. The procedure for one focus is completed when the desk indication changes from ADAP to TEST. The WAITING message disappears from the PC screen together with a PC beep, followed by the screen: BEFORE CONTINUING THE GENERATOR MUST BE RESET.

- RESET the generator.
- Run the adaptation for each focus (small and large) and tube.

Note

As there is no tube type with a physical third (middle) focus yet, the third focus can not be adapted. VARIOFOCUS values are calculated by adapted small and large focus. APR programs using VARIOFOCUS can only be selected until small and large focus are both adapted.

 Set RGDV(s) programming to original status according to table "RGDV programming" 2Z-2.x in the end of this chapter.

2-52 OPTIMUS R/F (c/02.1)OPTIMUS_RF_2_c021

8.3.3. Final tube adjustment work

- 1. BuckyDiagnost TH with CAN interface, DigitalDiagnost, Thoravision:
 - Switch OFF the generator.
 - Reestablish signal bus connector EZX23.
 - Reestablish CAN connectors EZX42-1 and EZX43-1.
 - Switch ON the generator.
- 2. All other systems:
 - · Reset the generator.

8.3.4. Problems during adaptation - Symptoms and solutions

1. A warning cannot be displayed on the control desk, the WAITING screen on the PC is flickering instead during this event and logged in the error log index.

2. If the tube has already been at a high temperature level (but the tube load indication still indicates green or green-yellow for 100% power) it might happen that the load indication changes straight to red and the adaptation is on hold. WAITING is flickering on the PC.

Solution: Keep the handswitch pushed, once the temperature is down adaptation continues automatically.

Note

An increment of one of the temperature levels inhibits the 100% power condition. This event is always logged as warning 00BV in the error log index.

- 3. An error message just flashes for a very short moment and is instantly covered by ADAP on the desk afterwards. WAITING is flickering on the PC.
- 4. All buttons at the control desk including the RESET button are inactive during adaptation. The only way to RESET an error is to release the PREP switch which causes an interrupt similar to the RESET command.
- 5. After letting go of the PREP switch wait until the desk indicates READY. If READY does not appear at least after 20 seconds run a warmstart of the generator by pushing the RESET button on CU EZ139.
- 6. If adaptation seems to do nothing for more than 30 seconds let go of the PREP switch. Wait until the desk indicates READY. If READY does not appear at least after 20 seconds run a warmstart of the generator by pushing the RESET button on CU EZ139.
- 7. If a constant READY indication appears for more than 2 seconds while PREP and EXP is activated by the handswitch during adaptation let go of the handswitch. Wait until the desk indicates READY. If READY does not appear at least after 20 seconds run a warmstart of the generator by pushing the RESET button on CU EZ139.
- 8. If adaptation does not carry on with or without READY indication check whether one of the function units indicates a FATAL error by turning on the red LED. Let go of the handswitch and run a warmstart of the generator by pushing the RESET button on CU EZ139.
- 9. If adaptation has been interrupted by a generator warmstart check the error log index before restarting adaptation:
 - kV errors 02WG and/or 02WH indicate tube arcing. Run conditioning of the tube as described in chapter 8.3.1 and/or reduce the max. kV value to the required application value.

2-54 (c/02.1)OPTIMUS R/F OPTIMUS_RF_2_c021

8.4. Dose rate control

8.4.1. AMPLIMAT sensitivity

· Select menu:

PROGRAM/ DOSE RATE CONTROL AMPLIMAT/ SENSITIVITY

• Depending on HW programming of jumper EZ150:W4. W4 programs sensitivity accordingly:

```
high = \times 4 = EZ150:W4 in position 1
===> All screen/film combinations with a system speed > 200.
low = \times 1 = EZ150:W4 in position 3
===> At least one screen/film combination with a system speed \leq 200.
```

8.4.2. Screen/film combinations

5 screen/film combinations can be programmed for each of the 5 measuring chambers:

· Select menu:

PROGRAM/ DOSE RATE CONTROL/ AMPLIMAT CHAMBER 1 ... 5/ DATA SET 1 ... 5

The number of the chamber corresponds to the specified unit number of the dose measuring unit.

The choice between automatic and manual DRC processing is possible when an authorized hardware key is inserted in the PC.

Automatic is selected as default and must be used for the initial programming. Data sets of adjacent rooms can be copied but have to be aligned afterwards.

Access manual DRC processing by pushing the <ESC> key.

The manual mode is suitable for:

- Copying complete programming to other measuring chambers.
- Setting the basic density.
- Changing the desk-displayed names of the programmed screen/film combinations.
- Creating backups of the DRC programmings.

8.4.2.1. Automatic DRC processing

Select the desired data from the files offered for the following programming steps.

The files are part of the installation software.

- Select the programming field with the cursor and enter < **RETURN>**.
- Enter the desired file from the list offered.
- Select the desired data as required.

CORRECTION FACTOR:

FILM: File FILM.TDL Film types according to description of the manufacturer. File FILM_BL/ _GR/ _UV.TDL General classification of the film according to color, sensitivity S and RLF compensation. SCREEN: File SCREEN.TDL Screen types according to description of the manufacturer. File **LUMAT LG.TDL** Screen types according to luminous matter. - Imaging plates. **CHAMBER:** File CHAMBER.TDL Different types of measuring chambers. File CASSETTE.TDL **CASSETTE:** Different types of cassettes. SYSTEM CORRECTION: File SYSCOR.TDL Select no corr. (ISO 9236-1)

Based on the combination of the components entered, the processor calculates the switch-OFF dose, kV correction and RLF compensation. The name for the screen/film combination, e.g. "B400", is taken from the "screen" default data set.

Correction factor for switch-OFF dose.

Dose Rate Control setting Optimus for Computed Radiography (PCR or other imaging plates)

The following example is for a 400 speed system, determined by the selection of the LG06 400 speed type from file LUMAT LG.TDL (luminous groups).

Ignore the violet screen colour of LG06, the data set just requires its kV characteristic.

Default 1.00

FILM: File FILM.TDL X-CONSTANT RLF=1

SCREEN: File **LUMAT_LG.TDL** LG06 S400 vi

CHAMBER: File CHAMBER.TDL the installed chamber type CASSETTE: File CASSETTE.TDL normal cassette (def) File SYSCOR.TDL **SYSTEM CORRECTION:** no corr. (ISO9236-1)

CORRECTION FACTOR: 1.00

Note

Film, screen, etc. data selected are not directly stored in the generator. It is recommended that they be entered in the table "Data sets of chambers" 2Z-4 at the end of this section.

RESET the generator.

Color and sensitivity class of the screen/film combination are displayed on the desk, e.g.: "B400". Other screen/film combinations (data set 1 ... 5) for the chamber can be selected by the \pm buttons.

8.4.2.2. Manual DRC processing

The current data set of the screen/film combination is displayed.

* Abbreviation: Abbreviation for the screen/film combination.

Example: B400 = blue, speedclass 400.

Dose Request Chamber: Senstitivity of the measuring chamber type in $[\mu Gy/V]$.

* Dose of FSC: Switch-OFF dose of the screen/film combination in $[\mu Gy]$.

Linear ratio with respect to the film density.

kV70-Char. U 0 ... 9: Checkpoints for kV-dependent density correction.

kV70-Char. Drel 0 ... 9: Relative correction value for the dose.

RLF t_0 ... 9: Checkpoints for time-dependent density correction.

(RLF = Reciprocity Law Failure).

RLF Drel 0 ... 9: Relative correction value for the dose.

* = Only these fields may be changed according to the system requirements.

All other fields must not be changed.

• If required, change the data and the abbreviation name.

Usually no value except the basic density "Dose of FSC" must be changed (see next page).

- · Transmit the data set with <F2>.
- · RESET the generator.

The **SAVE** <**F3**> and **LOAD** <**F4**> functions of **XRGSCOPE** permit straightforward copying of the measuring chamber programmings.

8.4.2.3. Density correction for AEC technique (option)

Basic density per screen/film combination:

A hardware key is required at the PC for direct access to the switch-OFF dose.

- · Make a test exposure for each screen/film combination. To do so, set the density correction = 0.
- · Determine the density of the test exposures.
- · Select menu: PROGRAM/ DOSE RATE CONTROL/ AMPLIMAT/ CHAMBER 1 ... 5/ DATA SET 1 ... 5
- Select manual DRC programming by pushing the <ESC> key.
- Correct the switch-OFF dose = "Dose of FSC" according to formular below:

- Transmit the data set by pushing the <F2> key.
- · Repeat the procedure for each screen/film combination at each chamber.
- · RESET the generator.

The switch-OFF dose can be set on the PC even without a hardware key.

To do so, call up the automatic DRC programming, repeat all the selections and change the correction factor for switch-OFF dose accordingly. Each time this programming is called up all the selections must be repeated.

2-58 OPTIMUS R/F (c/02.1)OPTIMUS_RF_2_c021

8.4.3. Image intensifier (II)

- II lead time
- Density voltage correction: Default = 1V
- Dynamic factors

At auxiliaries which are using DSI tomo with TDC (Tomo Density Control) take care that the following settings are present to get a linear density voltage of 1V.

· Select menu:

PROGRAM/ REGISTRATION DEVICES/ RGDV x/ DATA SET A

Dose measurement input: EZX41

Dose measurement sensor type: Photo sensor / Amplimat input

· Select menu:

PROGRAM/ DOSE RATE CONTROL/ AMPLIMAT/ CHAMBER 5/ DATA SET 1

• Press <ESC> button.

Checkpoint	Correction factor	Remark
Abbreviation	[def1] A	A
Dose Request Chamber [μGy/V]	[6.40]	△ Note Do not care about these correction factors.
Dose of FSC [μGy]	[2.14]	
kV70-Char. U_0 [kV]	[40]	
kV70-Char. Drel_0	[1.00]	
kV70-Char. U_1 [kV]	[40]	
kV70-Char. Drel_1	[1.00]	
kV70-Char. U_2 [kV]	[50]	
kV70-Char. Drel_2	[1.00]	
kV70-Char. U_3 [kV]	[60]	
kV70-Char. Drel_3	[1.00] *	
kV70-Char. U_4 [kV]	[70]	
kV70-Char. Drel_4	[1.00]	
kV70-Char. U_5 [kV]	[80]	These fields of the kV dependent correction factors must always be 1.00.
kV70-Char. Drel_5	[1.00] *	
kV70-Char. U_6 [kV]	[90]	
kV70-Char. Drel_6	[1.00] *	
kV70-Char. U_7 [kV]	[110]	
kV70-Char. Drel_7	[1.00] '	
kV70-Char. U_8 [kV]	[130]	
kV70-Char. Drel_8	[1.00]	
kV70-Char. U_9 [kV]	[150]	
kV70-Char. Drel_9	[1.00] *	

Checkpoint	Correction factor		Remark
RLF t_0 [ms]	[0]	Δ	
RLF Drel_0	1.000]	Δ	
RLF t_1 [ms]	[20]	Δ	
RLF Drel_1	[1.000]	Δ	
RLF t_2 [ms]	[60]	Δ	
RLF Drel_2	[1.000]	Δ	
RLF t_3 [ms]	[100]	Δ	
RLF Drel_3	[1.000]	Δ	
RLF t_4 [ms]	[500]	Δ	Δ Note
RLF Drel_4	[1.000]	Δ	Do not care about these RLF correction factors.
RLF t_5 [ms]	[1000]	Δ	
RLF Drel_5	[1.000]	Δ	
RLF t_6 [ms]	[1500]	Δ	
RLF Drel_6	[1.000]	Δ	
RLF t_7 [ms]	[2000]	Δ	
RLF Drel_7	[1.000]	Δ	
RLF t_8 [ms]	[3000]	Δ	
RLF Drel_8	[1.000]	Δ	
RLF t_9 [ms]	[4000]	Δ	
RLF Drel_9	[1.000]	Δ	

8.4.4. Fault exposure detection

Fault exposure detection is switched ON as a default for AEC and TDC. If in the initial phase of an exposure too little dose is measured, the exposure is aborted to protect the patient.

- Time of control measurement: 10% of backup time,

min. 250ms at TDC

- Dose minimum: 4% of set density voltage at AEC,

4 ... 10% at TDC

- Backup time AEC: Calculated time from 9.5 times mAs of the respective 2-factor technique,

 Backup time TDC: Exposure time set 0.3 ... 6s

This additional precaution can be switched OFF for both techniques individually in the menu: PROGRAM/ DOSE RATE CONTROL/ FAULT EXPOSURE DETECTION/ AEC OR TDC

For details see section "FAULT FINDING", chapter "Optimus AEC switch-OFF philosophy".

This monitoring does not take effect in the following cases, irrespective of programming:

- Using screen/film combinations with high speed in AEC technique.
- Exposure time in TDC technique is shorter than 1s.

8.4.5. Continuous fluoroscopy

Factors for continuous fluoroscopy can be modified: (see also section 4 FAULT FINDING, "Explanation of programming")

scantime TV [ms]: = 20.00ms (default) for scantime of the TV system ≤ 20ms (50Hz and 60Hz)

= Scantime of the TV system for scantime > 20ms

- scantime TV valid : = YES default

= NO no function vet

- P max EDL [W] : = 250W (default) (EDL = Entrance Dose Limiter)

(range = 0 ... 9000W) The maximum output of the tube during continuous fluoroscopy and with the SID

signal prevailing is limited to the programmed wattage.

P max EDL [W] is only active if the DRLM (Dose Rate Limiter) signal on 1 WBX12:1

is active (see also Z1-15.2)

For adjustment see section 6 ADJUSTMENTS.

- TV pos limit [V] : = 6.8V (default) Adjustment necessary if the fluoroscopic image does not come up

properly.

See section 6 ADJUSTMENTS. (range = +3 ... +7.5V)

TV neg limit [V] : = -6.8V (default) (range = -7.5 ... -3V)

Application limits 8.5.

8.5.1. X-mode limits

Using the menu:

PROGRAM/ APPLICATION LIMITS/ X-MODE LIMITS

Limit values can be defined for all available techniques. Some values look as if they are out of limit which they are indeed, but there are additional basic limit values programmed in the generator firmware. These are exposure technique dependent.

As an example the field of the "Falling Load" technique:

X-ray Mode:	AEC falling load kV
Min. Time Limit [ms]:	[1.00]
Max. Time Limit [ms]:	[60000.00]
Min. Current Time Product Limit [mAs]:	[0.001]
Max. Current Time Product Limit [mAs]:	[580.000]

Min. Time Limit [ms]: Is always 1ms for all non-AEC (Automatic Exposure Control) techniques. Exposures with AEC might be switched shorter than 1ms.

Max. Time Limit [ms]: Basic limits are technique dependent and can not be changed or increased:

AEC falling load	kV	4000ms
AEC fixed current	kV-mA	4000ms
TDC (Tomo Density Control)		6000ms
Γ	kV-mA-ms	16000ms
free techniques	kV-mAs	16000ms
	kV-mAs-ms	16000ms

Min. Current Time Product Limit [mAs]: The smallest mAs - product is 0.5mAs.

AEC exposures with less than 0.5mAs are possible.

Max. Current Time Product Limit [mAs]: The default mAs - product is 580mAs for all AEC-techniques.

850mAs is the absolute limit the generator terminates.

Note

Local limits have to be taken into consideration.

8.5.2. Thoravison limits

The kV-dependent mAs limits can be accessed via the menu: PROGRAM/ APPLICATION LIMITS/ **THORAVISION LIMITS**

They are activated only in conjunction with an online Thoravision unit. A change may only be made if instructed to do so by the service center.

Reference files on floppy disk: - ref limx.tdl X-ray limits

Reference files on system disk Thoravision: - TH_LI128.tdl

- TH_LI64.tdl

- TH LI32.tdl = ref limt.tdl

8.5.3. Overload-dependent limits

- max. current for continuous fluro [mA]: Range = 3 ... 30mA

- max current during overload for continuous fluro [mA]: Range = 1 ... 3mA

8.6. **Human interface**

Note

The Optimus generator of this system might has been delivered without APR

Before spending time trying to load APR check:

XRGSCOPE

>> Faultfind

>> Power ON Results

>> APR disabled: [No] = APR option present APR disabled: [Yes] = No APR possible

If APR shall be possible order

MGR 0011 (change of an existing configuration)

MGR 2181 **APR** with the generator serial number.

A maximum of up to 1024 APRs can be stored in the generator.

On a single RGDV button either up to 80 APRs can be programmed directly (10 pages of 8 each) or up to 400 APRs via menus.

The initial data sets are called ### APR name ### and they all have the same exposure parameters.

They can be directly assigned or via menu and submenu levels to registration devices RGDV 1 ... 8.

In case TEST APR is displayed after selection of a registration device, at least this particular registration device has not been assigned to any APRs.

8.6.1. Additional program settings for generators without option APR

In case of a generator without option APR 9890 000 0282x, the menu/APR display field must be dark until selection of all exam./aux. units.

If it does not work AR_EMPTY.tdl must be loaded from the service disk for each exam./aux. unit.

In case TEST- APR appears in the menu/APR display field upon selection of the exam./aux. units, the AR EMPTY.tdl can be loaded immediately because no APR/T is assigned to the exam./aux. unit.

In case another APR or menu appears on the desk it must be deleted before the AR_EMPTY.tdl can be loaded:

• Delete the APR/T under each exam./aux. unit.

Select menu:

PROGRAM/ HUMAN INTERFACE/ RGDV RELATED ASSIGNMENT/ RGDV 1 ... 8/ MANUAL ASSIGNMENT/ **DELETE MENU**

After a RESET of the generator TEST APR appears in the menu/APR display field for each exam./aux. unit.

Now the AR EMPTY.tdl can be loaded from the service disk for all exam./aux. units.

Load AR_EMPTY.tdl for all exam./aux. units.

Select menu:

PROGRAM/ HUMAN INTERFACE/ RGDV RELATED ASSIGNMENT/ RGDV 1 ... 8/

PREDEFINED ASSIGNMENTS

A list of possible TDL-files appears on the screen.

· Select file "AR EMPTY.tdl" from the list.

Only "none APR" appears.

Extended APR data table																		
No	APR Name	Tomo no.	Foc	V [%]	kV	mA	ms	mAs	Dose [μGy]	Dens.	Fields I m r	Technique pref background		Spectral Filter	FSC	Fluo. Curve	Menu Name	Submenu Name
1	buckytable	1	1		77	600.0	53.3	32.0	2.50	0	m	RAEC	RUQ	none	400	1		
2	tomography	1	1		77	16.0	2000.0	32.0	2.50	0	m	RAEC	RUQt	none	400	1		
3	wallstand	1	1		77	600.0	53.3	32.0	2.50	0	m	RAEC	RUQ	none	400	1		
4	free cassette	1	1		77	600.0	53.3	32.0	2.50	0	m	RAEC	RUQ	none	400	1		

Select it for each active exam./aux. unit.

After a RESET of the generator the menu/APR display field is dark for all exam./aux. units.

8.6.2. Language

· Assign the desired language.

Select menu:

PROGRAM/ HUMAN INTERFACE/ SELECT LANGUAGE

A language menu appears:

- English
- German
- French
- Spanish
- · Select the desired language.
- · RESET the generator.

A table lists which characters can be displayed on the control desk and how they can be indicated/entered at the service PC, e.g. for APR names, see drawing 2Z-3 "List of characters".

Certain characters can be generated at the PC only via the decimal code. To do so, push the <ALT> key on the PC and enter the numerical code.

8.6.3. Automatic programming of APRs (option)

Note

Pay attention in the rules of "PC and generator settings to avoid problems during up/downloading of CU complete files" described in chapter 7.1 of this section.

The installation disk contains data files for a complete, typical APR programming in different languages. Standard APR programs for each application can easily and quickly be loaded for each registration device.

· Select menu:

PROGRAM/ HUMAN INTERFACE/ RGDV RELATED ASSIGNMENT/ RGDV 1 ... 8/ **PREDEFINED ASSIGNMENTS**

All files consist of the the same format "A***my**.tdl"

• Select with <TAB> and <UP-/DOWN-ARROWS> one of the files listed:

for generators without APR option 9890 000 0282x AR EMPTY A***** APR data file = ***my** Release month / year *TS..V Bucky TS with VARIOFOCUS RO 1648, RO 1750, SRO 2550 SRO 33 100 = *T..V9 Bucky TS with VARIOFOCUS SRO 0951 = *CO Cosmos = *D DigitalDiagnost *R Radiography *DS DSI, EasyDiagnost *D15/76/96 R/F systems *SI EasyDiagnost, D15 *F Fluoro stands, DSI spotfilm *TR Trauma *TS **BuckyTS** *V VARIOFOCUS RO 1648, RO 1750, SRO 2550, SRO 33 100 = *V9 VARIOFOCUS SRO 0951 Mono focus tube *M *P **Pediatrics** = ***DE(U) German version ***ES(P) Spanish version ***EN(G) English version = ***FR(A) French version Example: A D 91 V 9 EN.tdl **English version** Tube SRO 0951 | VARIOFOCUS Release: September 2001 DigitalDiagnost APR data file

Some files contain data sets in a window:

Available Examination Unit Type

· Select one of the files listed and load it:

BUCKY GR Bucky table APR group layer (menu technique)

BUCKY PA Bucky table APR paging (scrolling through pages technique)

WALLSTAND GR Bucky wallstand group WALLSTAND PA Bucky wallstand paging =

FREE GR Free technique (non-AEC) group FREE PA Free technique (non-AEC) paging =

TOMO LT/HDH Linear tomography with units HDH, BTS2, BTS4 group TOMO LIN/PA Linear tomography with units HDH, BTS2, BTS4 paging

TOMO BTC GR Tomography with unit BuckyDiagnost TC group = Tomography with unit BuckyDiagnost TC paging TOMO BTC PA TOMO BTH GR Tomography with unit BuckyDiagnost TH group TOMO BTH PA Tomography with unit BuckyDiagnost TH paging

- · Repeat this procedure for each registration device.
- · RESET the generator.

Now all APR programs which have been loaded are displayed on the desk.

If required:

- Change the menus and the assignments of the APRs according to chapter 8.6.4.
- Change the names and the contents of the APRs according to chapter 8.6.5.

Note

In case a complete APR program setting has to be replaced by another program setting, all other APR program settings under the registration device concerned must be deleted first. Use APRMAN.

For deleting a program setting call menu:

PROGRAM/ HUMAN INTERFACE/ RGDV RELATED ASSIGNMENT/ RGDV 1 ... 8/ MANUAL ASSIGNMENTS/ **DELETE MENU**


and select the blank line.

8.6.4. Manual programming of APRs (option)

Use APRMAN.

Select service menu:

PROGRAM/ HUMAN INTERFACE/ RGDV RELATED ASSIGNMENT/ RGDV 1 ... 8/ MANUAL ASSIGNMENT/ . . .



8.6.4.1. Creating menus

· Assign the desired menu.

Select menu:

PROGRAM/ HUMAN INTERFACE/ RGDV RELATED ASSIGNMENT/ RGDV 1 ... 8/ MANUAL ASSIGNMENT/ **ASSIGN MENU**

- Enter the first menu name, e.g. <Body region 1>.
- If required, change the location suggested in the display. Otherwise the next vacant location is assigned.
- Enter the second menu name, e.g. <Body region 2>.

Etc.

8.6.4.2. Creating sub menus

· Assign the desired sub menu.

Select menu:

PROGRAM/ HUMAN INTERFACE/ RGDV RELATED ASSIGNMENT/ RGDV 1 ... 8/ MANUAL ASSIGNMENT/ MENU ASSIGNMENT/ SELECT MENU

- Select with the <CURSOR> from one of the windows a menu to be assigned with submenus.
- · Select menu:

PROGRAM/ HUMAN INTERFACE/ RGDV RELATED ASSIGNMENT/ RGDV 1 ... 8/ MANUAL ASSIGNMENT/ MENU ASSIGNMENT/ ASSIGN MENU

- Enter the first submenu name, e.g. <Left side>.
- If required, change the location suggested in the display. Otherwise the next vacant location is assigned.
- Enter the second submenu name, e.g. <Right side>.

Etc.

8.6.4.3. Assigning APRs

· Assign desired APRs.

Select menu:

PROGRAM/ HUMAN INTERFACE/ RGDV RELATED ASSIGNMENT/ RGDV 1 ... 8/ MANUAL ASSIGNMENT/ APR ASSIGNMENT/ SELECT MENU

• Select with the **<CURSOR>** from one of the windows a menu or submenu. If no menu layer is desired, proceed.

· Select menu:

PROGRAM/ HUMAN INTERFACE/ RGDV RELATED ASSIGNMENT/ RGDV 1 ... 8/ MANUAL ASSIGNMENT/ APR ASSIGNMENT/ ASSIGN APR

· Select from one of the windows an initial APR or an APR which has not been assigned.

Initial APR: ### APR name ###
Not assigned APR: e.g. <Thorax ap>

- If required, change the location suggested in the display. Otherwise the next vacant location is assigned.
- · Assign the next APR.

Etc.

· RESET the generator.

Only after a RESET of the generator the menus, submenus and APRs are displayed on the desk.

INSTALLATION OPTIMUS R/F

8.6.5. Changing of APR radiation parameters

- Select the APR to be changed on the desk, e.g. ### APR name ###.
- · Select menu:

PROGRAM/ HUMAN INTERFACE/ APR DATA SET/ SELECT APR DATA SET

The number of the APR selected on the desk is displayed. Transmit data by pushing <F2>.

· Select menu:

PROGRAM/ HUMAN INTERFACE/ APR DATA SET/ CHANGE APR DATA SET

- Change the contents of the APR, e.g. <NAME>, <kV> etc.
- Select the kV value according to the dose-equivalent series: 40 - 41 - 42 - 44 - 46 - 48 - 50 - 52 - 55 - 57 - 60 - 63 - 66 - 70 - 73 - 77 - 81 - 85 - 90 - 96 - 102 - 109 - 117 - 125 - 133 - 125 -141-150 kV.


Parameter	Value	Remark
APR number	1 1024	
APR name	up to 16 characters	
Focus	– Δ small – Δ vario – Δ large	
VARIOFOCUS ratio [%]	20%, 35%, 50%, 65%, 80% of small focus	
Dose measurement field (left)	Δ ON/OFF	AEC operable
Dose measurement field (middle)	Δ ON/OFF	only if at least
Dose measurement field (right)	Δ ON/OFF	one field is ON
Preferred technique	-non automatic -automatic	
AEC technique	- AEC falling load kV - AEC fixed current kV-mA - TDC	
No AEC technique	- kV-mA-ms - kV-mAs - kV-mAs-ms	
Tube current max. factor [%]	1 100	
PSC U thin (dose equiv. steps)	0 5	
PSC U thick (dose equiv. steps)	0 5	
PSC Q thin (6% steps)	0 10	
PSC Q thick (6% steps)	0 10	
PSC dens. thin (6% steps)	0 10 (max. limit = 8)	
PSC dens. thick (6% steps)	0 10 (max. limit = 8)	
Exposure data U [kV]	Δ 40 150	
Exposure data I [mA]	Δ 0.1 2000	
Exposure data Q [mAs]	Δ 0,5 850	
Exposure time [ms]	Δ 1 16000 (60000)	
Exposure data density (6% steps)	Δ -40 +40 (max. limit = ±8)	
Screen/film combination	Δ RGDV dependent	Default = Data Set 1
Tomo no.	1 16 assignment of a tomographic figure	
Spectral filter	- none - 2mm Al - 0.1mm Cu + 1mm Al - 0.2mm Cu + 1mm Al	Default = none
Fluoroscopy curve (for details see figure on next page)	1-20 Select kV/mA = curve 1, 2, 3 Curve 1: 70kV> 1.7mA 110kV 6mA Curve 2: 70kV> 0.8mA 110kV 6mA Curve 3: 70kV> 1.2mA 110kV 3mA	Default = curve 1 Curves loaded from the disk.

The basic setting of these data can also be performed from the desk (RESET + APR). Refer to the operator's manual.

ALL RIGHTS RESERVED

INSTALLATION OPTIMUS R/F

kV	Curve 1 mA	Curve 2 mA	Curve 3 mA
40	0.25	0.25	0.22
50	0.5	0.4	0.42
60	0.95	0.6	0.70
70	1.7	8.0	1.20
80	3	2	2.15
90	5	4	2.70
100	5.5	5	2.95
110	6	6	3.00



kV/mA - curves: Fluoroscopy Optimus R/F

414H01

OPTIMUS R/F INSTALLATION

The following parameters must also be taken into account for AEC techniques:

 mAs: basis for calculating the backup time for AEC, the tube current for AEC fixed current (kV-mA) and the initial mA value for TDC

- t: exposure time for TDC and AEC fixed current

For details see chapters 8.6.5.2 and 8.7.

If AEC FIXED CURRENT kV-mA or TDC is programmed as the preferred technique, the kV/mA/s or kV/mAs/s technique must be selected under NO AEC TECHNIQUE.

- Transmit data by pushing <F2>.
- · Select the next APR on the desk, select it in the programming menu and change it.

Etc.

· RESET the generator.

Organ-dependent correction:

- Select the APR to be changed on the desk.
- · Select menu:

PROGRAM/ HUMAN INTERFACE/ APR DATA SET/ CHANGE APR DATA SET

Transmit the displayed APR number by pushing the <F2> key.

Select menu:

PROGRAM/ HUMAN INTERFACE/ APR DATA SET/ CHANGE APR DATA SET

• Exposure data density: -40 ... +40 = correction in steps of 6%. Max. limit = ±8.

The number of correction steps must match the programmed step rate.

Example:

The desk display has been programmed to 12% step rate in the menu: PROGRAM/ REGISTRATION DEVICE/ RGDV 1 ... 8/ DATA SET B/ **DENSITY STEPS**

To be able to display a density correction of +1 for a certain APR, two corrections steps (2 x 6% = 12%) must be programmed under this APR.

· Select the next APR on the desk, select it in the programming menu and change it.

Etc.

RESET the generator.

Correction for each RGDV 1 ... 8:

This correction is possible but for reasons of clarity it should not be used.

· Select menu:

PROGRAM/ REGISTRATION DEVICE/ RGDV 1 ... 8/ DATA SET B

Density correction: -8 ... +8 = correction in steps of 6%.

INSTALLATION OPTIMUS R/F

8.6.5.1. Density correction for non-AEC techniques

The supplied APR standard sets are mostly based on a screen/film combination with a speed class of 400. APRs for extremities and some other applications are based on a 100 or 200-type system. Depending on the local situation the mAs or s parameters of all the relevant APRs must be adapted.

Example:

The customer uses a 200-type system. To change from the existing 400 values the relevant APRs must be reprogrammed to double the mAs products or to double the exposure time (400 divided by 200 = 2).

- · Select the relevant APR at the control desk.
- Set the new parameters at the control desk.
- Save the new parameters as default values.

To do this push the RESET button and the corresponding APR button. The asterisk in the APR name as an indication of overwritten data disappears. Use APRMAN.

8.6.5.2. AEC fixed current (kV-mA)

For this exposure technique the APRs must have the following programmings:

 Dose measurement field: ON (at least 1 field must be set to ON)

- Preferred technique: automatic

AEC fixed current kV-mA AEC technique:

- No AEC technique: kV-mAs-ms technique (RUQT) or kV-mA-ms (RUIT)

- Exposure data U: = anatomical kV value

- Exposure data Q/I: anatomical mAs/mA product based on the screen / film combination used

- Exposure time t: anatomical exposure time

The mA value is calculated automatically.

In the APR standard files supplied the following APR's is programmed to kV-mA technique:

```
Dens axis F
                 dens axis F
                                   dens axis F
                                                     atlas F
              /
  German
                    English
                                      French
                                                     Spanish)
```

The APR is marked with F.

The mAs value is based on a 400-type screen/film combination and must be adapted to the combinations actually

E. g.: If a 200-type combination is used, the mAs value must be doubled.

If the TDC option is installed, the preferable technique for all exposures is the one where the exposure time is the determining factor. TDC is not restricted to tomography applications only.

ALL RIGHTS RESERVED

OPTIMUS R/F INSTALLATION

8.6.6. Moving / copying of an APR data set

Determination of the number of APR data set X to where APR data set Y is to be moved / copied.

- Select APR data set X on the desk.
- Select menu:
 PROGRAM/ HUMAN INTERFACE/ APR DATA SET/ SELECT APR DATA SET
- Note down the number of APR data set X, e.g. 100.

Changing of the number of APR data set Y to be moved / copied to the number of APR data set X

- · Select APR data set Y on the desk.
- · Select menu:

PROGRAM/ HUMAN INTERFACE/ APR DATA SET/ SELECT APR DATA SET

Transmit by pushing <F2>.

· Select menu:

PROGRAM/ HUMAN INTERFACE/ APR DATA SET/ CHANGE APR DATA SET

- Replace the number of APR data set **Y** with the number of APR data set **X** in the input mask, e.g. nnn ⇒ 100.
- Transmit this number by pushing <F2>.
- · RESET the generator.

APR data set **Y** is displayed in place of the old APR data set **X** on the desk.

In case APR data set **Y** is merely moved and not copied to the location of APR data set **X**, the original APR data set **Y** must be deleted at the end of programming.

INSTALLATION OPTIMUS R/F

8.6.7. Deleting of APRs

· Select menu:

PROGRAM/ HUMAN INTERFACE/ RGDV RELATED ASSIGNMENT/ RGDV 1 ... 8/ MANUAL ASSIGNMENT/ APR DEASSIGN/ **SELECT MENU**

Required only if the APR is assigned to a menu or submenu.

- Select the associated menu or submenu from one of the windows.
- · Select menu:

PROGRAM/ HUMAN INTERFACE/ RGDV RELATED ASSIGNMENT/ RGDV 1 ... 8/ MANUAL ASSIGNMENT/ APR DEASSIGN/ **DEASSIGN APR**

- Select the APR to be deleted from one of the windows.
- · RESET the generator.

An APR which has been deleted is no longer displayed on the desk but remains stored in the generator. It can be re-activated according to chapter 8.6.4.

OPTIMUS_RF_2_c021

OPTIMUS R/F INSTALLATION

8.6.8. Manipulating menus

Deleting:

· Select menu:

PROGRAM/ HUMAN INTERFACE/ RGDV RELATED ASSIGNMENT/ RGDV 1 ... 8/ MANUAL ASSIGNMENT/ **DELETE MENU**

- Select the menu or submenu to be deleted from one of the windows.
 For deleting all APRs of an RGDV program select the blank line.
- · RESET the generator.

Shifting:

· Select menu:

PROGRAM/ HUMAN INTERFACE/ RGDV RELATED ASSIGNMENT/ RGDV 1 ... 8/ MANUAL ASSIGNMENT/ MOVE MENU

- Select the menu or submenu from one of the windows.
- Enter the new positions.
- · RESET the generator.

Renaming:

· Select menu:

PROGRAM/ HUMAN INTERFACE/ RGDV RELATED ASSIGNMENT/ RGDV 1 ... 8/ MANUAL ASSIGNMENT/ RENAME MENU

- Select the menu or submenu from one of the windows.
- Enter the name in the lowermost line.
- · RESET the generator.

INSTALLATION OPTIMUS R/F

8.6.9. External APR assignment

A maximum of 2 external APRT modules can be connected to the Optimus R/F via adapter 4 auxil. units 1WA, 2WA and WB.

Refer Z1-15.1 "Adapter 4 auxil. units".

- Define the external APRT modules by placing links W1 through W3 on adapters WA/WB. Also see table "Assignment of the external APRT modules" below and Z1-15.1 "Adapter 4 auxil. units".
- · Program the desired APRT in the generator.
- Select the desired APRT numbers for the external APRT modules from the RGDVs: PROGRAM/ HUMAN INTERFACE/ APR DATA SET/ SELECT APR DATA SET
- Program the RGDVs that are to be switched over by the external APRT modules. Select menu:

PROGRAM/ HUMAN INTERFACE/ RGDV RELATED ASSIGNMENT/ RGDV 1 ... 8/ EXTERNAL APR ASSIGNMENT/ ...

.... DEVICE INTERFACE 1

or

.... DEVICE INTERFACE 3

Only two RGDVs can be programmed and, respectively, switched over. Once programmed these two RGDVs are valid irrespective of the selected APRs of the other buttons of the module.

• Program the APRs that are to be selected by the external APR/T modules. Select menu:

PROGRAM/ HUMAN INTERFACE/ RGDV RELATED ASSIGNMENT/ RGDV 1 ... 8/ EXTERNAL APR ASSIGNMENT/...

.... **DEVICE INTERFACE 1**: 6 APR keys can be assigned for each RGDV


and/or

.... DEVICE INTERFACE 2 : 8 APR keys can be assigned for each RGDVand/or

and/or

.... DEVICE INTERFACE 3 : 8 APR keys can be assigned for each RGDV

• Enter the RGDVs and APRs programmed for the 3 device interfaces in table 2Z-5 "Programming of device interfaces".



ALL RIGHTS RESERVED

OPTIMUS R/F **INSTALLATION**

Assignment of the external APRT modules:

Device interface	Adapter	Programming W1/W2/W3	Hardware connection	Remark
1	WA or 1 WA	W1-W3 open	WA X23 1 RGDVSL (n) 2 RGDVSL (m) 3 SLPG 1 4 SLPG 2 5 SLPG 3 6 SLPG 4 7 SLPG 5 8 SLPG 6 9 PO26V	2 selectable RGDVs can be programmed. They are independent of the APRs programmed. Maximum number of APRs: 6 x 8 (RGDV) = 48 APRs Recommendation: 2 APRs assigned to each button.
2	2 WA	only W1 closed	2 WA X23 \$\frac{1}{2} \text{ SLPG 1} 2 \text{ SLPG 2} 3 \text{ SLPG 3} 4 \text{ SLPG 4} 5 \text{ SLPG 5} 6 \text{ SLPG 6} 7 \text{ SLPG 7} 8 \text{ SLPG 8} PO26V	RGDV cannot be selected. Maximum number of APRs: 8 x 8 (RGDV) = 64 APRs Recommendation: 2 APRs assigned to each button. Operation as 2 WA is possible also without 1 WA. Refer to drawing Z2-15.1 Cabinet wiring: Decade adapter WA
3	WB	only W2 closed	WB X21 RGDVSL 1 RGDVSL 2 WB X23 SLPG 1 SLPG 2 SLPG 3 SLPG 3 SLPG 4 SLPG 4 SLPG 5 SLPG 5 SLPG 5 SLPG 6 SLPG 7 SLPG 7 SLPG 7 SLPG 7 SLPG 7 SLPG 8 PO26V	2 selectable RGDVs can be programmed. They are independent of the APRs programmed. Maximum number of APRs: 8 x 8 (RGDV) = 64 APRs Recommendation: 2 APRs assigned to each button

INSTALLATION OPTIMUS R/F

8.7. **Option: Tomo Density Control TDC**

For this exposure technique the APRs must have the following programming:

- Dose measurement field: ON (at least 1 field must be set to ON)

- Preferred technique: automatic

- AEC technique: TDC

- No AEC technique: kV-mAs-ms technique (RUQT) or kV-mA-ms (RUIT)

 Exposure data U: = anatomical kV

 Exposure data Q/I: = anatomical mAs/mA product based on the screen/film combination used

- Exposure time t: = anatomical exposure time

The mAs product is used to calculate the start current, indicated under "Exposure data I". In the APR files supplied all the APRs for tomography applications are programmed to TDC. If there is no TDC option installed, the manual technique is selected as the preferred technique automatically.

TDC is not restricted to tomography applications only so it can be preferred for all exposures where exposure time is the determining factor.

The respective mAs product is generally based on a 400-type screen/film combination and must be adapted to the combinations actually used.

E. g.: If a 200-type combination is used, the mAs product must be doubled.

2-82 (c/02.1)OPTIMUS R/F OPTIMUS_RF_2_c021

OPTIMUS R/F INSTALLATION

8.8. Option: VARIOFOCUS

For the VARIOFOCUS option, special APR files have to be loaded. These are designated with a "V" in the file name and contain correspondingly defined APRs.

VARIOFOCUS is programmed as a percentage mix of the small focus with the large focus.

The following steps are possible:

20%, 35%, 50%, 65% and 80% of the small focus.

Default: 50%.

The percentage mix is not displayed directly on the control desk and can only be estimated indirectly via the small/large focus exposure time.

It is only possible to display and change the percentage mix via XRGSCOPE.

Select menu:

PROGRAM/ HUMAN INTERFACE/ APR DATA SET/ CHANGE APR DATA SET

However, VARIOFOCUS can be selected via the control desk and can be stored as default focus for an APR, too.

The default ratio is always 50%.

VARIOFOCUS is possible for tubes with superimposed focal spots only.

The following tubes are suitable for the application: RO 1648, RO 1750, SRO 0951, SRO 2550, SRO 33 100.

8.9. Option: Area Dose Calculator

This option operates in conjunction with a unit and a collimator which are CAN-controlled and supply information about SID, collimation and added filters only.

See section 6: ADJUSTMENTS chapter 1.x: Checking / Correction

8.10. Acceptance test

Execute the acceptance test.
 See section 7: ACCEPTANCE

· Observe all applicable national regulations.

For U.S.A applications check the H.H.S requirements!

After completition of setting-to-work, the system must be tested for H.H.S. compliance according the P.M.S.I. comprehensive compliance testing workbook: Numeric code 4535 800 2035x.

INSTALLATION OPTIMUS R/F

8.11. Interlock facility for APR modification

· Select menu:

PROGRAM/ HUMAN INTERFACE/ APR MODIFIABLE BY USER

It is possible to prevent a customer from being able to store APR modifications as default setting via the control desk.

Default: YES

8.12. Backup of all configuration data

A hardware key is required for the PC.

To save the configuration data use the CONFIGURATION BACKUP disk supplied.

• Save the complete SW programming of the generator on the floppy disk by using the menu: ACCEPT/ BACKUP/ ...

1. ... RGDV RELATED ASSIGNMENTS/ RGDV 1 ... 8/ APR ASSIGNMENTS

Use APRMAN.

After READING (turning bar) the PC offers the default backup name:

Backup filename: APR BAK1 ... 8.tdl

The name can be changed into any other file name.

The path (harddisk) is automatically taken into account.

It is also possible to type:

A:\"filename" <RETURN>

to load the backup files directly to the floppy disk.

2. ... CU complete

Note

Pay attention to the rules of "PC and generator settings to avoid problems during up/downloading of CU complete files" described in chapter 7.1 of this section.

A disk space of 700kByte is required.

It takes about 8min to save the data to the disk.

The default backup name:

Backup File Name: CUBACKUP.TDL

The name can be changed into any other filename.

The path (harddisk) is automatically taken into account.

It is also possible to type:

A:\"filename" <RETURN>

to load the backup files directly to the floppy disk.

ALL RIGHTS RESERVED

OPTIMUS R/F INSTALLATION

9. Labels

· Check the labelling according to the respective generator type.

See drawing 2Z-10 "Labelling".

All lables become visible by swiveling out the label bracket simply by hand and without any tool. The bracket is located at the top left corner of the front side of the cabinet, visibly marked by an "i" (for information) and text "Certified Component Labels Here". If the label bracket is swiveled 90 degrees to the right the following labels appear at its bottom side:

- X-ray control: - type designation

- serial No.

- name and address of manufacturer

- DHHS certification statement (if necessary)

- date of manufacture

- X-ray H.V. generator: - type designation

- serial No.

- name and address of manufacturer

- DHHS certification statement (if necessary)

- date of manufacture

- Technical data label with UL / CSA classification (if necessary)

INSTALLATION OPTIMUS R/F

10. Final installation work


- Mount the side panels of the generator cabinet.
- Take care that all cables inside the wall junction box are routed in **closed** loops without any kinks. Push the generator cabinet against the wall.



Warning!

Block the two front wheels of the cabinet with the locking screws to guarantee that unauthorized persons cannot accidentally touch parts of the generator which might be dangerous.

- If necessary, level the cabinet with the locking screws.
- Mount the front cover of the generator.



A4 01-10-04 We optimus_rt_zz-2_0_010

RGDV8 Tomo mode switch : X11:3 SL_XG_TO / Bucky RGDV : X11:1 Format + :10 Bucky ready / Tomo RGDV : X11:1 Format + :5 Tomo ready switch related X11:3 --> Bucky - Tomo remote switchover RGDVs RGDV7 Tomo time 5: Tomo time 7: Tomo time 8: Tomo time 6 Bucky / Tomo 1WA ∖ Tomo time [s] RGDV6 Time setting for input at WA X21:1 RGDV5 Tomo time 1: Tomo time 2: Tomo time 3: Tomo time 4: RGDV4 WAX12 RGDV3 WAX11 Bucky / Tomo 1WA: Decade Bucky 1/2 RGDV2 Bucky RGDV - switch related Tomo RGDV - switch related Tomo mode switch **Bucky RGDV** Bucky RGDV RGDV1 Name Exposure series/Tomo movement (No break after exposure end): Desk For WBX11:9--10 (ready) und 1--2 (format size correction contact)_ RGDV1[] RGDV2[] RGDV3[] RGDV4[] RGDV5[] RGDV6[] RGDV7[] RGDV8[] - Medium II format kV correction (dose equiv. steps) Bucky / Scopo 1WB / Decade Bucky 1 (WBX11) Small II format kV correction (dose equiv. steps) - Medium II format density correction (6% steps) - Medium II format mAs correction (6% steps) Small II format density correction (6% steps) Small II format mAs correction (6% steps) - Enable handswitch at generator desk - Mounted radiographical controller - Bucky format density correction : - Release circuit adaptation unit - Density correction (6% steps) - Dose measurement sensor - Mounted tomo extension : - Cone density correction: - Dose measurement input - Tube overload protectior Underexposure display Release circuit number Exposure switch type - Disable time override - Syncmaster present: Used for fluoroscopy - Tube power factor - Used for tomo - Release delay Density steps Data Set B: Bucky RGDV: CT add on : Data Set A - mAs steps mA steps : - time steps - kV steps : - Room : - Tube

RGDV programming

A4 01-10-04 We optimus_rf_2z-2_1_a010

Data Set A:									
	Desk :	RGDV1	RGDV2	RGDV3	RGDV4	RGDV4	RGDV6	RGDV7	RGDV8
- Room :		-	-	1	-	٦			
- Tube :		-	-	2	-	٦			
- Release circuit number :		1	2	ဇ	(4)	4			
- Enable handswitch at generator desk :		ou	ou	yes	yes	yes			
- Syncmaster present :		yes	yes	yes	OU				
- Exposure switch type :		double step	donple step	double step	double step	donple steb			
- Bucky format density correction :		0	0	0	0	0			
- Cone density correction :		0	0	0	0	0			
- Dose measurement input :		EZ X21	EZ X41	EZ X31	none	none			
- Dose measurement sensor :		Scopo amplimat	photo sensor/ amplimat input	Bucky amplimat	(Bucky amplimat)	lat) (Bucky amplimat)			
- Exposure series :		yes	yes	ou	ou	ou			
- Release delay :		enable	enable	enable	enable	enable			
- Mounted radiographical controller :		none	none	none	none	none			
- Release circuit adaptation unit :		1WB	1WB	1WB	none	1WB			
- Mounted tomo extension :		none	none	none	none	none			
- Medium II format kV correction (dose equiv. steps) :		0	0 1)	0	0	0			
- Medium II format density correction (6% steps) :		0	0 1)	0	0	0			
- Medium II format mAs correction (6% steps) :		0	0 1)	0	0	0			
- Small II format kV correction (dose equiv. steps) :		0	0 1)	0	0	0			
- Small II format density correction (6% steps) :		0	0 1)	0	0	0			
- Small II format mAs correction (6% steps) :		0	0 1)	0	0	0			
Data Set B:									
- Used for tomo :		no	OU	ou	OU	ou			
- Used for fluoroscopy :		yes	yes	ou	OU	ou			
- CT add on :		no	OU	ou	ou	ou			
- Disable time override :		no	ou	ou	ou	ou			
- Tube power factor :		100 %	100 %	100 %	100 %	100 %			
- kV steps :		Dose equiv. 1)	Dose equiv. 1)	Dose equiv. 1)	Dose equiv. 1)	1) Dose equiv. 1)			
- mAs steps :		25 % 1)	25 % 1)	25 % 1)	25 % 1)	25 % 1)			
- mA steps :		25 % 1)	25 % 1)	25 % 1)	25 % 1)	25 % 1)			
- time steps :		25 % 1)	25 % 1)	25 % 1)	25 % 1)	25 % 1)			
- Density steps :		12 % 1)	12 % 1)	12 % 1)	12 % 1)	12%			
- Density correction (6% steps) :		0	0	0	0	0			
- Underexposure display :		yes	yes	yes	yes	yes			
- Tube overload protection :		on	no no		o	oo			
1WB/D	Buc	ky / Tomo 1WA : I	Decade Bucky 1/2	WAX11	WAX12 B	Bucky / Tomo 1WA \ Tomo time [s]	no time [s] :		
Bucky RGDV: RGDV1 [x] RGDV2 [x] RGDV3 [x]	Tom	Tomo mode switch	-	:		Tomo time 1 :		Tomo time 5:	
	Pac	Bucky RGDV - switch related	related	:		Iomo time z :	mo!	Iomo time 6 :	!
	Buck	Bucky RGDV				Tomo time 3 :	Tom	Tomo time 7 :	
	Buck	Bucky RGDV				Tomo time 4 :	Tom	Tomo time 8 :	
	Tom	o RGDV - switch r	elated		L	ime setting for input at V	VA X21:18		
For WBX11:910 (ready) und 12 (format size correction	Tom	o mode switch : X1	11:3 SL_XG_TO / BL	icky RGDV : X11:1	Format + :10 Bud	Tomo mode switch : X11:3 SL_XG_TO / Bucky RGDV : X11:1 Format + :10 Bucky ready / Tomo RGDV : X11:1 Format + :5 Tomo ready	: X11:1 Format +	:5 Tomo ready	

RGDV programming: example 1

(a/01.0)

A4 01-10-04 We optimus_ri_2z-2_2010

optimus_rf_2z-2_2_a010										
- Bucky Diagnost TH - D76 : Exposure Sco	- Bucky Diagnost TH with Bucky-Controller - D76 : Exposure Scopo / BV-DSI	Name :	Bucky	Тото	Bucky wall stand	Free cassette	D 76 Scopo	D 76 BV - DSI		
Data Set A:		Desk:	RGDV1	RGDV2	RGDV3	RGDV4	RGDV5	RGDV6	RGDV7	RGDV8
- Room :			1	٦	-	-	-	-		
- Tube :			F	٦	F	-	2	2		
- Release circuit number	ier :		do not care	do not care	do not care	do not care	-	2		
- Enable handswitch at generator desk	t generator desk :		yes	səÁ	yes	səƙ	ou	ou		
- Syncmaster present :			yes	səÁ	yes	səƙ	yes	yes		
- Exposure switch type :			donple step	dets elduob	dets elduob	dets elduob	denple step	donple step		
- Bucky format density correction :	correction :		0	0	0	0	0	0		
- Cone density correction	ion :		0	0	0	0	0	0		
- Dose measurement input	nput :		EZ X21	none / [EZ X21]	EZ X31	euou	EZ X22	EZ X41		
- Dose measurement sensor :	sensor :		Bucky amplimat	Bucky amplimat	Bucky amplimat	(Bucky amplimat)	:) Scopo amplimat	photo sensor / amplimat input		
- Exposure series/Tomo movement	io movement :		no	ou	no	ou	yes	yes		
- Release delay :			enable	enable	enable	enable	enable	enable		
- Mounted radiographical controller :	cal controller :		Bucky contr. 1 / DigitalDiagnost	Bucky contr. 1 / DigitalDiagnost	Bucky contr. 1 / DigitalDiagnost	Bucky contr. 1 / DigitalDiagnost		none		
- Release circuit adaptation unit :	ation unit :		none	none	none	none	1WB	1WB		
- Mounted tomo extension	sion :		none	none	none	none	none	none		
- Medium II format kV	- Medium II format kV correction (dose equiv. steps) :		0	0	0	0	0	0 1)		
- Medium II format den	- Medium II format density correction (6% steps) :		0	0	0	0	0	0 1)		
- Medium II format mA	- Medium II format mAs correction (6% steps) :		0	0	0	0	0	0 1)		
- Small II format kV co.	- Small II format kV correction (dose equiv. steps) :		0	0	0	0	0	(10		
- Small II format densit	- Small II format density correction (6% steps) :		0	0	0	0	0	0 1)		
- Small II format mAs correction (6% steps)	correction (6% steps) :		0	0	0	0	0	0 1)		
Data Set B:										
- Used for tomo :			no	yes	no	ou	no	ou		
- Used for fluoroscopy			no	OU	no	ou	yes	yes		
- CT add on :			UO	ou	no	ou	ou	OU		
- Disable time override	::		no	ou	no	ou	ou	ou		
- Tube power factor:			100 %	100 %	100 %	100 %	100 %	100 %		
- kV steps :			Dose equiv. 1)	Dose equiv. ¹⁾	Dose equiv. 1)	Dose equiv. 1)	Dose equiv. 1)	Dose equiv. 1)		
- mAs steps :			25 % 1)	72 % 1)	25 % 1)	25 % 1)	25 % 1)	25 % 1)		
- mA steps :			25 % 1)	25 % 1)	25 % 1)	25 % 1)	25 % 1)	25 % 1)		
- time steps :			25 % 1)	25 % 1)	25 % 1)	25 % 1)	25 % 1)	25 % 1)		
- Density steps :			12 % 1)	15 % 1)	12 % 1)	15 %	12 % 1)	12 %		
- Density correction (6% steps) :	% steps) :		0	0	0	0	0	0		
- Underexposure display :	ay:		yes	yes	yes	yes	yes	yes		
- Tube overload protection	tion :		on	on		on	on	on		
Bucky / Scopo 1WB /	Bucky / Scopo 1WB / Decade Bucky 1 (WBX11):		Bucky / Tomo 1W	Bucky / Tomo 1WA : Decade Bucky 1/2	2 WAX11	WAX12	no 1WA ∖⊺	omo time [s]		
			Tomo mode switch			:	Tomo time 1:		Tomo time 5:	-
	KGDV1		Bucky RGDV - switch related	itch related			lomo time 2 :		lomo time 6 :	
Bucky RGDV:	RGDV5 [x] RGDV6 [x] RGDV7 [Bucky RGDV				Tomo time 3 :	Ton	Tomo time 7 :	
	RGDV8 []	,	Bucky RGDV				Tomo time 4 :	Ton	Tomo time 8 :	
			Tomo RGDV - switch related	ch related			Time setting for input at WA X21:18	t WA X21:18		
For WBX11 : 910 (re	For WBX11:910 (ready) und 12 (format size correction contact)		Tomo mode switch switch related X11.	: X11:3 SL_XG_TO / 3 > Bucky - Tomo	Bucky RGDV : X11	::1 Format + :10 B RGDVs	Tomo mode switch : X11:3 SL_XG_TO / Bucky RGDV : X11:1 Format + :10 Bucky ready / Tomo RGDV : X11:1 Format + :5 Tomo ready switch related X11:3 > Bucky - Tomo remote switchover RGDVs)V : X11:1 Format +	:5 Tomo ready	
Common			Which ichaca		Tellions of the second second	2,2				

RGDV programming: example 2

1) = has to be adjustet on site [] = TDC

A4 01-10-04 We optimus_rf_2z-2_3_a010

- Bucky wall stand / - Free cassette	- Dao : Cassette+DSI Will/Willout Tollio - Bucky wall stand / - Free cassette	e : D96 cassette	D96 Tomo	DSI	DSI Tomo	Free cassette	Wall stand		
Data Set A:	AsaC	RGDV1	RGDV2	RGDV3	RGDV4	RGDV5	RGDV6	RGDV7	RGDV8
			1	4	- T)))	5	5	5
- Hoom :		-	-	1	1	-	l		
- Tube :		1	1	1	1	2	2		
- Release circuit number	r:	1	2	3	4	1	2		
- Enable handswitch at generator desk	generator desk :	ou	OU	ou	OU	yes	səÁ		
- Syncmaster present :		yes	yes	yes	sək	sek	səƙ		
- Exposure switch type		donple step	double step	double step	double step	donple step	donple step		
- Bucky format density correction:	correction :	0	0	0	0	0	0		
- Cone density correction :	: u	0	0	0	0	0	0		
 Dose measurement input 	put :	EZX 21	none / [EZX 21]	EZX 41	none / [EZX 41]	none	EZX 31		
- Dose measurement sensor :	ensor :	Bucky amplimat	Bucky amplimat	photo sensor/ amplimat input	(photo sensor/ amplimat input)	(Bucky amplimat)	Bucky amplimat		
- Exposure series/Tomo movement	movement :	yes	yes	yes	yes	OU	ou		
- Release delay :		enable	enable	enable	enable	enable	enable		
- Mounted radiographical controller	al controller :	none	none	none	none	none	euou		
- Release circuit adaptation unit	tion unit :	1WB	1WB	1WB	1WB	1WA	1WA		
- Mounted tomo extension	on :	1WA	1WA	1WA	1WA	none	euou		
· Medium II format kV c	- Medium II format kV correction (dose equiv. steps) :	0	0	0 1)	0 1)	0	0		
· Medium II format dens	 Medium II format density correction (6% steps): 	0	0	0 1)	0 1)	0	0		
- Medium II format mAs correction (6% steps)	correction (6% steps) :	0	0	0 1)	0 1)	0	0		
Small II format kV corr	- Small II format kV correction (dose equiv. steps) :	0	0	0 1)	0 1)	0	0		
- Small II format density correction (6% steps)	correction (6% steps) :	0	0	0 1)	0 1)	0	0		
- Small II format mAs correction (6% steps)	prection (6% steps) :	0	0	0 1)	0 1)	0	0		
Data Set B:									
- Used for tomo :		ou	yes	no	yes	OU	ou		
- Used for fluoroscopy:		yes	yes	yes	sək	OU	ou		
- CT add on :		ou	ou	no	ou	OU	ou		
- Disable time override :		ou	ou	ou	ou	OU	ou		
- Tube power factor :		100 %	100 %	100 %	400 %	100 %	100 %		
- kV steps :		Dose equiv. 1)	Dose equiv. 1)	Dose equiv. 1)	Dose equiv. 1)	Dose equiv. 1)	Dose equiv. 1)		
- mAs steps :		25 % 1)	25 % 1)	25 % 1)	25 % 1)	25 % 1)	25 % 1)		
- mA steps :		25 % 1)	25 % 1)	25 % 1)	25 % 1)	25 % 1)	25 % 1)		
- time steps :		25 % 1)	25 % 1)	25 % 1)	22 % 1)	25 % 1)	25 % 1)		
- Density steps :		12 % 1)	12 % 1)	12%	12 %	12 %	12 % 1)		
- Density correction (6% steps)	steps):	0	0	0	0	0	0		
- Underexposure display	у:	yes	no	yes	ou	yes	yes		
- Tube overload protection	on :	uo	on		uo	on	uo		
3ucky / Scopo 1WB / I	Bucky / Scopo 1WB / Decade Bucky 1 (WBX11):	Bucky / Tomo 1	Bucky / Tomo 1WA : Decade Bucky 1/2	1/2 WAX11	WAX12	Bucky / Tomo 1W	Bucky / Tomo 1WA \ Tomo time [s]	3]:	
		Tomo mode switch	ત્	enable	enable	Tomo time 1:	9,0	Tomo time 5:	3,0
	RGDV1 [x] RGDV2 [x] RGDV3 [x]	Bucky RGDV - switch related	witch related	RGDV 1	RGDV3	Tomo time 2:	6,0	Tomo time 6:	1,5
Bucky RGDV:	RGDV4 [x] BGDV5 [1 BGDV6 [1 BGDV7 [1	Bucky RGDV		none	none	Tomo time 3:	1,5	Tomo time 7:	1,0
	RGDV8[]	Bucky RGDV		none	none	Tomo time 4:	2'0	Tomo time 8:	1,0
		Tomo RGDV - switch related	witch related	RGDV 2	RGDV 4	Time setting for in	Time setting for input at WA X21:18	8	
For WBX11 : 910 (rea	For WBX11:910 (ready) und 12 (format size correction	Tomo mode swit	Tomo mode switch: X11:3 SL_XG_TO / Bucky RGDV: X11:1 Format + :10 Bucky ready / Tomo RGDV: X11:1 Format + :5 Tomo ready switch related X11:3 Surky - Tomo remote switchover RGDVs	/ Bucky RGDV : X o remote switchow	(11:1 Format + :10	Bucky ready / Tomo	RGDV : X11:1 Fo	rmat + :5 Tomo rea	dy
(2000)			mer frame and		0.00				

A4 01-10-04 We optimus_ri_2z-2_1_a010

Data Set A: Desk : RGDV1 - Room : 1 - Tube : 1 - Enable handswitch at generator desk : yes - Syncmaster present : yes - Exposure switch type : yes - Exposure switch type : 0 - Bucky format density correction : 0 - Cone density correction : 0 - Cone density correction : 0 - Dose measurement input : 0 - Dose measurement sensor : EZ X21 - Dose measurement sensor : none - Exposure series / Tomo movement : no - Release delay : 1WA - Mounted radiographical controller : none - Mounted tomo extension : 0 - Medium II format kV correction (dose equiv. steps) : 0 - Medium II format was correction (dose equiv. steps) : 0 - Small II format density correction (dose equiv. steps) : 0 - Small II format density correction (dose equiv. steps) : 0 - Small II format density correction (6% steps) : 0 - Small II format mAs correction (6% steps) : <th> 1</th> <th>Mail stand RGDV3 1 1 3 3 yes yes yes double step 0 0 0 EZ X31 Bucky amplimat none 1WA none 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0</th> <th>HGDV4 1 (4) yes no double step 0 0 none none</th> <th> HGDV4</th> <th>RGDV6</th> <th>RGDV7</th> <th>RGDV8</th>	1	Mail stand RGDV3 1 1 3 3 yes yes yes double step 0 0 0 EZ X31 Bucky amplimat none 1WA none 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0	HGDV4 1 (4) yes no double step 0 0 none none	HGDV4	RGDV6	RGDV7	RGDV8
in number: witch at generator desk: esent : esent : ch type: ch t	+++++++++++++++++++++++++++++++++++++++	1 1 1 3 3 yes yes yes yes Couble step 0 0 0 0 EZ X31 Bucky amplimat none none 0 0 0 0 0 0 0 0 0 0 0 0	1 1 (4) yes no	 			
in number : witch at generator desk : esent : ch type : ch typ	+++++++++++++++++++++++++++++++++++++++	1 3 3 yes yes yes yes couble step 0 0 0 0 EZ X31 Bucky amplimat none none 1WA none 0 0 0 0 0 0	(4) yes no no double step 0 0 0 0 none none none none none 0 0 0 0 0 0 0 0 0 0 0 0	 			
inumber : sett at generator desk : ch type : ch type : checking correction : ment input : ment input : sy / Tomo movement : adapathical controller : charablical controller : at density correction (6% steps) : at the K correction (6% steps) : and K correction (6% steps) : and K correction (6% steps) : and K correction (6% steps) : checky	+++++++++++++++++++++++++++++++++++++++	9 9 8 9 9 8 9 9 8 9 9 8 9 9 9 9 9 9 9 9	(4) yes no no double step 0 0 0 0 0 none none none none none 0 0 0 0 0 0 0 0 0 0	 			
virth at generator desk : sesent : cesent : ces	 	yes yes yes double step 0 0 0 0 EZ X31 Bucky amplimat no enable none 1WA none 0 0 0 0	yes noble step double step 0 0 0 0 none none none none 0 0 0 0 0 0 0 0 0				
Second S	+++++++++++++++++++++++++++++++++++++++	yes yes yes double step 0 0 0 EZ X31 Bucky amplimat no enable none 1WA none 0 0 0 0	yes no double step 0 none (Bucky amplimat none none none 0 0 0 0 0 0 0	 			
Sesent : Good		yes double step 0 0 EZ X31 Bucky amplimat no enable none 1WA none 0 0	double step 0 0 0 0 none (Bucky amplimat none none 0 0 0 0 0 0 0	 			
About the content of the content o		double step 0 0 0 EZ X31 Bucky amplimat no enable none 1WA none 0 0 0	double step 0 0 0 0 none (Bucky amplimat none none 0 0 0 0 0 0	- 			
lensity correction : correction : ment input : ment sensor : sel / Tomo movement : graphical controller : radaptation unit : extension : at AV correction (6% steps) : that NA correction (6% steps) : at mAs correction (6% steps) : correction (6% steps) : mAs correction (6% steps) : correction (6% steps) : mAs correction (6% steps) : correction (6% steps) : mAs correction (6% steps) :		Ducky amplimat Bucky amplimat no enable none 1WA none 0 0 0	(Bucky amplimate) (Bucky amplimate) no none none none 0 0 0 0	 			
ment input: ment sensor: ment sensor: self-graphical controller: i adaptation unit: extension: at the correction (6% steps): at mAs correction (6% steps): where the controller (6% steps): mat mAs correction (6% steps): mat mAs correction (6% steps): controller: mat mat correction (6% steps): mat mat correction (6% steps): controller: contr		EZ X31 Bucky amplimat no enable none 1WA none 0 0 0	(Bucky amplimat none no no no no no no no none none no				
ment input: ment sensor: graphical controller: extension: adaptation unit: extension: at KV correction (6% steps): kV correction (6% steps): kV correction (6% steps): mat MAs correction (6% steps): mat mAs correction (6% steps): mAs correction (6% steps): mAs correction (6% steps): correction (6% steps): mAs correction (6% steps): correction (6% steps): scopy:	 	EZ X31 Bucky amplimat no enable none 1WA none 0 0 0	(Bucky amplimat no enable none none none none none none none 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0	 			
ment sensor: :: :: graphical controller: extension : at Av correction (6% steps) : the the correction (6% steps) :	 	Bucky amplimat no enable none 1WA none 0 0 0 0	(Bucky amplimat no enable none none none 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0	 			
ss / Tomo movement: :: graphical controller: extension atlablation unit : extension atl kV correction (6% steps) atl density correction (6% steps) kV correction (6% steps) kV correction (6% steps) can mak correction (6% steps) mak correction (6% steps)	yes enable none 1WA none (1WA) 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0	no enable none 1WA none 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0	no enable none none none 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0	100 enable none none none none none none none 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0			
graphical controller: adaptation unit: extension: att V correction (dose equiv. steps): att density correction (6% steps): kV correction (6% steps): kV correction (6% steps): central mAs correction (6% steps): mAs correction (6% steps): central correction (6% steps): central correction (6% steps): correction (6% steps):	enable none 1WA none (1WA) 0 0 0 0 0 0	enable none 1WA none 0 0 0 0 0	enable none none none 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0	enable none 1WA none 0 0 0 0 0 0			
graphical controller: adaptation unit: extension: att KV correction (dose equiv. steps): att density correction (6% steps): kV correction (6% steps): kV correction (6% steps): central mAs correction (6% steps): mAs correction (6% steps): central mAs correction (6% steps): central mAs correction (6% steps): central mAs correction (6% steps):	1WA 1WA none (1WA) 0 0 0 0 0 0 0	1WA 100e 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0	100ne 100ne 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0	1WA 1WA none 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0			
extension : extension : extension : at kV correction (dose equiv. steps) : at density correction (6% steps) : kV correction (6% steps) : kV correction (6% steps) : chensity correction (6% steps) : mAs correction (6% steps) : : : : : : : : : : : : : : : : : : :	1WA none (1WA) 0 0 0 0 0 0 0	1WA none 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0	100ne 100ne 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0	1WA none 0 0 0 0 0 0			
aut kV correction (dose equiv. steps): att density correction (6% steps): att density correction (6% steps): kV correction (6% steps): kV correction (6% steps): mAs correction (6% steps): : :: :: :::::::::::::::::::::::::::	0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0	0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0	0 0 0 0 0 0	0 0 0 0			
aut kV correction (dose equiv. steps) : att density correction (6% steps) : att mAs correction (6% steps) : kV correction (6% steps) : density correction (6% steps) : mAs correction (6% steps) : : : scopy :	00000	00000	00000	0 0 0			
at density correction (6% steps) : at mAs correction (6% steps) : kV correction (dose equiv. steps) : density correction (6% steps) : mAs correction (6% steps) : : scopy :	00000	0 0 0 0	0 0 0 0	0 0 0			
kV correction (6% steps): kV correction (dose equiv. steps): density correction (6% steps): mAs correction (6% steps):: scopy:	0 0 0 0	0 0 0	0 0 0	0 0			
kV correction (dose equiv. steps): density correction (6% steps): mAs correction (6% steps):: scopy:	0 0 0	0 0	0 0	0			
density correction (6% steps) : mAs correction (6% steps) : : scopy :	0 0	0	0	,			
mAs correction (6% steps) : : scopy :	0	•	•	0			
: :: scopy:	,	C	c	C			
: scopy:		•	,				
: scopy:	(==:) ==	9	9	4			
	no (yes)	UO	OLI	ОП			
	OU	no	OU	no			
	ou	no	ou	no			
	no	no	no	no			
er factor :	100 %	100 %	100 %	100 %			
- kV steps :	Dose equiv. 1)	Dose equiv. 1)	Dose equiv. 1)	Dose equiv. 1)			
	25 % 1)	25 % 1)	25 % 1)	25 % 1)			
	25 % 1)	25 % 1)	25 % 1)	25 % 1)			
	25 % 1)	25 % 1)	25 % 1)	25 % 1)			
12	12 % 1)	12 % 1)	12 %	12 %			
- Density correction (6% steps) :	0	0	0	0			
	ou	yes	yes	yes			
	no	on	on	on			
Bucky / Scope 1WB / Decade Bucky 1 (WBX11): Bucky / Tomo 1WA: Decade Bucky 1/2	: Decade Bucky 1/2	WAX11	WAX12 Bud	Bucky / Tomo 1WA \ Tomo time [s]	_		
		ele Se	disable Ton	Tomo time 1:	D,8	Tomo time 5 :	0,1
RGDV1 [] RGDV2 [] Bucky RGDV - switch related	ch related	RGDV 1	none Ton	Tomo time 2 :	3,2 T	Tomo time 6:	0,1
Bucky RGDV: RGDV4[]		RGDV 3	none Ton	Tomo time 3 :	D,8	Tomo time 7 :	0,1
		none	none Ton	Tomo time 4 :	3,2 T	Tomo time 8 :	0,1
-	h related	RGDV 2	none Tim	Time setting for input at WA X21:18	'A X21:18		
For WBX11:910 (ready) und 12 (format size correction Tomo mode switch: X1 contact)	Tomo mode switch : X11:3 SL_XG_TO / Bucky RGDV : X11:1 Format + :10 Bucky ready / Tomo RGDV : X11:1 Format + :5 Tomo ready switch related X11:3 > Bucky - Tomo remote switchover RGDVs	oky RGDV : X11:1	Format + :10 Buck)	ready / Tomo RGDV :	: X11:1 Forma	t + :5 Tomo ready	
1) = has to be adjusted on site							

A4 01-10-04 We optimus_rf_22-2_5,8010

			OILIO	wall stand		Free cassette			
Data Set A:	Desk :	RGDV1	RGDV2	RGDV3	RGDV4	V4 RGDV5	RGDV6	RGDV7	RGDV8
- Room :		1	٦	-	-				
- Tube:		-	-	-	-				
- Release circuit number :		do not care	do not care	do not care	do not care	care			
- Enable handswitch at generator desk :		yes	yes	yes	yes				
- Syncmaster present :		yes	yes	yes	yes				
- Exposure switch type :		double step	donple step	double step	donple step	step			
- Bucky format density correction :		0	0	0	0				
- Cone density correction :		0	0	0	0				
- Dose measurement input :		EZ X21	none / [EZ X21]	EZ X31	none	Ф			
- Dose measurement sensor :		Bucky amplimat	Bucky amplimat	Bucky amplimat	at (Bucky amplimat)	nplimat)			
- Exposure series / Tomo movement :		no	ou	ou	ou				
- Release delay :		enable	enable	enable		əle			
- Mounted radiographical controller :		Bucky contr. 1 / DigitalDiagnost	Bucky contr. 1 / DigitalDiagnost	Bucky contr. 1 / DigitalDiagnost	Bucky contr. 1 / St DigitalDiagnost	ontr. 1 /			
- Release circuit adaptation unit :		none	euou	none	-	9			
- Mounted tomo extension :		none	euou	none	none	Ф			
- Medium II format kV correction (dose equiv. steps) :		0	0	0	0				
- Medium II format density correction (6% steps) :		0	0	0	0				
- Medium II format mAs correction (6% steps) :		0	0	0	0				
- Small II format kV correction (dose equiv. steps) :		0	0	0	0				
- Small II format density correction (6% steps) :		0	0	0	0				
- Small II format mAs correction (6% steps) :		0	0	0	0				
Data Set B:									
- Used for tomo :		no	yes	ou	OU				
- Used for fluoroscopy :		no	ou	ou	ou				
- CT add on :		OU	ou	ou	ou				
- Disable time override :		no	OU	ou	ou				
- Tube power factor :		100 %	% 001	100 %	100 %	%			
- kV steps :		Dose equiv. 1)	Dose equiv. 1)	Dose equiv. 1)	Dose equiv. 1)	luiv. ¹⁾			
- mAs steps :		25 % 1)	72 % 1)	25 % 1)	25 % 1)	(1)			
- mA steps :		25 % 1)	72 % 1)	25 % 1)	25 % 1)	(1)			
- time steps :		25 % 1)	72 % 1)	25 % 1)	25 % 1)	(1)			
- Density steps :		12 % 1)	12 % 1)	15 % 1)	12 %	%			
- Density correction (6% steps) :		0	0	0	0				
- Underexposure display :		yes	yes	yes	yes	•			
- Tube overload protection :		on		on	0				
Bucky / Scope 1WB / Decade Bucky 1 (WBX11):	Bucky	•	: Decade Bucky 1/2	WAX11 IN	WAX12	Bucky / Tomo time			
	Tomo	Tomo mode switch				Tomo time 1:		Tomo time 5:	
RGDV1[] RGDV2[] RGDV3[]	Bucky	Bucky RGDV - switch related	ted	:	!	Tomo time 2 :		Tomo time 6:	
Bucky RGDV: RGDV4[]	Bucky	Bucky RGDV		:	!	Fomo time 3 :		Tomo time 7 :	
	Bucky	Bucky RGDV			!	Tomo time 4 :		Tomo time 8 :	:
	Tomo	Tomo RGDV - switch related	pe	!	-	Time setting for input at WA X21:18	at WA X21:1	8:	
For WBX11:910 (ready) und 12 (format size correction contact)	Tomo i switch	node switch : X11:3 related X11:3> E	SL_XG_TO / Buck Bucky - Tomo remo	y RGDV : X11:1 I te switchover RG	-ormat + :10 Bu DVs	Tomo mode switch : X11:3 SL_XG_TO / Bucky RGDV : X11:1 Format + :10 Bucky ready / Tomo RGDV : X11:1 Format + :5 Tomo ready switch related X11:3> Bucky - Tomo remote switchover RGDVs	DV : X11:1 Forn	nat + :5 Tomo ready	
otic actotalibe od ot sed = (1	-								

A4 01-10-04 We optimus_rf_zz-2_6_010

BuckyDIAGNOST TS	Name:	Bucky	Тото	wall stand	Free cassette	atte			
Data Set A:	Desk:	RGDV1	RGDV2	RGDV3	RGDV4	RGDV5	RGDV6	RGDV7	RGDV8
- Room :		-	-	-	-				
- Tube :		-	-	-	-				
- Release circuit number :		٦	-	2	(1)				
- Enable handswitch at generator desk :		yes	yes	yes	yes				
- Syncmaster present :		yes	yes	yes	ou				
- Exposure switch type :		donple step	double step	donple step	donple step	de			
- Bucky format density correction :		0	0	0	0				
- Cone density correction :		0	0	0	0				
- Dose measurement input :		EZ X21	none / [EZX21]	EZ X31	none				
- Dose measurement sensor :		Bucky amplimat	Bucky amplimat	Bucky amplimat	at Bucky amplimat	imat			
- Exposure series / Tomo movement :		ou	yes	OU	ou				
- Release delay :		enable	enable	enable	enable				
- Mounted radiographical controller :		none	none	none	none				
- Release circuit adaptation unit :		1WA	1WA	1WA	none				
- Mounted tomo extension :		none	1WA	none	none				
- Medium II format kV correction (dose equiv. steps) :		0	0	0	0				
- Medium II format density correction (6% steps) :		0	0	0	0				
- Medium II format mAs correction (6% steps) :		0	0	0	0				
- Small II format kV correction (dose equiv. steps) :		0	0	0	0				
- Small II format density correction (6% steps) :		0	0	0	0				
- Small II format mAs correction (6% steps) :		0	0	0	0				
Data Set B:									
- Used for tomo :		ou	yes	ou	ou				
- Used for fluoroscopy :		ou	ou	OU	OL				
- CT add on :		OU	ou	OU	OLI				
- Disable time override :		ou	ou	OU	ou				
- Tube power factor :		100 %	100 %	100 %	100 %				
- kV steps :		Dose equiv. 1)	Dose equiv. 1)	Dose equiv. ¹ ,	Dose eqiv. ¹⁾	(1)			
- mAs steps :		25 % 1)	25 % 1)	25 % 1)	25 % 1)				
- mA steps :		25 % 1)	25 % 1)	25 % 1)	25 % 1)				
- time steps :		25 % 1)	25 % 1)	25 % 1)	25 % 1)				
- Density steps :		12 % 1)	12 % 1)	12 % 1)	12 %				
- Density correction (6% steps) :		0	0	0	0				
- Underexposure display :		yes	ou	yes	yes				
- Tube overload protection :		uo		on	0				
Bucky / Scopo 1WB / Decade Bucky 1 (WBX11):	Buck	⋖	: Decade Bucky 1/2		WAX12 Bu	Bucky / Tomo time	,		
	Tomo	Tomo mode switch		enable	Tor	Tomo time 1 :	Tor	Tomo time 5:	
RGDV1[] RGDV2[] RGDV3[]	Bucky	Bucky RGDV - switch related	ated	RGDV1	Tor	Tomo time 2 :	Tor	Tomo time 6 :	
Bucky RGDV: RGDV4[]	Bucky	Bucky RGDV		none	Tor	Tomo time 3 :	Tor	Tomo time 7:	
	Bucky	Bucky RGDV		none	Tor	Tomo time 4 :	Tor	Tomo time 8:	
	Tomo	Tomo RGDV - switch related	ted	RGDV2	Tin	Time setting for input at WA X21:18	/A X21:18		
For WBX11:910 (ready) und 12 (format size correction contact)	Tomo switch	mode switch : X11:3 related X11:3> E	3 SL_XG_TO / Buci 3ucky - Tomo remo	ky RGDV : X11:1 F te switchover RGI	ormat + :10 Buck)Vs	Tomo mode switch : X11:3 SL_XG_TO / Bucky RGDV : X11:1 Format + :10 Bucky ready / Tomo RGDV : X11:1 Format + :5 Tomo ready switch related X11:3> Bucky - Tomo remote switchover RGDVs	: X11:1 Format	⊦ :5 Tomo ready	
1) = has to be adjustet on site									
1	OWIECI	l felateu Att.o A	SUCKY - IUIIIU IGIIII	אוניווטעפו דוכענ	JVS				

4 01-10-04 We

optimus_rf_2z-2_7_010

A4 01-10-04 We optimus_rf_2z-2_8_010

optimus_rf_2z-2_8_010										
BuckyDIAGNOST TH Generotor equipped Aux. for MCS (only) =	BuckyDIAGNOST TH any version with Bucky- Controller Generotor equipped with / without decade adapt. unit WA Aux. for MCS (only) = RGDV4 combined with free cassette	Name :	Bucky	Тото	Bucky wall stand	Free cassette	MCS			
Data Set A:		Desk:	RGDV1	RGDV2	RGDV3	RGDV4	RGDV4	RGDV6	RGDV7	RGDV8
- Room :			-	-	-	-	-			
- Tube :			٦	٦	٦	-	-			
- Release circuit number	per:		-	-	-	-	-			
- Enable handswitch at generator desk :	ւt generator desk ։		yes	yes	yes	yes	yes			
- Syncmaster present			yes	yes	yes	yes	yes			
- Exposure switch type :	0		donple step	donple step	donple step	donple step	donple step			
- Bucky format density correction	y correction :		0	0	0	0	0			
- Cone density correction	tion :		0	0	0	0	0			
- Dose measurement input	input :		EZ X21	none / [EZ X21]	EZ X31	none	EZ X22			
- Dose measurement sensor	sensor:		Bucky amplimat	Bucky amplimat	Bucky amplimat	t (Bucky amplimat)	t) Bucky amplimat	t t		
- Exposure series / Tomo movement	omo movement :		ou	ou	ou	OU	OLI			
- Release delay :			enable	enable			enable			
- Mounted radiographical controller :	ical controller :		Bucky contr. 1 / DigitalDiagnost	Bucky contr. 1 / DigitalDiagnost	Bucky contr. 1 / DigitalDiagnost	/ Bucky contr. 1 / DigitalDiagnost	/ Bucky contr. 1 / DigitalDiagnost	1		
- Release circuit adaptation unit :	itation unit :		none	none	none	none	none			
- Mounted tomo extension	sion :		none	euou	none	none	none			
- Medium II format kV	- Medium II format kV correction (dose equiv. steps) :		0	0	0	0	0			
- Medium II format der	- Medium II format density correction (6% steps) :		0	0	0	0	0			
- Medium II format m/	- Medium II format mAs correction (6% steps) :		0	0	0	0	0			
- Small II format kV cc	- Small II format kV correction (dose equiv. steps) :		0	0	0	0	0			
- Small II format densi	- Small II format density correction (6% steps) :		0	0	0	0	0			
- Small II format mAs correction (6% steps)	correction (6% steps) :		0	0	0	0	0			
Data Set B:										
- Used for tomo :			ou	yes	ou	OU	OL			
- Used for fluoroscopy	:/		ou	ou	ou	OU	ou			
- CT add on:			ou	ou	ou	OU	OU			
- Disable time override			ou	ou	ou	OU	OU			
- Tube power factor:			100 %	100 %	100 %	100 %	100 %			
- kV steps :			Dose equiv. 1)	Dose equiv. 1)	Dose equiv. 1)	Dose equiv. ¹⁾	Dose equiv. 1)			
- mAs steps :			25 % 1)	25 % 1)	25 % 1)	25 % 1)	25 % 1)			
- mA steps :			25 % 1)	25 % 1)	25 % 1)	25 % 1)	25 % 1)			
- time steps :			25 % 1)	25 % 1)	25 % 1)	25 % 1)	25 % 1)			
- Density steps :			12 % 1)	12 % 1)	12 % 1)	12 %	12 % 1)			
- Density correction (6% steps)	3% steps) :		0	0	0	0	0			
- Underexposure display	lay :		yes	yes	yes	yes	yes			
- Tube overload protection	ction :		on	no	on	0	on			
Bucky / Scopo 1WB,	Bucky / Scopo 1WB / Decade Bucky 1 (WBX11):	Buck	Bucky / Tomo 1WA : Decade Bucky 1/2		WAX11 W	WAX12 Bucky	Bucky / Tomo time			
		Tomo	Tomo mode switch			Tomo time 1	ime 1 :	Ton	Tomo time 5:	
	RGDV1[] RGDV2[] RGDV3[]	Bucky	Bucky RGDV - switch related	ited		Tomo time 2	ime 2 :	Ton	Tomo time 6 :	
Bucky RGDV:	RGDV4 RGDV5 1 RGDV6 1 RGDV7 1	Bucky	Bucky RGDV			Tomo time 3	ime 3 :	Ton	Tomo time 7 :	
		Buck	Bucky RGDV			Tomo time 4	ime 4 :	Ton	Tomo time 8 :	-
		Tomo	Tomo RGDV - switch related	peq		Time s	Time setting for input at WA X21:18	۸ X21:18		
For WBX11 : 910 (r contact)	For WBX11 : 910 (ready) und 12 (format size correction contact)	Tomo switch	Tomo mode switch : X11:3 SL_XG_TO / Bucky RGDV : X11:1 Form switch related X11:3> Bucky - Tomo remote switchover RGDVs	SL_XG_TO / Buc 3ucky - Tomo remo	ky RGDV : X11:1 Fe ote switchover RGD	ormat + :10 Bucky re Vs	.XG_TO / Bucky RGDV : X11:1 Format + :10 Bucky ready / Tomo RGDV : X11:1 Format + :5 Tomo ready y - Tomo remote switchover RGDVs	X11:1 Format +	.:5 Tomo ready	

A4 01-10-04 We optimus_rf_2z-2_9_010

eptimus_rt_2z-2_9.010 BuckyDIAGNOST TH any version	Name :	o o	L C	Bucky					O C N
deneroror equipped with without decade adapt, unit was Aux. for MCS (only) = any of RGDV5 8		Ducky		wall stand	ו בפ כמסספוופ	<u> </u>			200
Data Set A:	Desk:	RGDV1	RGDV2	RGDV3	RGDV4	RGDV5	RGDV6	RGDV7	RGDV8
- Room :		-	-	-	-				-
- Tube :		-	٠	-	-				-
- Release circuit number :		-	ŀ	-	-				-
- Enable handswitch at generator desk:		yes	yes	yes	yes				yes
- Syncmaster present :		yes	yes	yes	yes				ou
- Exposure switch type :		donple step	donple step	double step	donple step				double step
- Bucky format density correction :		0	0	0	0				0
- Cone density correction :		0	0	0	0				0
- Dose measurement input :		EZ X21	none / [EZ X21]	EZ X31	none				EZ X22
- Dose measurement sensor :		Bucky amplimat	Bucky amplimat	Bucky amplimat	t (Bucky amplimat)	nat)			Bucky amplimat
- Exposure series / Tomo movement :		ou	ou	ou	ou				ou
- Release delay :		enable	enable						enable
- Mounted radiographical controller :		Bucky contr. 1 / DigitalDiagnost	Bucky contr. 1 / DigitalDiagnost	Bucky contr. 1 / DigitalDiagnost	/ Bucky contr. 1 / t DigitalDiagnost	1 / st			none
- Release circuit adaptation unit :		none	none		none				none
- Mounted tomo extension :		none	none	none	none				none
- Medium II format kV correction (dose equiv. steps) :		0	0	0	0				0
- Medium II format density correction (6% steps) :		0	0	0	0				0
- Medium II format mAs correction (6% steps) :		0	0	0	0				0
- Small II format kV correction (dose equiv. steps) :		0	0	0	0				0
- Small II format density correction (6% steps) :		0	0	0	0				0
- Small II format mAs correction (6% steps) :		0	0	0	0				0
Data Set B:									
- Used for tomo :		ou	yes	ou	OU				ou
- Used for fluoroscopy :		ou	ou	ou	OU				ou
- CT add on :		ou	ou	ou	OU				ou
- Disable time override :		ou	ou	ou	OU				ou
- Tube power factor :		100 %	100 %	100 %	100 %				100 %
- kV steps :		Dose equiv. ¹⁾	Dose equiv. 1)	Dose equiv. 1)	Dose equiv. 1)	1)			Dose equiv. 1)
- mAs steps :		25 % 1)	25 % 1)	25 % 1)	25 % 1)				25 % 1)
- mA steps :		25 % 1)	25 % 1)	25 % 1)	25 % 1)				25 % 1)
- time steps :		25 % 1)	25 % 1)	25 % 1)	25 % 1)				25 % 1)
- Density steps :		12 % 1)	12 % 1)	12 % 1)	12 %				12 % 1)
- Density correction (6% steps) :		0	0	0	0				0
- Underexposure display :		yes	yes	yes	yes				yes
- Tube overload protection :		uo	uo	uo	uo				uo
Bucky / Scope 1WB / Decade Bucky 1 (WBX11) :	Bucky	Bucky / Tomo 1WA: Decade Bucky 1/2	-	WAX11 WAX	WAX12 Buc	Bucky / Tomo time			
	Tomo	Tomo mode switch			Tom	Tomo time 1 :		Tomo time 5:	
RGDV1[] RGDV2[] RGDV3[]	Bucky	Bucky RGDV - switch related	ated			Tomo time 2 :		Tomo time 6 :	
Bucky RGDV: RGDV4[]	Bucky	Bucky RGDV				Tomo time 3 :		Tomo time 7 :	* * * *
	Bucky	Bucky RGDV				Tomo time 4 :	*	Tomo time 8:	* *
	Tomo	Tomo RGDV - switch related	ted		Time	Time setting for input at WA X21:1	WA X21:18	8	
For WBX11:910 (ready) und 12 (format size correction contact)	Tomo r	node switch : X11:3 related X11:3 > F	3 SL_XG_TO / Buc Bucky - Tomo rem	ky RGDV : X11:1 Fe ote switchover RGE	ormat + :10 Bucky V/s	Tomo mode switch : X11:3 SL_XG_TO / Bucky RGDV : X11:1 Format + :10 Bucky ready / Tomo RGDV : X11:1 Format + :5 Tomo ready switch related X11:3 > Bucky - Tomo remote switchover RGDVs	V : X11:1 Form	at + :5 Tomo read	<u> </u>
connact		ופומיפים	deny reme	Jo Official Control of the Control o	64				

A4 01-10-04 We optimus_rf_2z-2_10_010

Data Cot A .	Aux. RGDV1 4 must not be used with a Bucky TH stystem via CAN			wall stand					Irauma
	Desk:	RGDV1	RGDV2	RGDV3	RGDV4	RGDV5	RGDV6	RGDV7	RGDV8
- Room :		1	٢	-	-				-
- Tube :		1	-	-	-				2
- Release circuit number :		1	ŀ	ļ	1				٠
- Enable handswitch at generator desk :		yes	yes	yes	yes				yes
- Syncmaster present :		yes	yes	yes	yes				OU
- Exposure switch type :		double step	double step	double step	dets elduop				double step
- Bucky format density correction :		0	0	0	0				0
- Cone density correction :		0	0	0	0				0
- Dose measurement input :		EZ X21	none / [EZ X21]	EZ X31	none				EZ X22
- Dose measurement sensor :	Ш	Bucky amplimat	Bucky amplimat	Bucky amplimat	(Bucky amplimat)				Scopo amplimat
- Exposure series / Tomo movement :		ou	OU	ou	ou				UO
- Release delay :		enable	enable	enable	enable				enable
- Mounted radiographical controller :		Bucky contr. 1 / DigitalDiagnost	Bucky contr. 1 / DigitalDiagnost	Bucky contr. 1 / DigitalDiagnost	Bucky contr. 1 / DigitalDiagnost				none
- Release circuit adaptation unit :		none	none	none	none				none
- Mounted tomo extension :		none	none	none	none				none
- Medium II format kV correction (dose equiv. steps) :		0	0	0	0				0
- Medium II format density correction (6% steps) :		0	0	0	0				0
- Medium II format mAs correction (6% steps) :		0	0	0	0				0
- Small II format kV correction (dose equiv. steps) :		0	0	0	0				0
- Small II format density correction (6% steps) :		0	0	0	0				0
- Small II format mAs correction (6% steps) :		0	0	0	0				0
Data Set B:									
- Used for tomo :		no	yes	ou	ou				OU
- Used for fluoroscopy :		ou	OU	ou	ou				yes
- CT add on :		no	ou	ou	ou				ou
- Disable time override :		no	ou	ou	no				no
- Tube power factor :		100%	100 %	100 %	100 %				100 %
- kV steps :		Dose equiv. ¹⁾	Dose equiv. 1)	Dose equiv. 1)	Dose equiv. ¹⁾				Dose equiv. 1)
- mAs steps :		25 % 1)	25 % 1)	25 % 1)	25 % 1)				25 % 1)
- mA steps :		25 % 1)	25 % 1)	25 % 1)	25 % 1)				25 % 1)
- time steps :		25 % 1)	25 % 1)	25 % 1)	25 % 1)				25 % 1)
- Density steps :		12 % 1)	12 % 1)	12 % 1)	12%				12 % 1)
- Density correction (6% steps) :		0	0	0	0				0
- Underexposure display :		yes	sek	sek	yes				yes
		on		on	0				on
Bucky / Scopo 1WB / Decade Bucky 1 (WBX11):	Bucky /	Bucky / Tomo 1WA: Dec	1/2	WAX11 WAX12	2 Bucky / Tomo time	omo time			
	Tomo mo	Tomo mode switch			Tomo time 1	1:	T	Tomo time 5:	
] RGDV2[] RGDV3[]	Bucky R(Bucky RGDV - switch related	ited		Tomo time 2	.2:	T	Tomo time 6:	
	Bucky RGDV	GDV			Tomo time 3	33:)	Tomo time 7:	* * *
HGDV8	Bucky RGDV	GDV			Tomo time 4	: 4 :)	Tomo time 8:	:
,	Tomo RG	Tomo RGDV - switch related	pel		Time settir	Time setting for input at WA X21:18	VA X21:18		
For WBX11: 910 (ready) und 12 (format size correction	Tomo mo	ode switch : X11:3	SL_XG_TO / Bucky	RGDV : X11:1 Form	Tomo mode switch : X11:3 SL_XG_TO / Bucky RGDV : X11:1 Format + :10 Bucky ready / Tomo RGDV : X11:1 Format + :5 Tomo ready switch salated X11:3 Sucky - Tomo remote switchover RGDVs	/ / Tomo RGDV	': X11:1 Forma	t + :5 Tomo read	У
1) - has to be adjusted on site									

Ch	naract	er display o	on the contro	ol desk	Possible PC	Input
Engli	sh	German	French	Spanish	display (code 850)	at the PC
!		!	!	!	!	
	-					
#		#	£	£	#	
\$		\$	\$	\$	\$	
%		%	%	%	%	
- &		&	&	&	&	
,		,	,	,	,	
(((((
)		*	*) *)	
+		+	+	+	+	
,		,	,	,	, -	
-		-	-	-	-	
/		/	/	/	/	
0		0	0	0	0	
1		1	1	1	1	
2		2	2	2	2	
3		3	3	3	3	
4		4	4	4	4	
5		5	5	5	5	
6		6	6	6	6	
7		7	7	7	7	
8		8	8	8	8	
9		9	9	9	9	
:		:	:	:	:	
;		;	;	;	;	
<		<	<	<	<	
=		=	=	=	=	
>		>	>	>	>	
?		?	?	?	?	
@		§	à	§	@	
Α		Α	Α	Α	А	
В		В	В	В	В	
С		С	С	С	С	
D		D	D	D	D	
Е		Е	Е	Е	Е	
F		F	F	F	F	
G		G	G	G	G	
Н		Н	Н	Н	Н	
ı		I	I	I	1	
J		J	J	J	J	
K		K	К	K	K	
L		L	L	L	L	
M		M	M	M	M	
N		N	N	N	N	
0		0	0	0	0	
P		Р	P	P	P	
Q		Q	Q	Q	Q	
٠, ٠,		~	~	~	l G	

		on the contro		Possible PC display	Input at the PC	
English	German	French	Spanish	(code 850)		
R	R	R	R	R		
S	S	S	S	S		
Т	Т	T	T	Т		
U	U	U	U	U		
V	V	V	V	V		
W	W	W	W	W		
Х	Х	Х	Х	Х		
Υ	Υ	Υ	Υ	Y		
Z	Z	Z	Z	Z		
[Ä		i	[
\	Ö	ç	Ñ	\		
]	Ü	§	i]		
^	٨	^	٨	^		
_		_	_			
,	,	,	,	,		
а	а	a	a	a		
b	b	b	b	b		
С	С	С	С	С		
d	d	d	d	d		
e	е	e	e	e		
f	f	f	f	f		
g h	g h	g h	g h	g h		
i	i	i	i	i		
j	j	j	j	j k		
k	k	k	k	k		
I	I	I	I	1		
m	m	m	m	m		
n	n	n	n	n		
0	0	0	0	0		
р	р	р	р	р		
q	q	q	q	q		
r	r	r	r	r		
s	s	s	s	s		
t	t	t	t	t		
u	u	u	u	u		
V	V	V	V	V		
w	W	w	W	w		
х	х	х	х	х		
у	у	у	у	у		
z	Z	Z	Z	Z		
{	ä	é	•	{	Alt +123	
ł	ö	ù	ñ		Alt +124	
}	ü	é	Ç	}	Alt +125	
~	ß	-	~	~	Alt +126	
				Δ	Alt +127	
4		4	4	á	Alt +160	
<u>-</u>		-		í	Alt +161	
-		+	- -	ó	Alt +162	
		I		Ī		

Chara	cter display o	on the contro	ol desk	Possible PC	Input
English	German	French	Spanish	display (code 850)	at the PC
£	£	£	£	ú	Alt +163
	•	•	•	ñ	Alt +164
§	- §	- §	§	<u>0</u>	Alt +167
"	<u> </u>			į	Alt +168
				8 B	Alt +169
■	=	=	=		Alt +170
=	=	=	=	1/2	Alt +171
_	_	_		/2	Alt +172
					Alt +172
					Alt +174
					Alt +174
0	0	0	0	-5	Alt +176
			1	왕: 육:	Alt +177
± 2	± 2	± 2	± 2		Alt +177
À	À	À	À	L .	Alt +192
Á	Á	Á	Á	Τ	Alt +193
Â	Â	Â	Â	Τ .	Alt +194
Ã	Ã	Ã	Ã	F	Alt +195
Ä	Ä	Ä	Ä	_	Alt +196
Å	Å	Å	Å	+	Alt +197
Æ	Æ	Æ	Æ	ã	Alt +198
Ç	Ç	Ç	Ç	Â	Alt +199
È	È	È	È	L	Alt +200
É	É	É	É	F	Alt +201
Ê	Ê	Ê	Ê	<u>. II</u>	Alt +202
Ë	Ë	Ë	Ë	TI .	Alt +203
I	ı	ı	ı	i.	Alt +204
ĺ	ĺ	ĺ	ĺ	11	Alt +205
Î	Î	Î	Î	JL II	Alt +206
Ϊ	Ϊ	Ϊ	Ϊ	¤	Alt +207
				δ	Alt +208
Ñ	Ñ	Ñ	Ñ	Đ	Alt +209
Ò	Ò	Ò	Ò	Ê	Alt +210
Ó	Ó	Ó	Ó	Ë	Alt +211
Ô	Ô	Ô	Ô	È	Alt +212
Õ	Õ	Õ	Õ	I	Alt +213
Ö	Ö	Ö	Ö	í	Alt +214
				î	Alt +215
Ø	Ø	Ø	Ø	Ϊ	Alt +216
Ù	Ù	Ù	Ù	J	Alt +217
Ú	Ú	Ú	Ú	[Alt +218
Û	Û	Û	Û		Alt +219
Ü	Ü	Ü	Ü		Alt +220
Ý	Ý	Ý	Ý		Alt +221
				ì	Alt +222
В	В	В	В		Alt +223
à	à	à	à	Ó	Alt +224
á	á	á	á	β	Alt +225
â	â	â	â	Ô	Alt +226

Chara	aracter display on the control desk			Possible PC display	Input at the PC
English	German	French	Spanish	(code 850)	
ã	ã	ã	ã	Ò	Alt +227
ä	ä	ä	ä	õ	Alt +228
å	å	å	å	Õ	Alt +229
æ	æ	æ	æ	μ	Alt +230
ç	Ç	Ç	Ç	Þ	Alt +231
è	è	è	è	þ	Alt +232
é	é	é	é	Ú	Alt +233
ê	ê	ê	ê	Û	Alt +234
ë	ë	ë	ë	Ù	Alt +235
ì	ì	ì	ì	ý	Alt +236
í	í	í	í	Ý	Alt +237
î	î	î	î	•	Alt +238
ï	ï	ï	ï	,	Alt +239
					Alt +240
ñ	ñ	ñ	ñ	±	Alt +241
ò	ò	ò	ò	=	Alt +242
ó	ó	ó	ó	3/4	Alt +243
ô	ô	ô	ô	¶	Alt +244
Õ	Õ	Õ	Õ	§	Alt +245
Ö	ö	ö	ö	÷	Alt +246
				•	Alt +247
Ø	Ø	Ø	Ø	٥	Alt +248
ù	ù	ù	ù	••	Alt +249
ú	ú	ú	ú	•	Alt +250
û	û	û	û	I	Alt +251
ü	ü	ü	ü	3	Alt +252
ý	ý	ý	ý	2	Alt +253
				•	Alt +254
					Alt +255

List of characters

		Chamber 1	Chamber 2	Chamber 3	Chamber 4	Chamber 5
	Film:					
_	Screen:					
Set	Chamber:					
ta S	Cassette:					
Data	Sys.corr.:					
	Corr. factor:					
	Film:					
7	Screen:					
Set 2	Chamber:					
ta S	Cassette:					
Data	Sys.corr.:					
	Corr. factor:					
	Film:					
က	Screen:					
Set	Chamber:					
Data S	Cassette:					
Da	Sys.corr.:					
	Corr. factor:					
	Film:					
4	Screen:				•••••	
Set	Chamber:					
Data (Cassette:					
Da	Sys.corr.:					
	Corr. factor:					
	Film:					
2	Screen:					
	Chamber:					
Data Set	Cassette:					
Da	Sys.corr.:					
	Corr. factor:					

A4 99-02-11 Re optimus_r_f_2z-4

Device interface 1

RGDV key 1:	
RGDV key 2:	
RGDV	
APR 1	APR 2
APR 3	APR 4
APR 5	APR 6

RGDV key 1:	
RGDV key 2:	
RGDV	
APR 1	APR 2
APR 3	APR 4
APR 5	APR 6

Device interface 2

RGDV	
APR 1	APR 2
APR 3	APR 4
APR 5	APR 6
APR 7	APR 8

RGDV	
APR 1	APR 2
APR 3	APR 4
APR 5	APR 6
APR 7	APR 8

Device interface 3

RGDV key 1:	
RGDV key 2:	
RGDV	
APR 1	APR 2
APR 3	APR 4
APR 5	APR 6
APR 7	APR 8

RGDV key 1:	
RGDV key 2:	
RGDV	
APR 1	APR 2
APR 3	APR 4
APR 5	APR 6
APR 7	APR 8

RGDV key 1:	
RGDV key 2:	
RGDV	
APR 1	APR 2
APR 3	APR 4
APR 5	APR 6

RGDV key 1:	
RGDV key 2:	
RGDV	
APR 1	APR 2
APR 3	APR 4
APR 5	APR 6

RGDV	
APR 1	APR 2
APR 3	APR 4
APR 5	APR 6
APR 7	APR 8


RGDV	
APR 1	APR 2
APR 3	APR 4
APR 5	APR 6
APR 7	APR 8

RGDV key 1:	
RGDV key 2:	
RGDV	
APR 1	APR 2
APR 3	APR 4
APR 5	APR 6
APR 7	APR 8

RGDV key 1:	
RGDV key 2:	
RGDV	
APR 1	APR 2
APR 3	APR 4
APR 5	APR 6
APR 7	APR 8

A4 97-09-16 Gu optimus_r_f_2z-5

Programming of Device interfaces



OPTIMUS R/F

FAULT FINDING

TEXT

	Contents	3-0.1		
1.	Tools	3-1		
2.	Notes	3-2		
3.	Strategy	3-3		
4.	Service-PC	3-4		
4.1.	Connection	3-4		
4.2.	Operation	3-5		
4.3.	Menu "OPTIMUS" structure	3-6		
4.4.	Saving data on disk and restoring data	3-12		
4.4.1.	PC and generator settings to avoid problems	3-12		
4.4.2.	Preparation of the service PC to guarantee a safe loading process	3-12		
4.4.3.	Preparation of the generator	3-13		
4.4.4.	Saving of data	3-14		
4.4.5.	Restoring of data	3-15		
5.	Initialization phase of the generator	3-16		
5.1.	Start-up sequence	3-16		
5.2.	Program status displayed on the operating panel	3-17		
6.	Switch-ON problems	3-18		
6.1.	Switch-ON not possible	3-18		
6.2.	No start up	3-18		
7.	Error numbers	3-19		
7.1.	Error classification	3-19		
7.2.	Error list	3-20		
7.3.	Elimination of error numbers	3-42		
8.	Power supply	3-44		
9.	Functional description of function unit mA	3-46		
10.	CAN bus	3-49		
11.	Incorrect exposure indicator	3-51		
12.	Mnemonic and routing list	3-54		
13.	Optimus AEC switch-OFF philosophy	3-87		
14.	AEC fault exposure detection strategy	3-89		
15.	Explanations on programming the generator	3-92		
16.	Printed-circuit boards	3-93		
DRAWINGS				
	Central rack, service aid	3Z-1		
	Comparison release decades - CP generators <-> OPTIMUS			
	Companison release decades - OF generators <-> OF HIMOS	UZ-21		

OPTIMUS R/F FAULT FINDING

1. Tools

- Service engineer mechanical tool kit
- mAs meter
- Multimeter
- Digital oscilloscope with 2-beam memory
- PC incl. 3.5" FDD, HW-dongle, serial interface cable, free RAM \geq 590KB
- Service software "XRGSCOPE" Version 2.2 or higher
- Recommended PLCC extraction tool (AMP 822154-1) 2422 487 89772

FAULT FINDING OPTIMUS R/F

2. **Notes**



Warning!

After the generator has been switched OFF, hazardous voltages are still applied to the D.C. intermediate circuits of the converter, the rotor control and the mA control.

These voltages are usually discharged within 2 minutes to values which are no longer dangerous.

For that reason always wait for a minimum of 2 minutes before starting any electrical work always after the generator has been switched OFF.

Note

Permanently interested in quality improvement of PMS products we depend on information from the field. Therefore please send us the current generator logfile information:

Please download the generator errorlog logfile in zipped format as described in chapter 4.4 "Saving data on disk and restoring data".

The filename must express the generator release and generator serial number. E. g. "35012105.tdl" for Rel. "3.5" and serial number "012105".

Send this file containing the serial number of the generator and customer data attached to an E-Mail to:

Carsten Mais Service Innovation Generators PMS DMC Hamburg

E-Mail: Carsten.Mais@philips.com

Appreciate your help and many thanks in advance!

AIL RIGHTS RESERVED

OPTIMUS R/F FAULT FINDING

3. Strategy

There are 3 categories of errors:

1. The generator cannot be switched ON at all or only for a short time.

```
See \Rightarrow 5. "Initialization phase of the generator"
```

- ⇒ 6.1. "Switch ON not possible"
- 2. The generator can be switched ON but no error numbers are displayed on the operating desk. For fault finding use the service PC.

See \implies 4.1. "Connecting the service PC"

⇒ 5. "Initialization phase of the generator"

⇒ 7. "Error numbers"

3. Error messages are displayed on the desk.

For fault finding use the service PC.

See \implies 4.1. "Connecting the service PC"


⇒ 7. "Error numbers"

FAULT FINDING OPTIMUS R/F

4. Service-PC

4.1. Connection

- · Switch the generator ON.
- Provide the service PC with the hardware key and switch it ON. The hardware key provides access to special program settings and to menu "Faultfind". Standard programming is possible without a hardware key.
- Connect the PC to X5 on EZ139 CENTRAL UNIT CU via a serial data cable: (A 5m long data cable can be ordered via 12NC: 4512 130 56931)



AIL RIGHTS RESERVED

4.2. Operation


For installation of generator firmware and newest service tools see "REPLACEMENT" chapter "Exchange of firmware ...".

- Call the program with <**xrgscope**> or with <**xrgscope lcd**> for PCs with LCD screen.
- Enter your password.

The following menu line appears:




Current data files, for instance, for online help, tube types, APR programming are available in the PHILIPS-Intranet. Use path: *http://technet.best.ms.philips.com*/ and pull down menu as shown below.





If you call the installation program with **<xrgscope** ?> the possible starting parameters for the service program are listed.

4.3. Menu "OPTIMUS" structure

OPTIMUS (Release 3.5)









- Accept


```
- RGDV related Assignments
        RGDV 1 ... 8
                        —— APR Assignments
     CU Complete
Restore [ = ]
     RGDV related Assignments
          ----- RGDV 1 ... 8
                      APR Assignments
     CU Complete
Inspect
    - Tube 1 Statistic -
     Tube 2 Statistic

    Show Tube Statistic

     - Tube 3 Statistic –
                             - Reset Tube Statistic [<del>- •</del> ]
     Generator Statistic (not used)
     Type of Tube 1 ... 3
```







[- A hardware key is required

Saving data on disk and restoring data 4.4.

All configuration data and logging tables are stored in battery-buffered CMOS areas of the CU board.

Therefore, these data have to be saved on disk as a backup.

In case data get lost they can easily be restored in the CMOS areas after the error source has been eliminated.

4.4.1. PC and generator settings to avoid problems

Optimus RAD release 3.x CMOS data are up/downloaded in one string without handshake.

Any kind of interruption can cause the loading process to fail.

Problems occur mainly during the download to the PC.

A download file which is not complete cannot be used as a safety backup file.

4.4.2. Preparation of the service PC to guarantee a safe loading process

Start XRGSCOPE always from DOS if possible.

When using any WINDOWS version:

- · Switch OFF all screensavers.
- · Do not run other programs.
- Do not insert any CD in the drive.

Any kind of power management of the PC hardware (BIOS) as well as the windows power management should be switched OFF.

If the PC is connected to mains power some of these might be automatically OFF.

3-12 OPTIMUS R/F (b/02.0)OPTIMUS_RF_3_b020

4.4.3. Preparation of the generator

Preparation of generators without a CAN interface:

· Switch ON the generator.

The loading process can be started once relay ENK1 has been energized.

Preparation of generators which are connected via a CAN interface:

- BuckyDiagnost TH and TH2
- DigitalDiagnost
- Thoravision
- EasyDiagnost with bucky unit
- Switch OFF the generator.
- · Disconnect the following plugs:

System	Connector			
	EZX23 signal bus	EZX42 or EZX42-1 system CAN	EZX43 or EZX43-1 system CAN	
BukyDiagnost TH / TH2	Χ		X	
DigitalDiagnost	Х	Х	Х	
Thoravision	Х	Х	Х	
EasyDiagnost with Bucky unit	X	Х	X	

· Switch ON the generator.

Note

The download procedure must not be started before relay ENK1 has been energized at least 2 minutes after the generator has been switched ON.

Note

In case of mailing CU complete download files or any other xxx.tdl files use a zipped file format. These files are ASCII files which might be destroyed while being mailed.

4.4.4. Saving of data

· Select menu:

ACCEPT/ BACKUP/ CU COMPLETE

• Store the data on floppy disk "Generator configuration data" found in the service documentation:

Default file name: cubackup.tdl

Recommended file name: s/n of the generator, e.g. cu012105.tdl

File size: approx. 500-700kB

Transfer time: approx. 8min.

· Recommendation:

Additionally save the APR programming individually for each RGDV via the menu:

ACCEPT/ BACKUP/ RGDV RELATED ASSIGNMENTS/ RGDV 1 ... 8/ APR ASSIGNMENT

Store them with the SAVE function <F3>-key on floppy disk.

File name: apr bak#.tdl #= RGDV - number

Since release 3 assignments of film/screen combinations to the individual APRs are saved in this procedure!

· Recommendation:

Additionally save the programmings for the film/screen combinations via the menu:

PROGRAM/ DOSE RATE CONTROL/ AMPLIMAT/ CHAMBER 1 ... 5/ DATA SET 1 ... 5 (manual processing)

Store them with the SAVE function <F3>-key on floppy disk.

Recommended file name: drc##.tdl ## = chamber and data set number

4.4.5. Restoring of data

- Select menu: ACCEPT/ RESTORE/ CU COMPLETE
- Restore the data from floppy disk.
 Transfer time approx. 15min ... 50min.
- · Reset the generator.
- · Program date and time.

Most of the programmings and logging tables can also be stored via the SAVE-function <F3>-key of XRGSCOPE. Some programmings can be restored via the LOAD-function <F4>-key.

- For service use, only keep the latest version of the backup.
- Never use a complete backup for a different generator, only if the hardware, firmware and option configuration are identical.
- If backups are used from other generators it is very important to carry out the adjustment of the "factor for duty cycle" as described in section 6 "ADJUSTMENTS", chapter 2.5.
- APR backups can also be loaded into other generators but take care about the release. Use APR manager to translate APR files into the required release format.
 Because specific kV and mA reductions are also transferred, load APR backups only in generators of the same or a lower power class.

Initialization phase of the generator 5.

5.1. Start-up sequence

```
Switch ON of the generator.
Ń
Pulling-up of ENK2.
Selftest of ...
                         ---> all display elements are switched ON for a short moment
| ... control desk C
    central unit EZ139
     kV control EZ130 ---> voltage E is measured in the D.C.. intermediate circuit
     mA control EZ119
     basic interface EZ150
... rotor control EY
      universal I/O EWA/B 102
         Indicating device: The red status LED of the associated printed-circuit board or assembly is illuminated.
After successful selftest the status LEDs blink.
The central unit establishes connection to each functional unit via the CAN bus.
         Indicating device: The red status LED of the associated printed-circuit board or assembly grows dark.
ENK1 is switched ON.
The generator is internally ready.
٧
The external ready circuits are checked ---> unit ready, door contact closed, thermal contact of the tube
                                               closed, tube not overloaded
```

The green READY lamp in the operating desk is illuminated. ===> **The generator is in the READY state.**

Program status displayed on the operating panel 5.2.

PHILIPS OPTIMUS	3		 No tube data loaded yet. No RGDVs programmed yet. No communication between desk and CU. Possible error entries: 00B3, 00B6, 00BA F, 00B0, 00BT, 00BX, 00CJ, 00L1, 00PE, 00XB, 00XL, 03FD
70kV	32.0mAs	Test	Tube data loaded.Selected focus not adapted.
70kV	32.0mAs	Adap	- Status after calling up the adaptation mode.
40kV	00.0mAs	Adap	 Start phase of adaptation mode. After the READY signal appears the adaptation can be started up with the release switch. Possible error entries after adaptation: 00BU, 00BV, 00X6
70kV	320mAs	100ms	 Selected focus is adapted. AEC / TDC technique: For the selected RGDV no measuring unit has been assigned yet.
70kV	0 🔏	def1	 For the selected RGDV no film/screen combination has been programmed yet.
Test APR			 No APR data have yet been loaded onto the selected RGDV.
81kV skull axial Schädel ax.	0 crâne axia cráneo ax		- Ready status. An APR with AEC technique has been selected.

6. Switch-ON problems

6.1. Switch-ON not possible

Z1-2.1 / 2.2 / 2.3 See drawings:

Z2-2

H1 on PCB EN100 is not illuminated

Error sources: - ENF1 was released.

For fault finding look in the error buffer.

- ENF1 is not switched ON.

- Mains voltage, especially phase L3, is not present.

- ENF2 was released.

Check: Low-voltage supply

Filament circuit Tube extension Rotor control

External current consumers

- ENF2 is not switched ON.

- PCB EN100 or its connections are not okay.

H1 on PCB EN100 is illuminated

- The EMERGENCY OFF circuit is open. Error sources:

- The operating desk is not connected.

6.2. No start up

Error sources: - EN100 V1 is defective.

The generator receives a continuous reset via signal: reset sw/.

All red LEDs of the generator are illuminated.

Also see Z1-2.1.

Test: Remove link EZX44:14 --- EZX44:6.

- No boot PROM present: EZ139 D3 (see 5Z-1).
- Flash PROMs EZ139 D4/D5 not correctly loaded.

7. Error numbers

7.1. Error classification

Errors:

- Errors are indicated by 4 digits.
- The first two digits indicate the Functional Unit FU reporting the error.

Example:

00xx = CU-functional unit is concerned 02xx = kV-functional unit is concerned 03xx = mA-functional unit is concerned

- The last two digits indicate the error symptom.

Displayed errors (Errors and Fatal errors):

- These errors are indicated on the display of the operating desk for the customer.
 Not all fatal errors come up on desk, use PC.
- The customer must call the service.
 The customer can inform the service about the respective error number and the service can order the spare parts needed at an early stage of the maintenance procedure.

Not displayed errors (Warnings):

- These errors are not relevant for the customer.
- In case an error of this category occurs frequently within a certain period of time, a displayed error can be generated.

7.2. **Error list**

Sources of error codes indicated in the first two digits decimally (hexadecimally):

Error code dec (hex)	FU (Function Unit)	Description
00xx	CU	central unit EZ139
01xx	FU_DRC	 dose rate control, control physically located on CU EZ139 parts of basic interface FU_CIE EZ150 also involved (Amplimat) FU_DRC also handles fluoro kV control EZ130
02xx	FU_kV	kV control EZ130
03xx	FU_mA_a	- 1st mA control EZ119 - handles 2 filaments
07xx	FU_CIE	central interface extension EZ150 basic interface
08xx	FU_HI_a	human interface 1: C300
10xx (0Axx)	FU_RC_a	1st rotor control high speed EY100
13xx (0Dxx)	FU_ADAP_a	adapter decade cable for 4 aux. units RAD WA/1WA or WA/1WA102
14xx (0Exx)	FU_ADAP_b	adapter decade cable for 4 aux. units RAD 2WA or 2WA102
15xx (0Fxx)	FU_ADAP_c	adapter decade cable for 4 aux. units R/F WB/1WB or WB/1WB102

Error classes: Fatal error, Error, Warning

Error code	Error class	Error text	FU
00B0	ERROR	CPU: Error in application data service interface	CU
00B1	WARNING	CPU: IIM was not expected by gen_order_list	CU
00B2	WARNING	CPU: HI order is not expected - NO Member in display_tab	CU
00B3	ERROR	NVRAM: Data language selector is invalid	CU
00B4	WARNING	CPU: Message invalid in ADopmes	CU
00B5	WARNING	CPU: Inputparameter out of range in ADsynta	CU
00B6	ERROR	NVRAM: FU adap data for DI are invalid	CU
00B7	ERROR	CPU: Message cannot be sent	CU
00B8	ERROR	NVRAM: Tomo mode switch cannot be enabled	CU
00BA	WARNING	NVRAM: Data of RGDV are invalid	CU
00BB	WARNING	NVRAM: Basedata of RGU are invalid	CU
00BC	WARNING	NVRAM: Statedata of RGU are invalid	CU
00BD	WARNING	NVRAM: Data of APR are invalid	CU
00BE	WARNING	NVRAM: Data of active RGU are invalid	CU
00BF	WARNING	NVRAM: Data of RGKeys are invalid	CU
00BG	ERROR	APR: No more lowest level menus available	CU
00BH	ERROR	APR: Display position collision	CU
00BI	ERROR	APR: Menu / APR mismatch in same level	CU
00BJ	ERROR	APR: Menu name not found	CU
00BK	ERROR	APR: APR is assigned to a different RGDV	CU
00BL	ERROR	APR: Menu name already exists	CU
00BM	ERROR	APR: Max. display position reached	CU
00BN	ERROR	APR: APR not found in this menu	CU
00BO	WARNING	NVRAM: Data of menu tree are invalid	CU
00BQ	ERROR	CPU: APR cannot be modified	CU
00BR	ERROR	CPU: APR is not assigned to an RGDV	CU
00BS	ERROR	APR: The RGDV of the APR is not ready for operation	CU
00BT	WARNING	NVRAM: Data of APR characteristics are invalid	CU
00BU	WARNING	Adaptation paused due to missing load	CU
00BV	WARNING	CPU: TTS status message during adaptation	CU

Error code	Error class	Error text	FU
00BW	ERROR	APR: APR not accepted by general calculation	CU
00BX	WARNING	NVRAM: Variofocus allowed invalid	CU
00BY	WARNING	RGDV order without active RGDV	CU
00CB	WARNING	CONF: Received IIM #1#2H unknown	CU
00CC	WARNING	CAN: Frame-repeat-counter overflow (IIM #1#2H)	CU
00CD	WARNING	CAN: FU #1H not addressable	CU
00CE	WARNING	CAN: rx-signal conflict (FU #1H)	CU
00CF	WARNING	CAN: No RTR from FU #1H	CU
00CG	WARNING	CPU: Domain tx response, Mailbox type wrong	CU
00CH	WARNING	CPU: Invalid tbdor-Parameter FU_type	CU
00CI	NOT USED	CAN: No FU acknowledges	CU
00CJ	WARNING	CAN auto configuration successfull (#1H)	CU
00CK	WARNING	CAN auto configuration without success (#1H)	CU
00CL	WARNING	CAN: FU #1H not addressable	CU
00CM	WARNING	CAN: FU #1H sent event and did not answer RTR	CU
00CP	WARNING	CAN: Max FU count exceeded	CU
00CQ	NOT USED	SYSCAN: Radiography system is not responding	CU
00CR	WARNING	SYSCAN: Guarded connection failed	CU
00CX	WARNING	CAN: Last-only-repeat-counter overflow (IIM #1#2H)	CU
00CY	WARNING	CAN: Abort of rx of IIM #1#2H (unexp frame)	CU
00CZ	WARNING	CAN: Unexpected frame received after IIM #1#2H	CU
00DA	WARNING ERROR	No CPU-Access to CAN-chip	CU
00DB	NOT USED	CAN-chip reset not acknowledged	CU
00DC	NOT USED	CAN-chip reset release not acknowledged	CU
00DD	WARNING	CAN-chip DPRAM check failed	CU
00DE	WARNING ERROR	Unexpected CAN-chip int-pointer	CU
00DF	NOT USED	CAN-chip state undefined	CU
00DG	WARNING	CAN-chip error-active after passive #1H	CU
00DH	WARNING	CAN-chip state error-passive #1H	CU
00DI	ERROR	CAN-chip state bus OFF #1H	CU
00DJ	NOT USED	CAN-chip state DPRAM-error	CU
00DK	NOT USED	CAN-chip state DPRAM-error & passive	CU

Error code	Error class	Error text	FU
00DL	NOT USED	Unexpected CAN-chip interrupt	CU
00DM	WARNING	CAN frame error (code #1H)	CU
00E0	ERROR	iRMX exception #2#1H occurred	CU
00G0	WARNING ERROR	Variable in case statement has undefined value	CU
00G1	ERROR	Condition_code <> OK after CALL to send	CU
00G2	WARNING	Condition_code <> OK after CALL to init	CU
0011	NOT USED	CPU index to I/O-table is wrong	CU
0012	WARNING	No interrupt reason on sig-bus	CU
0013	WARNING	No interrupt reason on XS-bus	CU
0014	ERROR	One FU has a WD-error, scantime_TV is not programmed correctly See: XRGSCOPE> OPTIMUS> Program> Dose rate control> CONT: scantime_TV = 20.00ms	CU
00L1	ERROR	GC: Checksum error	CU
00L2	ERROR	GC: Data access error	CU
00L3	ERROR	GC: Limit data error	CU
00L4	WARNING ERROR	GC: Limits inconsistent	CU
00L5	ERROR	GC: Calculation error	CU
00L6	ERROR	GC: Function not implemented	CU
00M0	ERROR	Unable to initialise FU(s) #1H, #2H, #3H, #4H, #5H, #6H	CU
00M1	ERROR	Configuration key is missing or defective	CU
00M2	NOT USED	Unable to initialise the FU mA	CU
00M3	ERROR	No response at all from FU(s) #1H, #2H, #3H, #4H, #5H, #6H	CU
00PA	WARNING	CPU: IIM/MSC number unknown	CU
00PB	WARNING	CPU: Technique mode unknown	CU
00PC	WARNING	CPU: Value limit overflow	CU
00PD	ERROR	PC comm: Unknown TDL proc ID	CU
00PE	WARNING	NVRAM: DRC NV checksum error	CU
00PF	WARNING	CPU: Equal kV-sets from CU come twice	CU
00PG	WARNING	CPU: kV sequence does not increase	CU
00PH	WARNING	CPU: EDL is not possible, min_mA limit	CU
00PI	WARNING	CPU: DCALC Dr_curve has only one element	CU
00PJ	WARNING	CPU: DCALC Dr_curve has strange values	CU

Error code	Error class	Error text	FU
00PK	WARNING	CPU: Equal kV-sets from CU with equal mA	CU
00PL	WARNING	CPU: Dose digits disturbance	CU
00S*	SERVICE	PCcomm: Service access trace	CU
00S?	WARNING ERROR	PCcomm: Unexpected error	CU
00\$0	ERROR	PCcomm: Tube programming error	CU
00SA	ERROR	PCcomm: Not enough space at destination segment	CU
00SB	NOT USED	PCcomm: Base out of range	CU
00SC	ERROR	PCcomm: Value too large	CU
00SD	ERROR	PCcomm: Terminator not found	CU
00SE	ERROR	PCcomm: Error in description	CU
00SF	ERROR	PCcomm: Item type unknown	CU
00SG	ERROR	PCcomm: Internal type unknown	CU
00SH	ERROR	PCcomm: Value negative	CU
00SI	NOT USED	PCcomm: Not enough space at destination buffer	CU
00SJ	ERROR	PCcomm: Syntax wrong	CU
00SK	ERROR	PCcomm: String too long	CU
00SL	WARNING	PCcomm: String truncated	CU
00SM	WARNING	PCcomm: TDL segment overflow	CU
00SN	ERROR	PCcomm: FU Reference Table full	CU
00SO	ERROR	PCcomm: Node ID unknown	CU
00SP	ERROR	PCcomm: FU Code unknown	CU
00SQ	ERROR	PCcomm: Syntax error in node ID	CU
00SR	WARNING	PCcomm: No node ID found	CU
00SS	ERROR	PCcomm: Request not performed	CU
00ST	ERROR	PCcomm: RMX error	CU
00SU	WARNING	PCcomm: Enumeration element not found	CU
00SV	ERROR	PCcomm: Mail corrupted	CU
00SW	ERROR	PCcomm: Procedure ID unknown	CU
00SX	ERROR	PCcomm: FU mA incompatible	CU
00SY	ERROR	PCcomm: FU Off request failed	CU
00SZ	ERROR	PCcomm: Wrong response	CU
00T?	ERROR	TTS: Unexpected error	CU
00TA	ERROR	TTS: Received message unknown	CU

Error code	Error class	Error text	FU
00TB	ERROR	TTS: Tube supervision error from FU kV. Thermal switch of tube housing okay?	CU
00TC	ERROR	TTS: Internal TTS error	CU
00TD	ERROR	TTS: Tube number unknown	CU
00TE	ERROR	TTS: NVRAM checksum error	CU
00TF	ERROR	TTS: NVRAM unavailable	CU
00TG	ERROR	TTS: Tube overheated	CU
00TH	WARNING	TTS: Load data supply inconsistent	CU
00X0	ERROR	CPU wrong timer ID	CU
00X1	ERROR	CPU wrong timer mode	CU
00X2	ERROR	CPU wrong message type	CU
00X3	WARNING	CPU DWORD does not fit into BYTE3	CU
00X4	WARNING ERROR	Timeout of X-ray backup timer	CU
00X5	WARNING	Timeout of X-ray rotation timer	CU
00X6	WARNING	Timeout setting FUs, response missing	CU
00X7	WARNING	CPU curve token is NO_TOKEN	CU
00XA	NOT USED	NVRAM switch table invalid	CU
00XB	WARNING	NVRAM tube data rotation invalid	CU
00XC	WARNING	NVRAM watch dog invalid	CU
00XD	WARNING	NVRAM configuration table invalid	CU
00XE	WARNING	NVRAM test data invalid	CU
00XF	WARNING	NVRAM RoCo data invalid	CU
00XG	NOT USED	CPU received IIM is unknown	CU
00XH	NOT USED	CPU received FU-type is unknown	CU
00XI	ERROR	Init with FU-RoCo not OK	CU
00XJ	WARNING ERROR	Exposure time too short	CU
00XK	WARNING	CPU FUmA refuses set data	CU
00XL	WARNING	NVRAM tube yield table invalid	CU
00XM	WARNING	NVRAM add filter corr table invalid	CU
00XN	WARNING	NVRAM wedge filter corr table invalid	CU
00XO	ERROR	Exposure time too long	CU
00XP	WARNING	Exposure time too long	CU
00XQ	WARNING	NVRAM tube statistic data invalid	CU

Error code	Error class	Error text	FU
00XR	WARNING	NVRAM gsta data invalid	CU
00XS	WARNING	Tube no. in CU and FUkV different	CU
00XT	WARNING	Rotation in CU and FURoCo FUCIE different	CU
00XU	ERROR	Transition endless loop	CU
00XV	WARNING	NVRAM HW test flags invalid	CU
00XW	ERROR	EN_X active in startup	CU
00XX	ERROR	RD_PR_X stays active after prep	CU
02AB	WARNING	Procedure called with wrong parameter	FU_kV
02AC	ERROR	Wrong index for table access	FU_kV
02AD	ERROR	Wrong do case entry	FU_kV
02AE	WARNING	Unknown IIM received	FU_kV
02AF	WARNING	IIM parameter out of range	FU_kV
02CA	WARNING	Error in CASE selector	FU_kV
02CB	WARNING	A CAN message with wrong IIM-no (no recipient defined) received	FU_kV
02CC	WARNING	Multiple receiption of the same CAN frame (transmitter ill)	FU_kV
02CE	WARNING	Unexpected signal value in CAN rx task	FU_kV
02CF	WARNING	CAN bus timeout while domain transmission	FU_kV
02CG	WARNING	Token of CAN response mailbox is not a mailbox token	FU_kV
02CX	WARNING	Multiple rx of the same CAN last/only frame (transmitter ill)	FU_kV
02CY	WARNING	Aborted CAN domain receive (because of timeout or wrong signal)	FU_kV
02CZ	WARNING	Unexpected CAN domain frame received (outside IIM-reception)	FU_kV
02DA	WARNING	No CPU access to the CAN controller	FU_kV
02DB	WARNING	Reset or release of the CAN controller was not acknowledged	FU_kV
02DD	WARNING	Check of the DPRAM of the CAN controller failed	FU_kV
02DE	WARNING	Unexpected interrupt pointer in the CAN controller	FU_kV
02DF	WARNING	CAN controller state undefined	FU_kV
02DG	WARNING	CAN controller state ERROR ACTIVE after ERROR PASSIVE	FU_kV
02DH	WARNING	CAN controller state ERROR PASSIVE	FU_kV
02DI	WARNING	CAN controller state BUS OFF	FU_kV
02DJ	WARNING	CAN controller state DPRAM ERROR	FU_kV
02DK	WARNING	CAN controller state DPRAM ERROR and ERROR PASSIVE	FU_kV
02EA	ERROR	Interrupt 0: Divide by zero	FU_kV
02EB	ERROR	Interrupt 1: Single step	FU_kV

Error code	Error class	Error text	FU
02EC	ERROR	Interrupt 2: NMI	FU_kV
02ED	ERROR	Interrupt 3: Breakpoint	FU_kV
02EE	ERROR	Interrupt 4: Overflow exception	FU_kV
02EF	ERROR	Interrupt 5: Array bounds exception	FU_kV
02EG	ERROR	Interrupt 6: Unused opcode	FU_kV
02EH	ERROR	Interrupt 7: ESC opcode	FU_kV
02EI	ERROR	CAN connection to CU lost	FU_kV
02GA	WARNING	Interpolation not possible	FU_kV
02HA	WARNING	kV nominal value out of range: ± (4% + 1kV); 3 detections within 30ms	FU_kV
02HB	ERROR	kV nominal value out of range: 0kV > U > 170kV	FU_kV
02HC	WARNING	Z nominal value out of range: $\pm1\%\pm0.2$; 3 detections within 30ms; duty cycle range 3% 30%	FU_kV
02HD	ERROR	Z nominal value out of range: 0% > Z > 50%	FU_kV
02HE	WARNING	kV value during standby too large: > 3kV for > 400ms after PREP	FU_kV
02HF	ERROR	kV value during standby too large: > 4kV for > 400ms after PREP	FU_kV
02HG	WARNING	kV actual value out of range: ± (4% + 1kV); 2 detections within 20ms	FU_kV
02HH	ERROR	kV actual value out of range: 20kV > U > 170kV; 3 detections within 30ms	FU_kV
02HI	WARNING	E value during standby out of range: 470V > E > 780V; 3 detections within 30ms	FU_kV
02HJ	ERROR	E value during standby out of range: 450V > E > 800V; 3 detections within 30ms	FU_kV
02HK	WARNING	E value during high tension out of range: 400V > E > 780V; 3 detections within 30ms	FU_kV
02HL	ERROR	E value during high tension out of range: 350V > E > 800V; 3 detections within 30ms	FU_kV
02HM	WARNING	Converter 1 temperature out of range: 0°C > T > 85°C; 3 detections within 30ms	FU_kV
02HN	ERROR	Converter 1 temperature out of range: 0°C > T > 90°C; 3 detections within 30ms	FU_kV
02HO	WARNING	Converter 2 temperature out of range: 0°C > T > 85°C; 3 detections within 30ms	FU_kV
02HP	ERROR	Converter 2 temperature out of range: 0°C > T > 90°C; 3 detections within 30ms	FU_kV

Error code	Error class	Error text	FU
02HQ	WARNING	High tension tank temperature out of range: 0°C > T > 80°C; 3 detections within 30ms	FU_kV
02HR	ERROR	High tension tank temperature out of range: 0°C > T > 85°C; 3 detections within 30ms	FU_kV
02HS	WARNING	Divider test cathode out of range: 45.5kV > U > 50.5kV; 3 detections within 30ms	FU_kV
02HT	ERROR	Divider test cathode out of range: 43kV ≥ U > 53kV; 3 detections within 30ms	FU_kV
02HU	WARNING	Divider test anode out of range: 45.5kV > U > 50.5kV; 3 detections within 30ms	FU_kV
02HV	ERROR	Divider test anode out of range: 43kV ≥ U > 53kV; 3 detections within 30ms	FU_kV
02HW	WARNING	kV asymmetrical: ± 15%; 2 detections within 20ms	FU_kV
02HX	ERROR	kV asymmetrical: ± 15%; 3 detections within 30ms	FU_kV
02MA	ERROR	State request not accepted because of grid mode	FU_kV
02MB	ERROR	State request not accepted because of error state	FU_kV
02MC	WARNING	State requested by CU unknown	FU_kV
02OA	ERROR	RMX error: Timeout	FU_kV
02OB	ERROR	RMX error: Memory	FU_kV
02OC	ERROR	RMX error: Busy	FU_kV
020E	ERROR	RMX error: Limit	FU_kV
02OF	ERROR	RMX error: Context	FU_kV
02OG	ERROR	RMX error: Exist	FU_kV
02OH	ERROR	RMX error: State	FU_kV
0201	ERROR	RMX error: Not configured	FU_kV
02OJ	ERROR	RMX error: Interrupt saturation	FU_kV
02OK	ERROR	RMX error: Interrupt overflow	FU_kV
02OL	ERROR	RMX error: Transmission	FU_kV
02OM	ERROR	RMX error: Divide by zero	FU_kV
02ON	ERROR	RMX error: Overflow	FU_kV
0200	ERROR	RMX error: Type	FU_kV
02OP	ERROR	RMX error: Parameter	FU_kV
02OQ	ERROR	RMX error: Bad call	FU_kV
02OR	ERROR	RMX error: Array bound	FU_kV
02OS	ERROR	RMX error: NDP error	FU_kV

Error code	Error class	Error text	FU
02OT	ERROR	RMX error: Illegal opcode	FU_kV
02OU	ERROR	RMX error: Emulator trap	FU_kV
02OV	ERROR	RMX error: Interrupt table limit	FU_kV
02OW	ERROR	RMX error: CPU xfer data limit	FU_kV
02OX	ERROR	RMX error: Wrap around	FU_kV
02OY	ERROR	RMX error: Check exception	FU_kV
02OZ	ERROR	RMX error: Unknown	FU_kV
02RA	WARNING	Grid mode changeover requested during prep	FU_kV
02RB	WARNING	Tube switch requested during preparation	FU_kV
02RC	WARNING	Requested P out of range	FU_kV
02SA	WARNING	Not enough space at the destination	FU_kV
02SB	WARNING	Base out of range	FU_kV
02SC	WARNING	PC comm.: Value too large	FU_kV
02SD	WARNING	Terminator not found	FU_kV
02SE	WARNING	PC comm.: Error in description	FU_kV
02SF	WARNING	PC comm.: Item type unknown	FU_kV
02SG	WARNING	PC comm.: Internal type unknown	FU_kV
02SH	WARNING	PC comm.: Value negative	FU_kV
02SI	WARNING	PC comm.: No space at dest. buffer	FU_kV
02SJ	WARNING	PC comm.: Syntax wrong	FU_kV
02SK	WARNING	PC comm.: String too long	FU_kV
02SL	WARNING	PC comm.: String truncated	FU_kV
02SO	WARNING	PC comm.: Unknown table ID received	FU_kV
02SP	WARNING	PC comm.: Access level to low	FU_kV
02SQ	WARNING	PC comm.: Unknown action requested	FU_kV
02SR	WARNING	PC comm.: Routing or message corrupt	FU_kV
02SS	WARNING	Source buffer to small for incoming message	FU_kV
02ST	WARNING	CAN buffer to small for outgoing message	FU_kV
02SU	WARNING	PC comm.: Access level is N/A	FU_kV
02UA	ERROR	HW configuration identifier wrong	FU_kV
02UB	WARNING	Set Up request received during preparation	FU_kV
02WA	WARNING	Wrong tube selected	FU_kV
02WB	ERROR	Wrong tube selected	FU_kV

Error code	Error class	Error text	FU
02WC	WARNING	EN X C signal faulty	FU_kV
02WD	ERROR	EN X C signal faulty	FU_kV
02WE	WARNING	Wrong grid mode selected	FU_kV
02WF	ERROR	Wrong grid mode selected	FU_kV
02WG	WARNING	Tube arcing detected	FU_kV
02WH	ERROR	Tube arcing detected	FU_kV
02WI	WARNING	kV over voltage detected	FU_kV
02WJ	ERROR	kV over voltage detected	FU_kV
02WK	WARNING	Measuring not stable	FU_kV
02WL	ERROR	Tube supervision error	FU_kV
02WM	ERROR	Tube supervision error	FU_kV
03AA	WARNING	Internal parameter error	FU_mA_a
03AB	WARNING	Wrong parameter from CU	FU_mA_a
03AC	WARNING	Ie-regulation active on two filaments; only in case of VARIOFOCUS	FU_mA_a
03AI	WARNING	Wrong IIM received	FU_mA_a
03BA	WARNING	Coordinates not monotonous	FU_mA_a
03BB	WARNING	No measurement values for adap. found	FU_mA_a
03CA	WARNING	Error in CASE selector	FU_mA_a
03CB	WARNING	A CAN message with wrong IIM-no (no recipient defined) received	FU_mA_a
03CC	WARNING	Multiple receiption of the same CAN frame (transmitter ill)	FU_mA_a
03CE	WARNING	Unexpected signal value in CAN rx task	FU_mA_a
03CF	WARNING	CAN bus timeout while domain transmission	FU_mA_a
03CG	WARNING	Token of CAN response mailbox is not a mailbox token	FU_mA_a
03CX	WARNING	Multiple rx of the same CAN last/only frame (transmitter ill)	FU_mA_a
03CY	WARNING	Aborted CAN domain receive (because of timeout or wrong signal)	FU_mA_a
03CZ	WARNING	Unexpected CAN domain frame received (outside IIM-reception)	FU_mA_a
03DA	WARNING	No CPU access to the CAN controller	FU_mA_a
03DB	WARNING	Reset or release of the CAN controller was not acknowledged	FU_mA_a
03DD	WARNING	Check of the DPRAM of the CAN controller failed	FU_mA_a
03DE	WARNING	Unexpected interrupt pointer in the CAN controller	FU_mA_a
03DF	WARNING	CAN controller state undefined	FU_mA_a
03DG	WARNING	CAN controller state ERROR ACTIVE after ERROR PASSIVE	FU_mA_a
03DH	WARNING	CAN controller state ERROR PASSIVE	FU_mA_a

Error code	Error class	Error text	FU
03DI	WARNING	CAN controller state BUS OFF	FU_mA_a
03DJ	WARNING	CAN controller state DPRAM ERROR	FU_mA_a
03DK	WARNING	CAN controller state DPRAM ERROR and ERROR PASSIVE	FU_mA_a
03EA	ERROR	CPU interrupt 0	FU_mA_a
03EB	ERROR	CPU interrupt 1	FU_mA_a
03ED	ERROR	CPU interrupt 3	FU_mA_a
03EE	ERROR	CPU interrupt 4	FU_mA_a
03EF	ERROR	CPU interrupt 5	FU_mA_a
03EG	ERROR	CPU interrupt 6	FU_mA_a
03EH	ERROR	CPU interrupt 7	FU_mA_a
03EI	ERROR	CAN is unable to send an error to CU	FU_mA_a
03FA	WARNING	NVRAM: Invalid checksum	FU_mA_a
03FB	WARNING	NVRAM: Standby filament not found	FU_mA_a
03FC	ERROR	No NVRAM plugged in	FU_mA_a
03FD	WARNING	NVRAM empty; battery?	FU_mA_a
03GA	ERROR	Limit error	FU_mA_a
03GB	WARNING	Real math. error: Real underflow	FU_mA_a
03GC	WARNING	Real math. error: Real overflow	FU_mA_a
03GD	WARNING	Real math. error: Dword overflow	FU_mA_a
03GE	WARNING	Real math. error: Integer overflow	FU_mA_a
03GF	WARNING	Real math. error: Word overflow	FU_mA_a
03GG	WARNING	Singular matrix	FU_mA_a
03НА	ERROR	Unknown hardware	FU_mA_a
03НВ	WARNING ERROR	Intermediate circuit voltage < 200V	FU_mA_a
03HF	WARNING	Undefined analog input channel	FU_mA_a
03HG	WARNING	If-actual out of tolerance	FU_mA_a
03HH	ERROR	If setpoint to large	FU_mA_a
03HI	ERROR	If-actual out of tolerance	FU_mA_a
03HJ	ERROR	If-actual out of tolerance	FU_mA_a
03HK	WARNING	If-nominal out of tolerance	FU_mA_a
03HL	ERROR	If-nominal out of tolerance	FU_mA_a
ознм	ERROR	If-nominal out of tolerance	FU_mA_a
03HN	ERROR	No retrigger received from CU	FU_mA_a

Error code	Error class	Error text	FU
03IA	WARNING	Adaptation cannot be completed	FU_mA_a
03IC	WARNING	No Ie-adaptation measurement values	FU_mA_a
03ID	WARNING	Ie-adaptation values not evaluable	FU_mA_a
03KA	WARNING	CondiX-Ray mode without mAs parameter	FU_mA_a
03МА	WARNING	Undefined status	FU_mA_a
03MB	WARNING	Status change not allowed	FU_mA_a
03МС	WARNING	FU init data not expected	FU_mA_a
03OA	ERROR	RMX exception: E\$TIME	FU_mA_a
03OB	ERROR	RMX exception: E\$MEM	FU_mA_a
03OC	ERROR	RMX exception: E\$BUSY	FU_mA_a
03OD	ERROR	RMX exception: E\$LIMIT	FU_mA_a
03OE	ERROR	RMX exception: E\$CONTEXT	FU_mA_a
03OF	ERROR	RMX exception: E\$EXIST	FU_mA_a
03OG	ERROR	RMX exception: E\$STATE	FU_mA_a
03OH	ERROR	RMX exception: E\$NOT\$CONFIGURED	FU_mA_a
03OI	ERROR	RMX exception: E\$INTERRUPT\$SATURATION	FU_mA_a
03OJ	ERROR	RMX exception: E\$INTERRUPT\$OVERFLOW	FU_mA_a
03OK	ERROR	RMX exception: E\$TRANSMISSION	FU_mA_a
03OL	ERROR	RMX exception: E\$ZERO\$DIVIDE	FU_mA_a
03OM	ERROR	RMX exception: E\$OVERFLOW	FU_mA_a
03ON	ERROR	RMX exception: E\$TYPE	FU_mA_a
0300	ERROR	RMX exception: E\$PARAM	FU_mA_a
03OP	ERROR	RMX exception: E\$BAD\$CALL	FU_mA_a
03OQ	ERROR	RMX exception: E\$ARRAY\$BOUND	FU_mA_a
03OR	ERROR	RMX exception: E\$NDP\$ERROR	FU_mA_a
03OS	ERROR	RMX exception: E\$ILLEGAL\$OPCODE	FU_mA_a
03OT	ERROR	RMX exception: E\$EMULATOR\$TRAP	FU_mA_a
03OU	ERROR	RMX exception: E\$INTERRUPT\$TABLE\$LIMIT	FU_mA_a
03OV	ERROR	RMX exception: E\$CPUXFER\$DATA\$LIMIT	FU_mA_a
03OW	ERROR	RMX exception: E\$SEG\$WRAP\$AROUND	FU_mA_a
03OX	ERROR	RMX exception: E\$CHECK\$EXCEPTION	FU_mA_a
03OY	ERROR	Unknown RMX exception	FU_mA_a
03PA	ERROR	le zero measured	FU_mA_a

Error code	Error class	Error text	FU
03PB	WARNING	le out of tolerance: $\pm 10\%$ (le > 5mA, exp. time ≤ 44 ms) or $\pm 3\%$ (le > 5mA, exp. time > 44ms)	FU_mA_a
03PC	ERROR	le out of tolerance: $\pm 30\%$ (le > 5mA, exp. time > 44ms)	FU_mA_a
03PD	WARNING	Setpoint for le-regulation incorrect	FU_mA_a
03PE	ERROR	Emergency OFF! Grid not closed!	FU_mA_a
03PF	ERROR	No kV discharged due to missing le	FU_mA_a
03SC	WARNING	PC comm.: Value too large	FU_mA_a
03SE	WARNING	PC comm.: Error in description	FU_mA_a
03SF	WARNING	PC comm.: Item type unknown	FU_mA_a
03SG	WARNING	PC comm.: Internal type unknown	FU_mA_a
03SH	WARNING	PC comm.: Value negative	FU_mA_a
03SI	WARNING	PC comm.: No space at dest. buffer	FU_mA_a
03SJ	WARNING	PC comm.: Syntax wrong	FU_mA_a
03SK	WARNING	PC comm.: String too long	FU_mA_a
03SL	WARNING	PC comm.: String truncated	FU_mA_a
03SO	WARNING	PC comm.: Unknown table ID received	FU_mA_a
03SP	WARNING	PC comm.: Access level to low	FU_mA_a
03SQ	WARNING	PC comm.: Unknown action requested	FU_mA_a
03SR	WARNING	PC comm.: Routing or message corrupt	FU_mA_a
03SU	WARNING	PC comm.: Access level is N/A	FU_mA_a
07CA	ERROR	CAN: Case-selector error	FU_CIE
07CB	WARNING	CAN: Invalid CAN ID %u	FU_CIE
07CC	ERROR	CAN: Frame rep. overflow IIM%u	FU_CIE
07CD	ERROR	CAN: No RTR from CU	FU_CIE
07CE	ERROR	CAN: rx signal conflict IIM%u	FU_CIE
07CF	ERROR	CAN: tx timeout	FU_CIE
07CI	WARNING	CAN: IMPOSSIBLE ERROR	FU_CIE
07CP	WARNING	CAN: CPU: PXerr %d %s(%d)	FU_CIE
07CR	WARNING	CAN: CPU: Message request fail	FU_CIE
07CS	WARNING	CAN: CPU: Message send error	FU_CIE
07CY	ERROR	CAN: rx abort IIM%u	FU_CIE
07CZ	WARNING	CAN: Unexpected frame (IIM%u)	FU_CIE
07DA	ERROR	CAN: Chip access error	FU_CIE

Error code	Error class	Error text	FU
07DB	ERROR	CAN: Chip reset error	FU_CIE
07DC	ERROR	CAN: Chip reset release error	FU_CIE
07DE	WARNING	CAN: Illegal interrupt pointer	FU_CIE
07DF	ERROR	CAN: Chip state undefined	FU_CIE
07DG	WARNING	CAN: Chip err act. after pass.	FU_CIE
07DH	WARNING	CAN: Chip state error passive	FU_CIE
07DI	WARNING	CAN: Chip state bus OFF	FU_CIE
07DJ	ERROR	CAN: Chip DPRAM error	FU_CIE
07DK	WARNING	CAN: Chip DPRAM error & passive	FU_CIE
07DL	WARNING	CAN: Unexpected interrupt	FU_CIE
07LA	WARNING	CV received IIM unknown	FU_CIE
07LB	WARNING	RC stator number out of range	FU_CIE
07LC	WARNING	RC stator not available	FU_CIE
07LD	ERROR	RC stator 1 readback failed	FU_CIE
07LE	ERROR	RC stator 2 readback failed	FU_CIE
07LF	ERROR	RC stator 3 readback failed	FU_CIE
07LG	WARNING	RC speed value out of range	FU_CIE
07LH	ERROR	RC speed set timeout	FU_CIE
07LI	WARNING	RC maximal stator load exceeded	FU_CIE
07LJ	ERROR	RC maximal rotation time exceeded	FU_CIE
07LK	WARNING	AM amplimat chamber number out of range	FU_CIE
07LL	WARNING	AM amplimat field number out of range	FU_CIE
07LM	WARNING	AM amplimat delay value out of range	FU_CIE
08CA	ERROR	CAN: Case-selector error	FU_HI_a
08CB	WARNING	CAN: Invalid CAN ID %u	FU_HI_a
08CC	ERROR	CAN: Frame rep. overflow IIM%u	FU_HI_a
08CD	ERROR	CAN: No RTR from CU	FU_HI_a
08CE	ERROR	CAN: rx signal conflict IIM%u	FU_HI_a
08CF	ERROR	CAN: tx timeout	FU_HI_a
08CI	WARNING	CAN: IMPOSSIBLE ERROR	FU_HI_a
08CP	WARNING	CAN: CPU: PXerr %d %s(%d)	FU_HI_a
08CR	WARNING	CAN: CPU: message request fail	FU_HI_a
08CS	WARNING	CAN: CPU: message send error	FU_HI_a

Error code	Error class	Error text	FU
08CY	ERROR	CAN: rx abort IIM%u	FU_HI_a
08CZ	WARNING	CAN: Unexpected frame (IIM%u)	FU_HI_a
08DA	ERROR	CAN: Chip access error	FU_HI_a
08DB	ERROR	CAN: Chip reset error	FU_HI_a
08DC	ERROR	CAN: Chip reset release error	FU_HI_a
08DE	WARNING	CAN: Illegal interrupt pointer	FU_HI_a
08DF	ERROR	CAN: Chip state undefined	FU_HI_a
08DG	WARNING	CAN: Chip err act. after pass.	FU_HI_a
08DH	WARNING	CAN: Chip state error passive	FU_HI_a
08DI	WARNING	CAN: Chip state bus OFF	FU_HI_a
08DJ	ERROR	CAN: Chip DPRAM error	FU_HI_a
08DK	WARNING	CAN: Chip DPRAM error & passive	FU_HI_a
08DL	WARNING	CAN: Unexpected interrupt	FU_HI_a
08HA	ERROR	No message receive displaytask	FU_HI_a
08HB	ERROR	No message release displaytask	FU_HI_a
08HC	ERROR	APR not found	FU_HI_a
08HD	ERROR	Offset in menu structure out of range	FU_HI_a
08HF	ERROR	No message request for test task	FU_HI_a
08HG	ERROR	No message send for test task	FU_HI_a
08НН	ERROR	APR buffer full	FU_HI_a
08HI	ERROR	No message send for test task	FU_HI_a
08HJ	ERROR	No send message to CU from ODD	FU_HI_a
08HK	ERROR	Data error in CAN message	FU_HI_a
08HL	ERROR	No message send for service task	FU_HI_a
08IE	ERROR	Wrong setup IIM	FU_HI_a
08SA	ERROR	No request domtxtask when scanning	FU_HI_a
08SB	ERROR	No request domtxtask when testing	FU_HI_a
08SC	ERROR	No send message to task2_sc	FU_HI_a
10CA	ERROR	CAN: Case-selector error	FU_RC_a
10CB	WARNING	CAN: Invalid CAN ID %u	FU_RC_a
10CC	ERROR	CAN: Frame rep. overflow IIM%u	FU_RC_a
10CD	ERROR	CAN: No RTR from CU	FU_RC_a
10CE	ERROR	CAN: rx signal conflict IIM%u	FU_RC_a

Error code	Error class	Error text	FU
10CF	ERROR	CAN: tx timeout	FU_RC_a
10Cl	WARNING	CAN: IMPOSSIBLE ERROR	FU_RC_a
10CP	WARNING	CAN: CPU: PXerr %d %s(%d)	FU_RC_a
10CR	WARNING	CAN: CPU: Message request fail	FU_RC_a
10CS	WARNING	CAN: CPU: Message send error	FU_RC_a
10CY	ERROR	CAN: rx abort IIM%u	FU_RC_a
10CZ	WARNING	CAN: Unexpected frame (IIM%u)	FU_RC_a
10DA	ERROR	CAN: Chip access error	FU_RC_a
10DB	ERROR	CAN: Chip reset error	FU_RC_a
10DC	ERROR	CAN: Chip reset release error	FU_RC_a
10DE	WARNING	CAN: Illegal interrupt pointer	FU_RC_a
10DF	ERROR	CAN: Chip state undefined	FU_RC_a
10DG	WARNING	CAN: Chip err act. after pass.	FU_RC_a
10DH	WARNING	CAN: Chip state error passive	FU_RC_a
10DI	WARNING	CAN: Chip state bus OFF	FU_RC_a
10DJ	ERROR	CAN: Chip DPRAM error	FU_RC_a
10DK	WARNING	CAN: Chip DPRAM error & passive	FU_RC_a
10DL	WARNING	CAN: Unexpected interrupt	FU_RC_a
10FB	ERROR	Short circuit detected	FU_RC_a
10FT	WARNING	Overcurrent detected	FU_RC_a
10IF	WARNING	Initialization failed	FU_RC_a
10LA	WARNING	Acceleration count limit exceeded	FU_RC_a
10LH	ERROR	Phase current %u mA (>%u)	FU_RC_a
10LL	ERROR	Phase current %u mA (<%u)	FU_RC_a
10LO	WARNING ERROR	Intermediate voltage %u V (>%u)	FU_RC_a
10LT	ERROR	Temperature limit exceeded	FU_RC_a
10LU	WARNING ERROR	Intermediate voltage %u V (<%u)	FU_RC_a
10LZ	ERROR	Temperature sensor failure	FU_RC_a
100E	WARNING	CPU: PXROS error %d	FU_RC_a
100F	WARNING	CPU: PXROS error %d %s(%d)	FU_RC_a
10RC	ERROR	Rotation check failed	FU_RC_a
10RI	ERROR	Invalid rotation request : %u	FU_RC_a

Error code	Error class	Error text	FU
10RT	ERROR	Rotation request timeout	FU_RC_a
10TD	ERROR	Invalid data for tube %u	FU_RC_a
10TE	ERROR	Stator %u hardware error	FU_RC_a
10TF	ERROR	Stator %u switching failed	FU_RC_a
10TI	ERROR	Invalid stator request : %u	FU_RC_a
10TR	ERROR	Stator change with rotating anode	FU_RC_a
10UI	WARNING	Unknown message from CU: IIM %u	FU_RC_a
10UM	WARNING	Unexpected message from CU: IIM %u	FU_RC_a
10WT	WARNING	CPU: Watchdog timeout	FU_RC_a
10XX	WARNING	IMPOSSIBLE ERROR	FU_RC_a
13CA	ERROR	CAN: Case-selector error	FU_AD_a
13CB	WARNING	CAN: Invalid CAN ID %u	FU_AD_a
13CC	ERROR	CAN: Frame rep. overflow IIM%u	FU_AD_a
13CD	ERROR	CAN: No RTR from CU	FU_AD_a
13CE	ERROR	CAN: rx signal conflict IIM%u	FU_AD_a
13CF	ERROR	CAN: tx timeout	FU_AD_a
13Cl	WARNING	CAN: IMPOSSIBLE ERROR	FU_AD_a
13CP	WARNING	CAN: CPU: PXerr %d %s(%d)	FU_AD_a
13CR	WARNING	CAN: CPU: Message request fail	FU_AD_a
13CS	WARNING	CAN: CPU: Message send error	FU_AD_a
13CY	ERROR	CAN: rx abort IIM%u	FU_AD_a
13CZ	WARNING	CAN: Unexpected frame (IIM%u)	FU_AD_a
13DA	ERROR	CAN: Chip access error	FU_AD_a
13DB	ERROR	CAN: Chip reset error	FU_AD_a
13DC	ERROR	CAN: Chip reset release error	FU_AD_a
13DE	WARNING	CAN: Illegal interrupt pointer	FU_AD_a
13DF	ERROR	CAN: Chip state undefined	FU_AD_a
13DG	WARNING	CAN: Chip err act. after pass.	FU_AD_a
13DH	WARNING	CAN: Chip state error passive	FU_AD_a
13DI	WARNING	CAN: Chip state bus OFF	FU_AD_a
13DJ	ERROR	CAN: Chip DPRAM error	FU_AD_a
13DK	WARNING	CAN: Chip DPRAM error & passive	FU_AD_a
13DL	WARNING	CAN: Unexpected interrupt	FU_AD_a

Error code	Error class	Error text	FU
13LA	WARNING	CV received IIM unknown	FU_AD_a
13LB	WARNING	IO wrong bidirectional lines output value	FU_AD_a
13LC	WARNING	TR TOMO value for K5 - K12 out of range	FU_AD_a
13LD	WARNING	TR RGDV value out of range	FU_AD_a
13LE	ERROR	TR RGDV readback failed	FU_AD_a
13LF	WARNING	TR wrong sync contact value	FU_AD_a
13LG	WARNING	TR wrong handswitch enable value	FU_AD_a
13LH	ERROR	PR S1/S2 switch active during startup	FU_AD_a
14CA	ERROR	CAN: Case-selector error	FU_AD_b
14CB	WARNING	CAN: Invalid CAN ID %u	FU_AD_b
14CC	ERROR	CAN: Frame rep. overflow IIM%u	FU_AD_b
14CD	ERROR	CAN: No RTR from CU	FU_AD_b
14CE	ERROR	CAN: rx signal conflict IIM%u	FU_AD_b
14CF	ERROR	CAN: tx timeout	FU_AD_b
14Cl	WARNING	CAN: IMPOSSIBLE ERROR	FU_AD_b
14CP	WARNING	CAN: CPU: PXerr %d %s(%d)	FU_AD_b
14CR	WARNING	CAN: CPU: Message request fail	FU_AD_b
14CS	WARNING	CAN: CPU: Message send error	FU_AD_b
14CY	ERROR	CAN: rx abort IIM%u	FU_AD_b
14CZ	WARNING	CAN: Unexpected frame (IIM%u)	FU_AD_b
14DA	ERROR	CAN: Chip access error	FU_AD_b
14DB	ERROR	CAN: Chip reset error	FU_AD_b
14DC	ERROR	CAN: Chip reset release error	FU_AD_b
14DE	WARNING	CAN: Illegal interrupt pointer	FU_AD_b
14DF	ERROR	CAN: Chip state undefined	FU_AD_b
14DG	WARNING	CAN: Chip err act. after pass.	FU_AD_b
14DH	WARNING	CAN: Chip state error passive	FU_AD_b
14DI	WARNING	CAN: Chip state bus OFF	FU_AD_b
14DJ	ERROR	CAN: Chip DPRAM error	FU_AD_b
14DK	WARNING	CAN: Chip DPRAM error & passive	FU_AD_b
14DL	WARNING	CAN: Unexpected interrupt	FU_AD_b
14LA	WARNING	CV received IIM unknown	FU_AD_b
14LB	WARNING	IO wrong bidirectional lines output value	FU_AD_b

Error code	Error class	Error text	FU
14LC	WARNING	TR TOMO value for K5 - K12 out of range	FU_AD_b
14LD	WARNING	TR RGDV value out of range	FU_AD_b
14LE	ERROR	TR RGDV readback failed	FU_AD_b
14LF	WARNING	TR wrong sync contact value	FU_AD_b
14LG	WARNING	TR wrong handswitch enable value	FU_AD_b
14LH	ERROR	PR S1/S2 switch active during startup	FU_AD_b
15CA	ERROR	CAN: Case-selector error	FU_AD_c
15CB	WARNING	CAN: Invalid CAN ID %u	FU_AD_c
15CC	ERROR	CAN: Frame rep. overflow IIM%u	FU_AD_c
15CD	ERROR	CAN: No RTR from CU	FU_AD_c
15CE	ERROR	CAN: rx signal conflict IIM%u	FU_AD_c
15CF	ERROR	CAN: tx timeout	FU_AD_c
15Cl	WARNING	CAN: IMPOSSIBLE ERROR	FU_AD_c
15CP	WARNING	CAN: CPU: PXerr %d %s(%d)	FU_AD_c
15CR	WARNING	CAN: CPU: Message request fail	FU_AD_c
15CS	WARNING	CAN: CPU: Message send error	FU_AD_c
15CY	ERROR	CAN: rx abort IIM%u	FU_AD_c
15CZ	WARNING	CAN: Unexpected frame (IIM%u)	FU_AD_c
15DA	ERROR	CAN: Chip access error	FU_AD_c
15DB	ERROR	CAN: Chip reset error	FU_AD_c
15DC	ERROR	CAN: Chip reset release error	FU_AD_c
15DE	WARNING	CAN: Illegal interrupt pointer	FU_AD_c
15DF	ERROR	CAN: Chip state undefined	FU_AD_c
15DG	WARNING	CAN: Chip err act. after pass.	FU_AD_c
15DH	WARNING	CAN: Chip state error passive	FU_AD_c
15DI	WARNING	CAN: Chip state bus OFF	FU_AD_c
15DJ	ERROR	CAN: Chip DPRAM error	FU_AD_c
15DK	WARNING	CAN: Chip DPRAM error & passive	FU_AD_c
15DL	WARNING	CAN: Unexpected interrupt	FU_AD_c
15LA	WARNING	CV received IIM unknown	FU_AD_c
15LB	WARNING	IO wrong bidirectional lines output value	FU_AD_c
15LC	WARNING	TR TOMO value for K5 - K12 out of range	FU_AD_c
15LD	WARNING	TR RGDV value out of range	FU_AD_c

Error code	Error class	Error text	FU
15LE	ERROR	TR RGDV readback failed	FU_AD_c
15LF	WARNING	TR wrong sync contact value	FU_AD_c
15LG	WARNING	TR wrong handswitch enable value	FU_AD_c
15LH	ERROR	PR S1/S2 switch active during startup	FU_AD_c
16CA	ERROR	CAN: Case-selector error	FU_AD_d
16CB	WARNING	CAN: Invalid CAN ID %u	FU_AD_d
16CC	ERROR	CAN: Frame rep. overflow IIM%u	FU_AD_d
16CD	ERROR	CAN: No RTR from CU	FU_AD_d
16CE	ERROR	CAN: rx signal conflict IIM%u	FU_AD_d
16CF	ERROR	CAN: tx timeout	FU_AD_d
16CI	WARNING	CAN: IMPOSSIBLE ERROR	FU_AD_d
16CP	WARNING	CAN: CPU: PXerr %d %s(%d)	FU_AD_d
16CR	WARNING	CAN: CPU: Message request fail	FU_AD_d
16CS	WARNING	CAN: CPU: Message send error	FU_AD_d
16CY	ERROR	CAN: rx abort IIM%u	FU_AD_d
16CZ	WARNING	CAN: Unexpected frame (IIM%u)	FU_AD_d
16DA	ERROR	CAN: Chip access error	FU_AD_d
16DB	ERROR	CAN: Chip reset error	FU_AD_d
16DC	ERROR	CAN: Chip reset release error	FU_AD_d
16DE	WARNING	CAN: Illegal interrupt pointer	FU_AD_d
16DF	ERROR	CAN: Chip state undefined	FU_AD_d
16DG	WARNING	CAN: Chip err actin after passive	FU_AD_d
16DH	WARNING	CAN: Chip state error passive	FU_AD_d
16DI	WARNING	CAN: Chip state bus OFF	FU_AD_d
16DJ	ERROR	CAN: Chip DPRAM error	FU_AD_d
16DK	WARNING	CAN: Chip DPRAM error & passive	FU_AD_d
16DL	WARNING	CAN: Unexpected interrupt	FU_AD_d
16LA	WARNING	CV received IIM unknown	FU_AD_d
16LB	WARNING	IO wrong bidirectional lines output value	FU_AD_d
16LC	WARNING	TR TOMO value for K5 - K12 out of range	FU_AD_d
16LD	WARNING	TR RGDV value out of range	FU_AD_d
16LE	ERROR	TR RGDV readback failed	FU_AD_d
16LF	WARNING	TR wrong sync contact value	FU_AD_d

Error code	Error class	Error text	FU
16LG	WARNING	TR wrong handswitch enable value	FU_AD_d
16LH	ERROR	PR S1/S2 switch active during startup	FU_AD_d

7.3. Elimination of error numbers

00PL:

The message 00PL (error of the AEC signal) may be a "warning" or an "error". It depends on the disturbance of the AEC signal.

The AEC signal can be measured at pin EZ150 X4 (signal) to EZ150 X3 (see also Z1 " Basic interface ").


When using measuring chambers there are three possibilities to get the error "00PL":

- 1. The shielding of the measuring chamber has a connection to system ground at the measuring chamber or interconnection.
- 2. In the cable to the measuring chamber is a missing ground connection. (This mistake is not possible with the ACL chamber type No. 9890 000 016xx).
- 3. The measuring chamber is defective.

Localization and elimination of the error source:

Re 1.)

- Remove the connector of the measuring chamber at the generator side.
- Measure connection:
 - shielding (Sub-D connector, 15 pins) to system ground or
 - pin M (3+ connector, 14 pins) to system ground ===> The connection must not be present!
- · Measure connection:
 - shielding (Sub-D connector, 15 pins) to chamber shielding
 - pin M (3+ connector, 14 pins) to chamber shielding
 ===> The connection must be present!



Re 2.)

The connector of the measuring chamber at the generator side has been removed.


 Measure the connection between pin 8 and pin 13 (Sub-D connector).

The connection must be present!

If the connection is not present, insert a link between pin 8 and 9 at the chamber connector at the chamber end as shown in the figure.
 In this case the system is most probably operated with an old hybrid measuring chamber 9803 509 xxxxx instead of an ACL measuring chamber 9890 000 016xx. In hybrid measuring chambers the connection between pin 8 <-> 9 is missing.

In case a 3+ connector is used, the connection pin ${\bf N} <-> {\bf J}$ is most probably missing because this connection is not present in hybrid measuring chambers.

To increase interference protection establish the above mentioned connection at the chamber connector of the chamber cable pin **8** <-> **9** in addition to the connection in the adapter for the AMPLIMAT cable (see Z1 "Basic interface").



Re 3.):

Use a test chamber and compare the function.

8. Power supply



Warning!

During fault finding within the power supply unit be very careful the unit is still connected to

Switch ON not possible:

- ENF1 not switched ON (visual check)

- ENF1 released

 check for damage before reactivating ENF1/2 (visual check, any smell?)

- ENF2 not switched ON (visual check)

- ENF2 released by

low-voltage supply filament circuit tube extension

• check for damage before reactivating ENF1/2 (visual check, any smell?)

external components supply

- ON circuit EN100 defective

Phase supervision:

1. Without mains adaptation transformer

Mains contactors ENK2 and ENK1 cannot be activated. - Phase L1 is missing:

- Phase L2 is missing: The generator can be switched ON but does not go into the READY state.

The filament-circuit supply is missing.

There is an error message from function unit kV.

- Phase L3 is missing: ON circuit without supply voltage.

Fault tracing:

Check leads and fuses up to the mains supply.

- 2. With mains adaptation transformer
 - In case at least one phase at the primary end is missing, the generator cannot be switched ON. If there is a problem concerning the leads at the secondary end, refer to section 1 "INTRODUCTION AND TECHNICAL DATA".

AIL RIGHTS RESERVED

After switch-ON or attempted switch-ON:

The generator cannot be brought into the READY state (e.g. no desk display).

- · Check the low-voltage supply.
- · Check for released ENF1:

Ground fault / short-circuit of one / several phase(s).

Check ENK2 and, if necessary, the contacts of ENK1.

Check the leads and the mains adaptation transformer.

Check visually whether the contacts of ENK2 or ENK1 have dropped out.

Check for missed voltage of intermediate-circuit:

The damping resistors are unsoldered which was caused by overcurrent during switch ON.

Cause: Short-circuit in the converter, defective charging capacitors, mains-filter capacitors or rectifiers.

Unsoldering happens about 45s after switch ON.

The damping resistors are unsoldered because the converter was active and ENK1 was not switched ON although activated by the software.

Probably termination of exposure.

This procedure can only happen once because the generator cannot go into STANDBY mode when intermediate-circuit voltage E is missing.

In case intermediate-circuit voltage E is present, ENK1 is activated by the software of the kV-control and remains activated for the complete time the unit is in operation.

In case of high impedance or the tolerance of the symmetry resistors of the intermediate-circuit capacitor battery is too large, capacitors may be destroyed by overvoltage. In case ENK1 has already been activated, ENF1 probably releases.

ENF3 is released by the rotor control units.

The release of ENF2 switches the generator OFF because the supply voltage for the ON circuit and, consequently, the supply voltage of contactors ENK2 and ENK1 is interrupted.

A converter test kit OPTIMUS is available to determine possible problems with the converter, the HT transformer or the tube.

Order no.: 4512 104 9168x

9. Functional description of function unit mA

Tube data must be loaded as a data set from floppy disk via PC and central unit CU into function unit mA.

The procedures described below cannot be carried out before the complete data set for the tube housing assembly is present in central unit CU.

Before the tube adaptation can be started, tube conditioning must be implemented as described in section 2, chapter 8.3.1.

Before adaptation is started, the mA offset value of the mA measuring circuit has to be determined.

This offset value consists of two components:

- 1. A current of 4mA is impressed upon the mA measuring circuit which is used for continuous calibration (during STANDBY about once per minute).
- 2. Additionally the kV measuring circuit delivers an offset current depending on the kV.

To measure this total value an exposure is released with 40kV and 500mA filament current. The emission current measured is the correction value for all standard exposures (4mA, measuring circuit current depending on the kV).

As opposed to the standby filament current value of the predecessor versions of generators, the standby filament current value of the Optimus generator is not fixed.

It is determined for each focus individually. A 40kV exposure is released with the focus to be measured while all other foci are switched OFF.

The filament current changes until an emission current of 100µA is obtained.

The associated filament current value is the individual standby filament current (1% to be substracted so that the fluoroscopic current of any of the other foci is not affected).

The following adaptation program takes place fully automatically.

Based on 120 single exposures for each focus a data field is created in the CMOS of function unit mA. The adjustments for all other exposures are interpolated from this data field during operation.

During the adaptation procedure all limit values such as maximum filament current, maximum kV, maximum tube load, maximum output, current of the generator etc. are taken into account.

3-46 (b/02.0)OPTIMUS R/F © 1997 Philips Medizin Systeme OPTIMUS RF 3 b020

Boost adaptation

Boost time determination: Positive boosting

With the predecessor versions of generators, a calculated boost current was added to the exposure filament current for a fixed time of 400ms.

With Optimus generators the boost current is always fixed but with a variable time.

The amount of the boost current is the sum of the maximum filament current (of the respective filament) plus 2000mA.

To determine the time values an exposure must be started at a kV stage from which on the filament current does not have to be increased anymore to obtain the max. kV dependent emission current.

As soon as the 100% kV value is reached, the filament current jumps from the STANDBY value to the maximum filament current plus 2000mA. The emission current is measured every 2ms until the maximum tube current or the maximum possible tube current is reached.

In case this procedure takes too long (warming up of the tube), the measurement is continued with a second exposure after a sufficient period of time has passed.

The measurement starts again at the value obtained last.

Boost time determination: Negative boosting

An innovation of the Optimus generator is the determination of the negative boosting (blanking of the filament current).

The measurement is started at the same kV stage as for the positive boost time but with maximum filament current. As soon as the 100% kV value is reached, the maximum filament current of the filament jumps down to 500mA. Every 2ms the emission current is measured until a value of 100μ A is obtained.

The values for the blanking times are required for techniques such as, for instance, cine.

A filament current value of 500mA must not be exceeded for otherwise the output to supply a gridswitch box (which might be present) is too low.

The following procedure takes place after the generator has been switched ON:

Function unit mA initializes itself and afterwards establishes connection with central unit CU via CAN.

For 3s every focus is boosted with the respective specified maximum filament current. Then blanking of the filament current (500mA) takes place for a variable period of time (derived from negative boost adaptation) to bring the filament current to the STANDBY value (large focus first followed by a smaller one).

The change of the filament current value upon a change of the focus which was the usual routine for the predecessor versions of the Optimus generator does no longer take place. All STANDBY values remain constant.

During operation the following procedure takes place after the release of PREPARATION:

The filament current rises from the individual STANDBY filament current to the boost current.

The switch-ON time of the boost current depends on the difference between STANDBY and the exposure (single boost) or intermediate filament (double boost) current.

Double boost

 The intermediate filament current is a calculated value. It is calculated in such a way that the filament current and thus the filament temperature is brought to exposure level when the boost current is switched ON for another 50ms by the exposure command.

- During exposure the filament current regulates as required.
- At the end of exposure the filament current is reduced to the minimum value of 500mA (negative boosting) for a calculated time to bring it from the exposure to the STANDBY value.
- In case PREPARATION is released, negative boosting takes place until heating can go on with the STANDBY filament current.

CAN bus 10.

All the intelligent assemblies and PCBs communicate via the CAN bus. There they are connected in parallel to the two lines CAN L (low) and CAN H (high).

The data are serially transmitted in the form of so-called frames.

Levels in quiescent status against chassis:

CAN L: 2.5V CAN_H: 2.5V

Levels during data transmission against chassis:

CAN L: 0.50 ... 2.25V Both levels are opposite. CAN H: 2.75 ... 4.50V J The difference must be > 1.5V!

Test points generator CAN		Test points system CAN	
CAN_L:	EZX71	S_CAN_L:	EZX42:2
CAN_H	EZX72	S_CAN_H:	EZX42:7
Chassis:	EZX5	Chassis:	EZX42:3

Z1-5.1, Z2-5.1/5.2 Reference:

Symptoms of errors:

- The generator is inoperable.
- The red LEDs of one or more of the assemblies or PCBs flash.
- Parameter settings at the control desk are accepted and displayed with a considerable delay.
- In the error memory are several entries which code begins with 00C (apart from 00CJ) or the error description contains a reference to signal conflicts.


Error localization:

- 1. Entries in the error memory clearly indicate that the assembly and PCB are not communicating properly or not at all.
- 2. Control measurement of CAN levels with an oscilloscope during data transmission and in the quiescent status: Data transmission is triggered by pressing any desk button.

If the levels are outside the tolerance or are not symmetrical, the CAN driver of an assembly or PCB is faulty. Because all the users are connected to the bus in parallel, the troublemaker can only be found by disconnecting one user after another.

Disconnection must only take place with the generator switched OFF.

AIL RIGHTS RESERVED



11. Incorrect exposure indicator

General causes:

An incorrect exposure is indicated on the control desk if an exposure cannot be terminated according to the parameters set. Frequent causes of underexposure are the following:

- 1. The operator lets go of the release switch prematurely.
- 2. The measuring chamber is incorrectly programmed, not connected or faulty.

Check the following: - RGDV programming

- Programming of Amplimat sensitivity

- Programming of EZ150 basic interface (gain, 15V/40V-supply)

- Programming of screen/film combination (data set 1 ... 5)

3. The APR selected does not match the technique used or the screen/film combination.

Check the following: - APR programming

The standard APRs supplied have parameters which generally match a 400-type screen/film combination. If the standard APRs are used, the exposure parameters have to be changed according to the speed of the screen/film combination actually used.

This also applies if an automatic technique is programmed as the preferred technique. In automatic techniques the mAs and ms-parameters are used for fault exposure detection.

Fault exposure detection AEC / TDC:

To protect patients there are 3 monitoring systems for automatic techniques:

1. Maximum mAs product: Can be set by XRGSCOPE

2. Maximum exposure time or backup time: Can be set by XRGSCOPE

3. Fault exposure detection: The fault exposure detection can be switched ON or OFF via

XRGSCOPE. Irrespective of this fault exposure detection does

not perform if levels fall below certain limits.

AEC / AECF limits:

- Maximum mAs product: 580mAs (default)

- Maximum exposure time: 4s (cannot be changed)

 Backup time AEC: Exposure time based on 9.5 times the mAs of the respective manual

technique (kV-mAs). 4s after overriding

- Backup time AECF: 9.5 times the exposure time of the respective manual technique (kV-mAs)

- Fault exposure detection: ≤ 4% dose at 10% backup time

Fault exposure detection is ignored under the following circumstances:

- Backup time: \leq 100ms (\leq 10ms at 10%)

 Switch OFF voltage (dose): \leq 610mV (\leq 24.4mV at 4%)

If there is a fault an exposure is aborted after about 10% of backup time. If the fault exposure detection fails to respond in the event of a fault, shutdown takes place after reaching backup time, max. exposure time or max. mAs product.

TDC limits:

 Maximum mAs product: 580mAs (default)

- Exposure time: 0.3 ... 6s

- Fault exposure detection: ≤ 10 ... 4% dose for 9.5 times the sample time

9.5 x sample time

----- x 40% nominal dose dose minimum =

exposure time (corr.)

- Backup time: Exposure time

- Sample time: 25 ... 60ms = 1% exposure time (corr.); min. 25ms

12 ... 100 - Sample steps:

Fault exposure detection is ignored under the following circumstances:

- Exposure time: < 1s

In the event of a fault the exposure is aborted after approx. 11 times sample time. If the fault exposure detection fails to respond in the event of a fault, shutdown takes place after reaching the backup time or the max. mAs product.

The switch OFF voltage should be at least 1.2V to guarantee good TDC regulation. Program the higher gain factor on EZ150 BASIC INTERFACE (≥ 4512 108 05964), if necessary.

AIL RIGHTS RESERVÉD

Programming possibilities:

· Select menu:

PROGRAM/ APPLICATION LIMITS/ X-MODE LIMITS

X-Ray Mode: AEC ... TDC Max. Current Time Product Limit: 580mAs

· Select menu:

PROGRAM/ DOSE RATE CONTROL/ FAULT EXPOSURE DETECTION/ AEC ... TDC

ON - OFF

Aids to fault finding:

· Select menu:

FAULTFIND/ LOGGING TABLE/ X-RAY LOG/ DOSE RATE CONTROL LOGGING/ ...

... READ ACTUAL STATUS Technique and parameters of the last exposure

... AEC/ AEC CALCULATION Data of the selected APR with AEC or AECF

... AEC/ **AEC TRACE** Control values of the last AEC exposure

... TDC/ TDC CALCULATION Data of the selected APR with TDC

... TDC/ **TDC TRACE** Control values of the last TDC exposure

Adjustment possibilities:

• Select menu:

ADJUST/ DOSE RATE CONTROL/ TDC AMPLIMAT

P gain factor (def. 50):

i gain factor (def. 8): Do not change any value here without order from DMC Hamburg!

d gain factor (def. 5):

min. sample time (def. 40) [ms]: 25 ... 65

12. Mnemonic and routing list

Explanation:

MNEMONIC
explanation
signal chain (-direct connection = connection via relay contact) all possible units mentioned
signal value / range / signal source
measuring point (in () at PCB front panel)
trigger point [preferred]
remarks
part of supply

AC 0V XG

AC mains supply 0 V X-ray generator

ENX1101/2-EZX13:2-EZ102X1:DBZ4-EZ119X1:DBZ24 Optimus RAD ENX3201-EZX13:2-EZ102X1:DBZ4-EZ119X1:DBZ24 Optimus R/F

neutral N of mains supply for EZ102 + EZ119

AC 230V L1

AC mains supply 230VAC phase 1

ENF2:L1-ENF2:T1-ENK2:2-ENK2:1-EZX13:1-EZX102X1:DBZ2

AC mains supply for low voltage power supply EZ102

AC 230V L2

mains supply 230V AC phase 2

ENF2:L2-ENF2:T2-ENK2:4-ENK2:3-EZX13:3-EZ119X1:DBZ26

AC mains supply for function unit mA control EZ119

AV_HT_AN

high tension actual value anode side $0V \dots +3.75V = 0 \dots 75kV \quad 1V = 20kV$ measuring point EZ130 (X4) [CRTL_X_C/ at EZX74]

AV_HT_CA

high tension actual value cathode side $0V \dots +3.75V = 0 \dots 75kV \quad 1V = 20kV$ measuring point EZ130 (X5) cathode value also positive! [CRTL X C/ at EZX74]

AV HT

high tension actual value

0 ... +7.5V = 0 ... 150kV 1V = 20kV

measuring point EZ130 (X3)

[CRTL X C/ at EZX74]

3-54 (b/02.0) OPTIMUS R/F

CAN H

generator CAN high active

-EZ119X2:C3-EZ130X2:C3-EZ139X2:C3-EZ150X2:C3-EZX44:10-EZX45:10-EZX46:10

-C300X1:10-EZX51:3-EZX151:3-EZX52:7-EZX72

-EWAX51:10-EWAX52:10-EWA100X2:C3-EWBX51:10-EWBX52:10-EWB100X2:C3

+2.5VDC standby, +3.2VDC during communication

F7X72

for communication of generator function units only part of: XRG bus

CAN L

generator CAN low active

-EZ119X2:A3-EZ130X2:A3-EZ139X2:A3-EZ150X2:A3-EZX44:2-EZX45:2-EZX46:2

-C300X1:2-EZX51:2-EZX151:2-EZX52:2-EZX71-EWAX51:2-EWAX52:2-EWA100X2:A3-EWBX51:2-EWBX52:2-EWB100X2:A3

+2.5VDC standby, +1.5VDC during communication

EZX71

for communication of generator function units only part of: XRG bus

CM EX SW 1

common for exposure switch of release decade 1

EWA100X1:C5-EWAX1:10 EWB100X1:C5-EWBX1:10

+26V non-active exposure request

CM EX SW 2

common for exposure switch of release decade 2

EWA100X1:C7-EWAX2:10 EWB100X1:C7-EWBX2:10

+26V non-active exposure request

CM_EX_SW_3

common for exposure switch of release decade 3

EWA100X1:C9-EWAX3:10 EWB100X1:C9-EWBX3:10

+26V non-active exposure request

CM EX SW 4

common for exposure switch of release decade 4

EWA100X1:C11-EWAX4:10 EWB100X1:C11-EWBX4:10

+26V non-active exposure request

CM SW

common for radiation indication

EZ150X1:C29-EZX1:6-EWGX1:6-EWGX2:6-EWGX3:6

partner of SW UN EX, potential free contact

CM TH

common for thermal sensor of tube housing

NTC temperature measurement in tube housing (not yet available)

EZ130X1:C12-EZX3:7-EWGX7:7-EWGX8:7-EWGX9:7 EZ130X1:C12-EZX3:4-EWGX7:4-EWGX8:4-EWGX9:4 backpanel 4512 108 05983

backpanel 4512 108 05983

backpanels 4512 108 05984 + 4512 108 09361/2

partner of TH OL

CM TH SW

common for tube housing temperature switch

EZ130X1:C11-EZX3:4-EWGX7:4-EWGX8:4-EWGX9:4

EZ130X1:C11-EZX3:7-EWGX7:7-EWGX8:7-EWGX9:7

backpanels 4512 108 05984 + 4512 108 09361/2

partner of TH_OL_SW/

COM EX CD

common for exposure end signal and other warning signals

EWB102X1:A12-EWBX22:6

partner of EX_CD + SW_XG_RD_1 + SW_PR_FL_1 + SW_WN_FL_1 + SW_UN_EX_1

CTRL_X/

control X-ray request command, system level or with decade adaptation units WA/WB

EZ139X1:A4-EZX23:4-EZX45:5-EWAX51:5-EWAX52:5-EWA100X2:C25-EWBX51:5-EWBX52:5

-EWB100X2:C25

0V active, +15V inactive

EZX85

part of: signal bus

CTRL_X_C/

control X-ray request command, internal generator signal

EZ119X2:C6-EZ130X2:C6-EZ139X2:C6-EZ150X2:C6-EZX52:8

0V active, +5V inactive

EZX74 as preferred trigger signal for kV measurement

final high tension on command if all conditions ready

part of: XRG bus, CAN/XS bus

CU CT1 1

cooling unit contact 1 1

EZ150X1:A22-EZX2:6-EWGX4:6=EGWX5:6=EWGX6:6

CU_CT1_2

cooling unit contact 1_2

EZ150X1:C22-EZX2:7-EWGX4:7=EWGX5:7=EWGX6:7

3-56 (b/02.0) OPTIMUS R/F

CU U

stator current U

high speed rotor control units 4512 104 71421/461

EY100 X15

9.3A/V

CU V

stator current V

high speed rotor control units 4512 104 71421/461

EY100 X16

9.3A/V

CU W

stator current W

high speed rotor control units 4512 104 71421/461

EY100 X17

9.3A/V

CV1 EN/

CV2 EN/

converter 1/2 enable

converter 1: EZ130X1:A9-EZX24:22-EQ100X1:22 converter 2: EZ130X1:A30-EZX34:22-E2Q100X1:22

not used, no function

CV1 GND

converter power part 1 ground

EZ130X1:AC8-EZX24:8/21-EQ100X1:8/21

in combination with: CV2_ID/ signal release 2 generators

in combination with: CV2_IDA/ and CV2_IDB/ release 3 generators

CV1_GND_OL

converter power part 1 ground overload (generator basic version ≥ 4512 104 70203/70602)

EZ130X1:A7-EZX24:20-EQ100X1:20

not used, no function

CV1_ID/

converter power part 1 identification EQ100X1:19-EZX24:19-EZ130X1:A6 open +5V, converter connected 0V in combination with: CV1_GND signal

release 2 generators only

CV1_IDA/

converter power part 1 identification A EQ100X1:19-EZX24:19-EZ130X1:A6 open +5V, converter connected +24mV in combination with: CV1_GND signal

release 3 generators only

CV1 IDB/

converter power part 1 identification B EQ100X1:21-EZX24:21-EZ130X1:C9 open +5V, converter connected +24mV in combination with: CV1 GND signal release 3 generators only

CV2 IDA/

converter power part 2 identification A E2Q100X1:19-EZX34:19-EZ130X1:A27 open +5V, converter connected +24mV in combination with: CV2_GND signal release 3 generators only

CV2 IDB/

converter power part 2 identification B E2Q100X1:21-EZX34:21-EZ130X1:C30 open +5V, converter connected +24mV in combination with: CV2 GND signal release 3 generators only

CV1 OL/

converter power part 1 overload EQ100X1:7-EZX24:7-EZ130X1:C7 not used, no function

CV1 TM

converter power part 1 temperature EQ100X1:6-EZX24:6-EZ130X1:C6 4.4V ... 1.5V = 20 ... 100 degrees C in combination with: CV1 GND signal

CV2_GND

converter power part 2 ground EZ130X1:AC29-EZX34:8/21-E2Q100X1:8/21

in combination with: CV2 ID/ signal release 2 generators

in combination with: CV2 IDA/ and CV2 IDB/ release 3 generators

CV2 GND OL

converter power part 2 ground overload (generator basic version ≥ 4512 104 70203/70602) EZ130X1:A28-EZX34:20-E2Q100X1:20

not used, no function

CV2 ID/

converter power part 2 identification E2Q100X1:19-EZX34:19-EZ130X1:A27 open +5V, converter connected 0V in combination with: CV2 GND signal

release 2 generators only

3-58 OPTIMUS R/F (b/02.0)OPTIMUS RF 3 b020

CV2 OL

converter power part 2 overload

E2Q100X1:7-EZX34:7-EZ130X1:C28

not used, no function

CV2_TM

converter power part 2 temperature

EZ130X1:C27-E2Q100X1:6-EZX34:6

4.4V...1.5V = 20...100 °C

in combination with: CV2 GND signal

DR BV 0V

dose rate (signal) reference of image intensifier

EZX61:3-EZ139X2:C18

negative potential of II unit, 0V +/-50mV against generator ground

differential signal with DR BV SG

not used, no function for generators release 2

DR BV NG

dose rate (signal) reference of image intensifier

EZX61:6-EZ139X2:C18

negative potential of II unit, 0V +/-50mV against generator ground

differential signal with DR BV SG

part of: dose rate control

backpanel 4512 108 09361/2 only

backpanel 4512 108 05983/4 only

DR_BV_SG

dose rate signal of image intensifier

EZX61:8-EZ139X2:A18

EZX61:4-EZ139X2:A18

positive potential, 0 ... 10V

differential signal with DR_BV_NG

no function for generators release 2

part of: dose rate control

backpanel 4512 108 05983/4 only backpanel 4512 108 09361/2 only

DR FL LO 1

dose rate fluoro lock-in 1

EWBX12:7-EWB100X1:A21

DR FQ NG

dose rate signal (pulses) negative

not used, no function

DR FQ PO

dose rate signal (pulses) positive

not used, no function

DR_LM

dose rate limiter

EWBX12:1-EWB100X1:A20

low_active if tubelift D76 / EZD on short SID (if tubelift option present)

OPTIMUS R/F (b/02.0) 3-59

OPTIMUS RF 3 b020 © 1997 Philips Medizin Systeme

DR TV NG

dose rate of TV chain signal negative, fluoro control

(II/TV adapter PCB X3:1-X2:8)-EZX61:8-EZ139X2:C19

+/-12V minus polarity

dual voltage differential signal

typically +6V in standby coming from TV chain

+V for more dose, -V for less dose, 0V stable image

part of: dose rate control

DR_TV_NT

dose rate of TV chain signal negative, fluoro control

EZX61:4-EZ139X2:C19

not used, no function

DR_TV_PO

dose rate of TV chain signal positive, fluoro control

(II/TV adapter PCB X3:3-X2:7)-EZX61:7-EZ139X2:A19

-/+12V positive polarity

dual voltage differential signal

typically -6V in standby coming from TV chain

-V for more dose, +V for less dose, 0V stable image

part of: dose rate control

DR TV PT

dose rate of TV chain signal positive, fluoro control

EZX61:9-EZ139X2:A19

not used, no function

DS BV NG

dose (signal ramp) reference of image intensifier

(II/TV_adapter_PCB_X1:P -X2:3)-EZX61:3-EZ139X2:C17

negative potential of II unit, 0V +/-50mV against generator ground

differential signal with DS_BV_SG

part of: dose rate control

DS BV 0V

dose (signal ramp) reference of image intensifier

EZX61:2-EZ139X2C17

not used, no function

backpanel 4512 108 05983/4

backpanel 4512 108 09361/2

backpanel 4512 108 05983/4

backpanel 4512 108 09361/2

backpanel 4512 108 05983/4

backpanel 4512 108 09361/2

DS BV SG

dose signal ramp of image intensifier signal

EZX61:7-EZ139X2:A17

(II/TV_adapter_PCB_X1:R-X2:2)-EZX61:2-EZ139X2:A17

0 ... 10V, polarity positive

differential signal with DS_BV_NG release 3 generators only

release 2 generators: not used, no function

part of: dose rate control

backpanel 4512 108 05983/4 backpanel 4512 108 09361/2

3-60 (b/02.0) OPTIMUS R/F
© 1997 Philips Medizin Systeme

DS MC 0V

dose (signal ramp) reference of selected measuring chamber

EZ150X2:C16-EZ139X2:C16

negative potential of selected measuring chamber, 0V +/-50mV against generator ground differential signal with DS MC SG

DS MC SG

dose signal ramp of selected measuring chamber

EZ150X2:A16-EZ139X2:A16

0 ... +12V

[EZ150 X4 against X5 ground] differential signal with DS_MC_0V

E NG CV1

E value converter DC supply negative

converter 1 (frontal 50/65/80kW): EQ100X1:5-EZX24:5-EZ130X1:C5

0 ... -12V = 0 ... -375V if converter is stand-alone (EQ100 X1 not connected)

if in normal operation: E PO + E NG >> 445VDC = 10V measuring input EZ130 X1:A5 - X1:C5

E_NG_CV2

E value converter DC supply negative

converter 2 (rear 65/80kW): E2Q100X1:5-EZX34:5-EZ130X1:C26

no input to EZ130 release 2 generators

release 3 generators only with 2 converters

0 ... -12V = 0 ... -375V if converter is stand-alone (E2Q100 X1 not connected)

if in normal operation: E PO + E NG >> 445VDC = 10V measuring input EZ130 X1:A26 - X1:C26

E PO CV1

E value converter DC supply positive

converter 1: EQ100X1:18-EZX24:18-EZ130X1:A5

0 ... +12V = 0 ... +375V if converter is stand-alone (EQ100 X1 not connected)

if in normal operation: E_PO + E_NG >> 445VDC = 10V measuring input EZ130 X1:A5 - X1:C5

E PO CV2

E value converter DC supply positive

converter 2: E2Q100X1:18-EZX34:18-EZ130X1:A26

no input to EZ130 version 4512 108 08661..4 release 2 generators

release 3 generators only with 2 converters EZ130 version 4512 108 09102 ... 4

0 ... +12V = 0 ... +375V if converter is stand-alone (E2Q100 X1 not connected)

if in normal operation: E PO + E NG >> 445VDC = 10V measuring input EZ130 X1:A26 - X1:C26

EN X/

enable X-ray, system level

preparation or fluoro request, only valid in combination with CAN message (RAD-R/F)

or hardware requests (Optimus C)

EZ139X1:C2-EZX10:1/3-EZX23:15-EZX45:11-EZX46:11-C300X1:11-EWAX51:11-EWAX52:11

-EWA100X2:C26-EWBX51:11-EWBX52:11-EWB100X2:C26

measuring point: EZX82, EZ139X9

part of: signal bus 0V/+15V low active

EX CD

exposure end signal contact to drive e.g. an external buzzer partner of COM EX CD

EN X C/

enable X-ray, internal generator signal preparation or fluoro request if confirmed by CAN message (RAD-R/F) or hardware requests (Optimus C) EZ119X2:C7-EZ130X1:C7-EZ130X2:C7-EZ139X2:C7-EZ150X2:C7-EZX52:9-EZX76 0V/+5V low active measuring point EZX76 driven by CU if EN X/ active (low)

EX ON

exposure on

part of: XS/XRG bus

EWA100X2:A9-EWAX14:7 EWB100X2:A9-EWBX14:7 potential free optocoupler driven signal

in combination with IT_0V supply: max 26V 10mA part of: EXON old world

FD C CH1

central field measuring chamber 1 EZ150X1:C4-EZX21:12 +15V, Ri of EZ150 = 220Ω

FD C CH2

central field measuring chamber 2 EZ150X1:A4-EZX22:12 +15V, Ri of EZ150 = 220Ω

FD C CH3

central field measuring chamber 3 EZ150X1:C10-EZX31:12 +15V, Ri of EZ150 = 220Ω

FD C_CH4

central field measuring chamber 4 EZ150X1:A10-EZX32:12 +15V, Ri of EZ150 = 220Ω

FD C CH5

central field measuring chamber 5 EZ150X1:C16-EZX41:12 +15V, Ri of EZ150 = 220Ω

3-62 (b/02.0) OPTIMUS R/F

FD_L_CH1 left field measuring chamber 1 EZ150X1:C3-EZX21:11 +15V, Ri of EZ150 = 220Ω

FD_L_CH2 left field measuring chamber 2 EZ150X1:A3-EZX22:11 +15V, Ri of EZ150 = 220Ω

FD_L_CH3 left field measuring chamber 3 EZ150X1:C9-EZX31:11 +15V, Ri of EZ150 = 220Ω

FD_L_CH4 left field measuring chamber 4 EZ150X1:A9-EZX32:11 +15V, Ri of EZ150 = 220Ω

FD_L_CH5 left field measuring chamber 5 EZ150X1:C15-EZX41:11 +15V, Ri of EZ150 = 220Ω

FD_R_CH1 right field measuring chamber 1 EZ150X1:C5-EZX21:3 +15V. Ri of EZ150 = 220Ω

FD_R_CH2 right field measuring chamber 2 EZ150X1:A5-EZX22:3 +15V. Ri of EZ150 = 220Ω

FD_R_CH3 right field measuring chamber 3 EZ150X1:C11-EZX31:3 +15V, Ri of EZ150 = 220Ω

FD_R_CH4 right field measuring chamber 4 EZ150X1:A11-EZX32:3 +15V, Ri of EZ150 = 220Ω FD_R_CH5 right field measuring chamber 5 EZ150X1:C17-EZX41:3 +15V, Ri of EZ150 = 220Ω

OPTIMUS R/F (b/02.0) 3-63

FI TF1 1

filament transformer 1 line 1

EZ119X1:DBZ4-EZX12:1-EG106X15:1

square pulses 100 ... 20kHz, amplitude ~ 300V

FI TF1 2

filament transformer 1 line 2

EZ119X1:DBZ6-EZX12:2-EG106X15:2

square pulses 100 ... 20kHz, amplitude ~ 300V

FI TF2 1

filament transformer 2 line 1

EZ119X1:DBZ8-EZX12:4-EG106X15:4

square pulses 100 ... 20kHz, amplitude ~ 300V

FI TF2 2

filament transformer 2 line 2

EZ119X1:DBZ10-EZX12:5-EG106X15:5

square pulses 100 ... 20kHz, amplitude ~ 300V

GND

ground

- -EZ102X1:DBZ6-EZ119X1:DBZ26-EZ102X2:DBZ8/10/12/14/16/18/20/26/30-EZ119X2:AC4/5/13/15/16/32
 - -EZ130X2:C16:AC4/5/13/15/32-EZ139X2:AC4/5/13/15/32-EZ150X2:AC4/5/13/15/32-EZX21:13
 - -EZX22:13-EZX31:13-EZX32:13-EZX41:13-EZX12:3/6-EZX51:11/12/13/14/15-EZX151:X11/12/13/14/15
 - -EZX44:1/7-EZX46:8/13-EZX1:9-EZX2:10-EZX3:10-EZX5-EZX6-EZX7:3-EZX8:3-EZX17:2-EZX18:2 -EZX19:2-EZX20:2
- -EWGX11:4-EWGX12:4-EWGX1:9-EWGX2:9-EWGX3:9-EWGX4:10-EWGX5:10
 - -EWGX6:10-EWGX7:10-EWGX8:10-EWGX9:10
- -EWAX41:2-EWAX42:2-EWAX51:15-EWAX52:15
 - -EWAX1:7-EWAX2:7-EWAX3:7-EWAX4:7-EWAX11:2-EWAX11:4-EWAX11:6-EWAX11:9-EWAX12:2
 - -EWAX12:4-EWAX12:6-EWAX12:9-EWAX13:9-EWAX14:9-EWAX21:10-EWAX23:10-EWAX24:1
 - -EWAX24:10
- -WA102X1AC2-WA102X2:AC15/28
- -EWBX41:2-EWBX42:2-EWBX51:15-EWBX52:15-EWBX1:7-EWBX2:7-EWBX3:7-EWB4:7
 - -EWBX11:2-EWBX11:9-EWBX12:10-EWBX13:4-EWBX13:6-EWBX21:6-EWBX22:10-EWBX23:10
 - -EWBX24:1-EWBX24:10
- -WB102X1AC2-WB102X2:AC15/28
- -EYAX1:15/16/17-EYAX2:1-EY100X1:11/12/13/14/15-EY100X13-EY100X41
- -C200X1:2-C200X2:17/18/19/20-X100X1:17/18/19/20-C100X10-C100X2:6/7/8/9/10-C300X4:6/7/8/9/10
 -C300X2:1/5
- -EZX87- (cannot be used as signal ground at Duo Diagnost, only Optimus RAD-R/F)

GND 15V

ground (+15V) for desk hand switch

C300X3:1/2/6

3-64 (b/02.0) OPTIMUS R/F

HT AN

high tension anode side actual value

EG100X14:2-EZX35:2-EZ130X1:C17

 $0 \dots +10V = 0 \dots +100 \text{ kV}$ measured at 10kOhm ($20\text{k}\Omega$ measuring circuit parallel to $20\text{k}\Omega$ kV control)

HT_AN_GND

high tension anode side ground

EG100X14:10-EZX35:10-EZ130X1:A17

0V

HT CA

high tension cathode side actual value

EG100X14:1-EZX35:1-EZ130X1:C16

 $0 \dots -10V = 0 \dots -100kV$ measured at 10kOhm ($20k\Omega$ measuring circuit parallel to $20k\Omega$ kV control)

HT CA GND

high tension cathode side ground

EG100X14:9-EZX35:9-EZ130X1:A16

0V

l1 1

partner of I1 1/ optocoupler signal IGBT1 power part 1

EQ100 = 4512 108 05882 release 2

EQ100 ≥ 4512 108 08621 * release 2

EQ100 ≥ 4512 108 09341 * release 3

EZ130X1:A1-EZX24:14-EQ100X1:14

end to X1 * EQ100 X6 measuring point: EQ100 R25 * = X10

value: on = 3.7V off = 1.2Vagainst ground

I1 1/

partner of I1 1 optocoupler signal IGBT1 power part 1

EQ100 = 4512 108 05882 release 2

EQ100 ≥ 4512 108 08621 * release 2

EQ100 ≥ 4512 108 09341 * release 3

EZ130X1:C1-EZX24:1-EQ100X1:1

```
l1 2
partner of I1_2/ optocoupler signal IGBT2 power part 1
                               release 2
EQ100 = 4512 108 05882
EQ100 ≥ 4512 108 08621 *
                               release 2
EQ100 ≥ 4512 108 09341 *
                               release 3
EZ130X1:A2-EZX24:15-EQ100X1:15
I1 2/
partner of I1 2 optocoupler signal IGBT2 power part 1
EQ100 = 4512 108 05882
                               release 2
EQ100 ≥ 4512 108 08621 *
                               release 2
EQ100 ≥ 4512 108 09341 *
                               release 3
EZ130X1:C2-EZX24:2-EQ100X1:2
                                             * EQ100 X7
measuring point: EQ100 R27
                            end to X1
value: ON = 3.7V OFF = 1.2V
                               against ground
                                                * = X10
I1 3
partner of I1 3/ optocoupler signal IGBT3 power part 1
EQ100 = 4512 108 05882
                               release 2
EQ100 ≥ 4512 108 08621 *
                               release 2
EQ100 ≥ 4512 108 09341 *
                               release 3
EZ130X1:A3-EZX24:16-EQ100X1:16
I1 3/
partner of I1 3 optocoupler signal IGBT3 power part 1
EQ100 = 4512 108 05882
                               release 2
EQ100 ≥ 4512 108 08621 *
                               release 2
EQ100 ≥ 4512 108 09341 *
                               release 3
EZ130X1:C3-EZX24:3-EQ100X1:3
                            end to X1
                                             * EQ100 X8
measuring point: EQ100 R29
value: ON = 3.7V OFF = 1.2V
                                                * = X10
                               against ground
I1 4
partner of I1_4/ optocoupler signal IGBT4 power part 1
                               release 2
EQ100 = 4512 108 05882
EQ100 > 4512 108 08621 *
                               release 2
EQ100 ≥ 4512 108 09341 *
                               release 3
EZ130X1:A4-EZX24:17-EQ100X1:17
I1 4/
partner of I1 4 optocoupler signal IGBT4 power part 1
EQ100 = 4512 108 05882
                               release 2
                               release 2
EQ100 ≥ 4512 108 08621 *
EQ100 ≥ 4512 108 09341 *
                               release 3
EZ130X1:C4-EZX24:4-EQ100X1:4
measuring point: EQ100 R31
                            end to X1
                                             * EQ100 X9
value: ON = 3.7V OFF = 1.2V
                               against ground
                                                * = X10
```

3-66 OPTIMUS R/F (b/02.0)© 1997 Philips Medizin Systeme

```
l2 1
partner of I2_1/ optocoupler signal IGBT1 power part 2
                               release 2
EQ100 = 4512 108 05882
EQ100 ≥ 4512 108 08621 *
                               release 2
EQ100 ≥ 4512 108 09341 *
                               release 3
EZ130X1:A22-EZX34:14-E2Q100X1:14
I2 1/
partner of I2 1 optocoupler signal IGBT1 power part 2
EQ100 = 4512 108 05882
                               release 2
EQ100 ≥ 4512 108 08621 *
                               release 2
EQ100 ≥ 4512 108 09341 *
                               release 3
EZ130X1:C22-EZX34:1-E2Q100X1:1
measuring point: EQ100 R25 end to X1
                                             * E2Q100 X6
value: ON = 3.7V OFF = 1.2V
                               against ground
                                                * = X10
12 2
partner of I2 2/ optocoupler signal IGBT2 power part 2
EQ100 = 4512 108 05882
                               release 2
                               release 2
EQ100 ≥ 4512 108 08621 *
EQ100 ≥ 4512 108 09341 *
                               release 3
EZ130X1:A23-EZX34:15-E2Q100X1:15
12 2/
partner of I2 2 optocoupler signal IGBT2 power part 2
EQ100 = 4512 108 05882
                               release 2
EQ100 ≥ 4512 108 08621 *
                               release 2
EQ100 ≥ 4512 108 09341 *
                               release 3
EZ130X1:C23-EZX34:2-E2Q100X1:2
measuring point: EQ100 R27
                            end to X1
                                             * E2Q100 X7
value: ON = 3.7V OFF = 1.2V
                               against ground
                                                * = X10
12_3
partner of I2_3/ optocoupler signal IGBT3 power part 2
EQ100 = 4512 108 05882
                               release 2
EQ100 > 4512 108 08621 *
                               release 2
EQ100 ≥ 4512 108 09341 *
                               release 3
EZ130X1:A24-EZX34:16-E2Q100X1:16
12 3/
partner of I2 3 optocoupler signal IGBT3 power part 2
EQ100 = 4512 108 05882
                               release 2
EQ100 ≥ 4512 108 08621 *
                               release 2
EQ100 ≥ 4512 108 09341 *
                               release 3
EZ130X1:C24-EZX34:3-E2Q100X1:3
measuring point: EQ100 R29 end to X1
                                             * E2Q100 X8
value: ON = 3.7V OFF = 1.2V
                               against ground
                                                * = X10
```

12 4 partner of I2_4/ optocoupler signal IGBT4 power part 2 EQ100 = 4512 108 05882 release 2 EQ100 ≥ 4512 108 08621 * release 2 EQ100 ≥ 4512 108 09341 * release 3 EZ130X1:A25-EZX34:17-E2Q100X1:17 12 4/ partner of I2 4 optocoupler signal IGBT4 power part 2 EQ100 = 4512 108 05882 release 2 EQ100 ≥ 4512 108 08621 * release 2 EQ100 ≥ 4512 108 09341 * release 3 EZ130X1:C25-EZX34:4-E2Q100X1:4 measuring point: EQ100 R31 end to X1 * E2Q100 X9 value: ON = 3.7V OFF = 1.2V against ground * = X10IT 0V emitter 0V exposure on signal EWA100X2:C9-EWAX14:9 EWB100X2:C9-EWBX14:9 potential free optocoupler driven signal in combination with EX ON part of: EXON old world lu stator current phase U of Low Speed Rotor Control measuring point EYAX22 10A/V lw stator current phase W of Low Speed Rotor Control measuring point EYAX21 10A/V MN EM OF mains power emergency off EZX4:1-EZX47:6-EN100X1:6 MN ON mains on Optimus RAD - R/F C300X1:6-EZX46:6-EZX47:2-EN100X1:2-EZX44:14 CB100X10:3-EZX46:6-EZX47:2-EN100X1:2-EZX44:14 Optimus C NG 15V -15V supply Vee EZ102X2:DBZ24-EZ119X2:AC12-EZ130X2:AC12-EZ139X2:AC12-EZ150X2:AC12-EZX21:6-EZX22:6

3-68 (b/02.0) OPTIMUS R/F

-EZX31:6-EZX32:6-EZX41:6-EZX35:15-EZX51:8-EZX151:8-EG100X14:15

-14.5V ... -15.5V

NR PR X/

not ready preparing for X-ray

EZ139X1:A3-EZX23:3-EZX45:4-EZX46:4-C300X1:4-EWAX51:4-EWAX52:4-EWA100X2:A24-EWBX51:4

-EWBX52:4-EWB100X2:A24 driven by CU and/or system controller

measuring point: EZX83 part of: signal bus 0V/+15V high active

PO 0V

signal bus ground GNDS

EZ139X1:AC1-EZX23:1/14-EZX44:15-EZX45:15-EWAX51:15-EWAX52:15-EWBX51:15-EWBX52:15 part of: signal bus, supply via X44 Optimus RAD+R/F, from Cockpit at Duo Diagnost systems

PO_12V

+12 V supply

EN100X1:1-EZX47:1-EZX46:7-C300X1:7

PO 15V

+15V supply Vdd

EZ102X2:DBZ22-EZ119X2:AC11-EZ130X2:AC11-EZ139X2:AC11

- -EZ150X2:AC11-EZX2:8/9-EZX35:7-EZX44:12/13-EZX46:5
- -EZX51:7-EG100X14:7-C300X1:5
- -EZX21/22/31/32/41:5 backpanel 4512 108 05983 only
- -EZX151:7 backpanels 4512 108 05984 + 4512 108 09361/2 only

+14.5V ... +15.5V

PO 15/40V

+15V or +40V supply for measuring chamber

EZ150X1:A20-EZX21/22/31/32/41:5EZ150

version ≥ 4512 108 05964 EZX21/22/31/32/41:5 via (15/40V Sub-D/3+ adapter) EZX21/22/31/32/41:L EZ150 version 4512 108 05963

PO_26V

+26V supply

EZ102X2:DBZ28-EZ119X2:AC14-EZ130X2:AC14-EZ139X2:AC14-EZ150X2:AC14-EZX1:5-EZX2:3

-EZX3:9-EZX11:1-EWGX11:1-EWGX12:1-EZX17:1-EZX18:1-EQ100X2:1-E2Q100X2:1

PO 26V 1

+26V supply options

EZ102X2:DBZ32-EZX19:1-EZX20:1-

- -EWAX1:4-EWAX2:4-EWAX3:4-EWAX4:4-EWAX41:1-EWAX42:1-EWAX23:9-EWAX24:5
 - -EWA100X2:AC14-EWBX1:4-EWBW2:4-EWBX3:4-EWBX4:4-EWBX41:1-EWAX42:1-EWBX21:9
 - -EWBX22:9- EWBX23:9-EWBX24:5-EWB100X2:AC14
 - -EZX8:1 backpanels 4512 108 05984 + 4512 108 09361/2

PO 26V RE

+26V reverse supply

EWAW11-EWAW12-EWAX1/2/3/4:4-EWAX42:1

if generator and system release voltages do not match

normal condition: PO 26V RE = +26V of generator against ground

(jumper WA W11 + W13 closed, W12 open)

special condition: PO 26V RE = 0V against -24V, supply from stand

(jumper WA W11 + W13 open, W12 closed)

PO 26V SW

+26V supply switched, for cooling fan low voltage power supply

EZ102X1:D32-EZX7:1-EM1

backpanels 4512 108 05984 + 4512 108 09361/2

PO 40V

+15V or + 40V supply for measuring chamber

EZ150X1:A20-EZX21/22/31/32/41:5 EZ150

version ≥ 4512 108 05964

EZX21/22/31/32/41:5 via (15/40V Sub-D/3+ adapter) EZX21/22/31/32/41:L EZ150

version 4512 108 05963

PO 400V

+400V supply measuring chamber

EZ150X1:AC1-EZX21/22/31/32/41:1

+400V, Ri of EZ150 = 100kOhms

PO₅V

+5V supply Vcc

EZ102X2:DBZ2/4/6-EZ119X2:AC1/2-EZ130X2:AC1/2-EZ139X2:AC1/2-EZ150X2:AC1/2-EZX46:9-C300X1:9-EZX51:4/5/6-EZX151:4/5/6

+4.74V ... +5.25V

PO V

signal bus supply

EZX23:13/25-EZX44:5-EZX45:7-EZ139X1:AC6

(V15S = -EWAX51:7-EWAX52:7-EWA100X2:AC27-EWBX51:7-EWBX52:7-EWB100X2:AC27)

+15V Vsqn, supply via X44 Optimus RAD+R/F, from Cockpit at DuoDiagnost systems

part of: signal bus

POWERFAIL/

power fail signal of low voltage power supply, initiates warm-boot if supply voltage phase L1 drops below 196VAC EZ102X1:D30-EZ139X1:A10

PW ON NG

relay power on negative, energizes ENK1 if generator ready

EZ130X1:A15-EZX47:9-EN100X1:9

partner of PW_ON_PO

0V/+15V (pulled up by relay coil EN100 K2), low active

3-70 (b/02.0) OPTIMUS R/F
© 1997 Phillips Medizin Systeme

PW ON PO

supply relay power on positive,

EZ130X1:C15-EZX47:4-EN100X1:4

partner of PW_ON_NG

+15V

RC_ON/

rotor control on, low speed rotor control only

EZ150X1:A25-EZX51:1 backpanel 4512 108 05983

EZ150X1:A25-EZX51:1-EZX151:1 backpanels 4512 108 05984 + 4512 108 09361/2

measuring point EYAX28

RC RD/

rotor control ready, low speed rotor control only

EYAX1:9-EXZ51:9-EZ150X1:C25 backpanel 4512 108 05983

EYAX1:9-EXZ51:9-EZX151:9-EZ150X1:C25 backpanels 4512 108 05984 + 4512 108 09361/2

measuring point EYAX25

RC ST 2/

rotor control stator 2

EZ150X1:A26-EZX16:1-EWGX14:1 low speed rotor control EY100X3:1-EWGX14:1 high speed rotor control

RC_ST_3/

rotor control stator 3

EZ150X1:C26-EZX16:2-EWGX14:2-EWGX15:1 low speed rotor control EY100X3:2-EWGX14:2-EWGX15:1 high speed rotor control

RD MN ON

ready mains power on

C100X2:50-C300X4:50-C300X1:14-EZX46:14-EZX47:7-EN100X1:7 Optimus RAD - R/F

CB100X10:4- EZX46:14-EZX47:7-EN100X1:7 Optimus C

RD PR X

NR PR X/

ReaDy preparing for X-ray or Not Ready preparing for X-ray

EZ139X1:A3-EZX23:3-EZX45:4-EZX46:4-C300X1:4- -EWAX51:4-EWAX52:4-EWA100X2:A24

driven by CU or other system components

measuring point: EZX83

part of: signal bus

0V/+15V high active signal

REL CH1

release (reset integrator) chamber 1

EZ150X1:C6-EZX21:4

0V/+15V, typically +13V, high active

REL CH2

release (reset integrator) chamber 2

EZ150X1:A6-EZX22:4

0V/+15V, typically +13V, high active

REL CH3

release (reset integrator) chamber 3

EZ150X1:C12-EZX31:4

0V/+15V, typically +13V, high active

REL CH4

release (reset integrator) chamber 4

EZ150X1:A12-EZX32:4

0V/+15V, typically +13V, high active

REL CH5

release (reset integrator) chamber 5

also used as EXON signal for DSI

EZ150X1:C18-EZX41:4

0V/+15V, typically +13V, high active

RESET 1

external reset

resets incorrect, exposure indication, 5min fluoro buzzer, errors

EWBX22:7-EWB100X1:C23

0V/+26V low active

RESET C/

internal RESET command for function units

EZ119X2:A6- EZ130X2:A6-EZ139X2:A6-EZ150X2:A6-EZX52:3-EZX45:3-EZX46:3-C300X1:3

-EZX51:10-EZX73-EWAX51:3-EWAX52:3-EWA100X1:A6-EWBX51:3-EWBX52:3-EWB100X1:A6-

-EZX151:10 backpanels 4512 108 05984 + 4512 108 09361/2

0V/+5V

measuring point EZX73

driven by CU, active (low) if: EZ139 S1 activated, RESET SW/ ON signal bus active,

threatening power supply drop in, watchdog alarm, switch ON or warm-start,

resets FU's

part of: XS/XRG bus

RESET SW/

signal bus reset, generator reset with turn-ON or push of turn-ON button as warm-start

EZX23:2-EZX44:6-EZ139X1:A2

0V/+15V low active

time constant \geq 200ms

resets CU only

measuring point: EZX81

part of: signal bus

3-72 (b/02.0) OPTIMUS R/F

RF 0V CH1

0V reference value measuring chamber 1

EZX21:8-EZ150X1:C8

differential signal with SIGN_CH1

RF_0V_CH2

0V reference value measuring chamber 2

EZX22:8-EZ150X1:A8

differential signal with SIGN_CH2

RF 0V CH3

0V reference value measuring chamber 3

EZX31:8-EZ150X1:C14

differential signal with SIGN CH3

RF 0V CH4

0V reference value measuring chamber 4

EZX32:8-EZ150X1:A14

differential signal with SIGN CH4

RF_0V_CH5

0V reference value measuring chamber 5

EZX41:8-EZ150X1:C20

differential signal with SIGN_CH5

RG_DV_1

registration device 1 selected

EWA100X1:C4-EWAX1:5

EWB100X1:C4-EWBX1:5

RG_DV_2

registration device 2 selected

EWA100X1:A7-EWAX2:5

EWB100X1:A7-EWBX2:5

RG_DV_3

registration device 3 selected

EWA100X1:A9-EWAX3:5

EWB100X1:A9-EWBX3:5

RG DV 4

OPTIMUS RF 3 b020

registration device 4 selected

EWA100X1:A11-EWAX4:5

EWB100X1:A11-EWBX4:5

RG DV SL 1

registration device selection 1

cassette / camera switchover signal

EWBX21:1-EWB100X1:C18

0V/+26V low active

partner of RG DV SL 2, only one of these should be low active at a time

RG DV SL 2

registration device selection 2

camera / cassette switchover signal

EWBX21:2-EWB100X1:A19

0V/+26V low active

partner of RG DV SL 1, only one of these should be low active at a time

RM_DR_0V

room door contact 0V

EZ150X1:C28-EZX1:10-EWGX1:10-EWGX2:10-EWGX3:10

release 2 generators only, not used release 3 RAD-R/F and Optimus C

partner of RM DR CT signal release 2 RAD generators only

0V/+26V low active, detects room door contact signal short circuit at release 2 RAD generators during turn-ON

RM DR CT

room door contact

EZ150X1:A28-EZX1:8-EWGX1:8=EWGX2:8=EWGX3:8 backpanels 4512 108 05983/4

EZ150X1:A28-EZX45:8-EWBX51:8-EWBX52:8-EWBX22:8-EZX1:8-EWGX1:8=EWGX2:8=EWGX3:8

backpanels 4512 108 09361/2

partner of RM DR 0V signal release 2 RAD generators only

0V/+26V low active = door closed

RQ M1 X/

request mode 1 (fluoro)

Optimus C only, not used

EZX23:9-EZ139X1:C4

RQ M2 X/

request mode 2 (exposure)

Optimus C only, not used

EZX23:22-EZ139X1:C5

RQ M3_X/

request mode 3

Optimus C only, not used

EZX23:10-EZ139X1:C7

3-74 OPTIMUS R/F (b/02.0)

not possible with WA

RQ SN X/

request synchronization of X-ray, exposure request signal

EZX23:16-EZX45:12-EZX46:12-C300X1:12-EZ139X1:C3-EWAX51:12-EWAX52:12-EWA100X2:A25

-EWBX51:12-EWBX52:12-EWB100X2:A25

measuring point: EZX84

0V/+15V

part of: signal bus

RQ XG EX

request X-ray generator for exposure

EWAX1:1- EWAX1:2- EWAX1:3- EWAX1:4-EWA100X1:A3

EWBX1:1- EWBX1:2- EWBX1:3- EWBX1:4-EWB100X1:A3

0V/+26V low active, high if waiting for sync contact

partner of XG RD EX for grid sync (20-21)

RQ XG FL

request X-ray generator for fluoroscopy

EWAX1:6-EWAX2:6-EWAX3:6-EWAX4:6-EWA100X1:A5

EWBX1:6-EWBX2:6-EWBX3:6-EWBX4:6-EWB100X1:A5

0V/+26V low active

RQ_XG_PR_1

request X-ray generator for preparation

EWAX1:3-EWA100X1:A4

EWBX1:3-EWB100X1:A4

0V/+26V low active

RQ XG PR 2

request X-ray generator for preparation

EWAX2:3-EWA100X1:C6

EWBX2:3-EWB100X1:C6

0V/+26V low active

RQ_XG_PR_3

request X-ray generator for preparation

EWAX3:3-EWA100X1:C8

EWBX3:3-EWB100X1:C8

0V/+26V low active

RQ XG PR 4

request X-ray generator for preparation

EWAX4:3-EWA100X1:C10

EWBX4:3-EWB100X1:C10

0V/+26V low active

RX CAN 1

system CAN 1 optional

EZX44:3-EZ139X1:C15

RX CAN 2

system CAN 2 optional

EZX43:1-EZX44:11

S CAN GND

system CAN bus ground

EZ139X1:C17-EZX42:3/6-EZX43:3/6-EZX44:9

-(EZX44:9- EZX44:1- to GND via function programming plug 4512 130 54441 Optimus RAD only)

part of: system CAN

S CAN L

system CAN low active

EZ139X1:C16-EZX42:2-EZX43:2

+2.5VDC standby, +1.5VDC during communication

part of: system CAN

S CAN H

system CAN high active

EZ139X1:A16-EZX42:7-EZX43:7

+2.5VDC standby, +3.2VDC during communication

part of: system CAN

S CAN PO

system CAN supply

EZX44:4-EZX42:9-EZX43:9-EZ139X1:A17

-(EZX44:12-EZX44:9 supply via function programming plug 4512 130 54441 Optimus RAD only)

typically +12V, Vcan

part of: system CAN

SI_PH/

single phase identifier

EN100X1:5-EZX47:5-EZ130X1:C14

SI PH ID

single phase identifier

EN100X1:5-EZX47:5-EZ130X1:C14

SIGN CH1

dose signal of measuring chamber 1

EZX21:7-EZ150X1:C7

0 ... 12V (24V out of range possible)

differential signal with RF 0V CH1

SIGN_CH2

dose signal of measuring chamber 2

EZX22:7-EZ150X1:A7

0 ... 12V (24V out of range possible)

differential signal with RF 0V CH2

3-76 (b/02.0) OPTIMUS R/F

SIGN CH3

dose signal of measuring chamber 3

EZX31:7-EZ150X1:C13

0 ... 12V (24V out of range possible) differential signal with RF 0V CH3

SIGN CH4

dose signal of measuring chamber 4

EZX32:7-EZ150X1:A13

0 ... 12V (24V out of range possible) differential signal with RF 0V CH4

SIGN_CH5

dose signal of measuring chamber 5

EZX41:7-EZ150X1:C19

0 ... 12V (24V out of range possible) differential signal with RF 0V CH5

SL CO 1

select correction 1

external patients size correction, slim patient

EWA100X1:A32-EWAX24:8 EWB100X1:A32-EWBX24:8

0V/+26V low active for selection or when selected from generator desk

SL CO 2

select correction 2

external patients size correction, stout patient

EWA100X1:C32-EWAX24:9 EWB100X1:C32-EWBX24:9

0V/+26V low active for selection or when selected from generator desk

SL_PG_1

select external APRT program 1 EWA100X1:A28-EWAX23:1 EWB100X1:A28-EWBX23:1

0V/+26V low active for selection or when selected from generator desk

SL PG 2

select external APRT program 2 EWA100X1:C28-EWAX23:2 EWB100X1:C28-EWBX23:2

0V/+26V low active for selection or when selected from generator desk

SL PG 3

select external APRT program 3 EWA100X1:A29-EWAX23:3 EWB100X1:A29-EWBX23:3

0V/+26V low active for selection or when selected from generator desk

AIL RIGHTS RESERVED

FAULT FINDING OPTIMUS R/F

SL PG 4

select external APRT program 4 EWA100X1:C29-EWAX23:4 EWB100X1:C29-EWBX23:4

0V/+26V low active for selection or when selected from generator desk

SL_PG_5

select external APRT program 5 EWA100X1:A30-EWAX23:5 EWB100X1:A30-EWBX23:5

0V/+26V low active for selection or when selected from generator desk

SL_PG_6

select external APRT program 6 EWA100X1:C30-EWAX23:6 EWB100X1:C30-EWBX23:6

0V/+26V low active for selection or when selected from generator desk

SL PG 7

select external APRT program 7 EWA100X1:A31-EWAX23:7 EWB100X1:A31-EWBX23:7

0V/+26V low active for selection or when selected from generator desk

SL PG 8

select external APRT program 8 EWA100X1:C31-EWAX23:8 EWB100X1:C31-EWBX23:8

0V/+26V low active for selection or when selected from generator desk

SL_TO_TM_1

select tomo time 1

tomo time input from stand EWAX21:1-EWA100X1:A24

0V/+26V low active

SL_TO_TM_2

select tomo time 2

tomo time input from stand

EWAX21:2-EWA100X1:C24

0V/+26V low active

SL TO TM 3

select tomo time 3

tomo time input from stand EWAX21:3-EWA100X1:A25

0V/+26V low active

3-78 (b/02.0) OPTIMUS R/F

SL_TO_TM_4

select tomo time 4

tomo time input from stand

EWAX21:4-EWA100X1:C25

0V/+26V low active

SL_TO_TM_5

select tomo time 5

tomo time input from stand

EWAX21:5-EWA100X1:A26

0V/+26V low active

SL TO TM 6

select tomo time 6

tomo time input from stand

EWAX21:6-EWA100X1:C26

0V/+26V low active

SL_TO_TM_7

select tomo time 7

tomo time input from stand

EWAX21:7-EWA100X1:A27

0V/+26V low active

SL TO TM 8

select tomo time 8

tomo time input from stand

EWAX21:8-EWA100X1:C27

0V/+26V low active

SL XG TO

select X-ray generator for tomography

EWAX11:3-EWAX12:3-EWA100X1:C18

0V/+26V, low active

STOP_X_C/

stop X-ray command, X-ray OFF from function units mA and dose rate control (on-board of CU)

EZ119X2:A7-EZ130X2:A7- EZ139X2:A7-EZ150X2:A7-EZX52:4

0V/5V

measuring point EZX75

inactivates CTRL X C/

EXOF exposure OFF command

part of: XS/XRG bus

FAULT FINDING OPTIMUS R/F

STU

stator phase U

EYAX2:2-EX1101 low speed rotor control single tube

EYAX2:2-EWGK11:1-EWGK12:1=EWGK11:2=EWGK12:2

EY100X46:2-EX1101

low speed two tubes high speed rotor control

vers. 4512 104 33791/2 or 71401..6 single tube

EY100X46:2-EWGK11:1-EWGK12:1=EWGK11:2=EWGK12:2 high speed rotor control

vers. 4512 104 33791/2 or 71401..6 two tubes

EY100X51-EX1101 high speed rotor control

vers. 4512 104 71421/61 single tube

EY100X51--EWGK11:1-EWGK12:1=EWGK11:2=EWGK12:2 high speed rotor control

vers. 4512 104 71421/61 two tubes

STV

stator phase V = common

EYAX2:3-EX1102

EYAX2:3-EWGK11:3-EWGK12:3=EWGK11:4=EWGK12:4

EY100X47:1-EX1102

EY100X52-EX1102

low speed rotor control single tube

low speed two tubes high speed rotor control

vers. 4512 104 33791/2 or 71401..6 single tube

EY100X47:1-EWGK11:3-EWGK12:3=EWGK11:4=EWGK12:4 high speed rotor control

vers. 4512 104 33791/2 or 71401..6 two tubes

high speed rotor control

vers. 4512 104 71421/61 single tube

EY100X52--EWGK11:3-EWGK12:3=EWGK11:4=EWGK12:4 high speed rotor control

vers. 4512 104 71421/61 two tubes

STW

stator phase W

EYAX2:4-EX1103

EYAX2:4-EWGK11:5-EWGK12:5=EWGK11:6=EWGK12:6

EY100X47:2-EX1103

low speed rotor control single tube

low speed two tubes

high speed rotor control

vers: 4512 104 33791/2 or 71401..6 single tube

EY100X47:2-EWGK11:5-EWGK12:5=EWGK11:6=EWGK12:6 high speed rotor control

vers. 4512 104 33791/2 or 71401..6 two tubes

EY100X53-EX1103 high speed rotor control

vers. 4512 104 71421/61 single tube

EY100X53--EWGK11:5-EWGK12:5=EWGK11:6=EWGK12:6

high speed rotor control

vers. 4512 104 71421/61 two tubes

SW BU 1

switch bucky 1 ready (WA + WB)

EWAX11:10-EWA100X1:C19

EWBX11:10-EWB100X1:C19

part of: bucky ready contact

0V/+26V low active

SW BU 2

switch bucky 2 ready (WA only)

EWAX12:10-EWA100X1:A21

part of: bucky ready contact

0V/+26V low active

3-80 (b/02.0)OPTIMUS R/F © 1997 Philips Medizin Systeme

SW OF FD 1

switch OFF field 1

format size correction < 14cm or if cone in use serial changer chamber

EWBX13:5-EWB100X1:C21

0V/+26V low active

SW ON FD 3

switch ON field 3

format size correction > 24x24cm serial changer chamber

EWBX13:7-EWB100X1:A22

0V/+26V low active

SW PR FL 1

switch preparation or fluoro 1

contact to drive an external prep or fluoro indication lamp

EWBX22:2-EWB100X1:C13

partner of COM_EX_CD

SW SF CF 1

switch side field to central field bucky measuring chamber (WA + WB)

EWAX11:1-EWA100X1:A18

EWBX11:1-EWB100X1:A18

cassettes < 23cm

0V/+26V low active

SW SF CF 2

switch side field to central field bucky measuring chamber 2 (WA only)

EWAX12:1-EWA100X1:A20

cassettes < 23cm

0V/+26V low active

SW_TO_1

switch tomography 1 ready

EWAX11:5-EWA100X1:A19

part of: tomo ready contact

0V/+26V low active

SW TO 2

switch tomography 2 ready

EWAX12:5-EWA100X1:C20

part of: tomo ready contact

0V/+26V low active

SW_UN_EX

radiation indication

EZ150X1:A29-EZX1:4-EWGX1:4

partner of CM_SW, potential free contact

FAULT FINDING OPTIMUS R/F

SW UN EX 1

radiation indication

(EWGX1:4)=EWGX2:4

partner of CM SW, potential free contact

SW_UN_EX_1

switch radiation indication 1

contact to drive an external X-ray indication lamp

EWBX22:4-EWB100X1:C14 partner of COM EX CD

SW UN EX 2

radiation indication

(EWGX1:4)=EWGX3:4

partner of CM_SW, potential free contact

SW XG RD 1

switch generator ready 1

contact to drive an external ready indication lamp

EWBX22:1-EWB100X1:A13 partner of COM EX CD

SW_WN_FL_1

switch warning fluoro 1

contact to drive an external fluoro warning indication lamp (> 5 minutes)

EWBX22:3-EWB100X1:A14 partner of COM_EX_CD

TB 2/

tube 2 selected

EZ130X1:A13-EZX11:2-EWGX11:2

0V/15V, low active

TB 2 RT

tube 2 return signal, tube selection check

EWGX11:3-EZX11:3-EZ130X1:A10

0V/5V, low active

TB 3/

tube 3 selected

EZ130X1:C13-EZX11:5-EWGX11:5-EWGX12:2

0V/15V, low active

TB_3_RT

tube 3 return signal, tube selection check

E2WGX11:3-E1WGX12:3-E1WGX11:6-EZX11:6-EZ130X1:C10

0V/5V, low active

3-82 (b/02.0) OPTIMUS R/F

TB CU FR NG

tube current frequency negative

EG100X14:14-EZX35:14-EZ119X1:BZ32

-14V against ground, frequency: 1 kHz = 2mA, 0 ... 1500mA 500kHz/A

differential signal with TB_CU_FR_PO

TB CU FR PO

tube current frequency positive

EG100X16:6-EZX35:6-EZ119X1:BZ30

-14V against ground, frequency: 1 kHz = 2mA, 0 ... 1500mA 500kHz/A

differential signal with TB CU FR NG

TH OL

tube housing overload

NTC temperature measurement in tube housing (not yet available)

EZ130X1:A12-EZX3:6-EWGX7:6-EWGX8:6-EWGX9:6

backpanel 4512 108 05983

EZ130X1:A12-EZX3:3-EWGX7:3-EWGX8:3-EWGX9:3

backpanels 4512 108 05984 + 4512 108 09361/2

4.4V ... 1.5V = 20 ... 100 degrees C

partner of CM TH

TH_OL_SW/

tube housing overload switch

EZ130X1:A11-EZX3:3-EWGX7:3-EWGX8:3-EWGX9:3

backpanel 4512 108 05983

EZ130X1:A11-EZX3:6-EWGX7:6-EWGX8:6-EWGX9:6

backpanels 4512 108 05984 + 4512 108 09361/2

0V ... 1.7V = short circuit, 1.7V ... 3.3V = closed, >3.3V open

partner of CM_TH_SW

TOMO PG

tomo mode programmed

EWA100X1:A17-EWAX22:9

common line for tomo trajectory selection TO PG 1 ... 8 to stand, potential free

TO PG 1

tomo program 1

EWA100X1:A13-EWAX22:1

tomo trajectory selection, potential free contact with TOMO PG

TO PG 2

tomo program 2

EWA100X1:C13-EWAX22:2

tomo trajectory selection, potential free contact with TOMO_PG

TO_PG_3

tomo program 3

EWA100X1:A14-EWAX22:3

tomo trajectory selection, potential free contact with TOMO_PG

FAULT FINDING OPTIMUS R/F

TO PG 4

tomo program 4

EWA100X1:C14-EWAX22:4

tomo trajectory selection, potential free contact with TOMO_PG

TO_PG_5

tomo program 5

EWA100X1:A15-EWAX22:5

tomo trajectory selection, potential free contact with TOMO PG

TO PG 6

tomo program 6

EWA100X1:C15-EWAX22:6

tomo trajectory selection, potential free contact with TOMO_PG

TO PG 7

tomo program 7

EWA100X1:A16-EWAX22:7

tomo trajectory selection, potential free contact with TOMO_PG

TO PG 8

tomo program 8

EWA100X1:C16-EWAX22:8

tomo trajectory selection, potential free contact with TOMO PG

TO_PG_SL

tomo program selected

EWA100X1:C17-EWAX22:10

tomo APR selected = closed, overriding = open, potential free contact with TOMO PG

TP_HT_GND

temperature high tension tank ground

EZ130X1:A19-EZX35:12-EG100X14:4

partner of TP_HT_SG

TP HT SG

temperature signal high tension tank

NTC in high tension tank oil

EG100X14:12-EZX35:4-EZ130X1:C19

4.4V ... 1.5V = 20 ... 100°C

+25 C(12kW) ... +100 C(950W)

partner of TP HT GND

V15C

(S_CAN_PO)

system CAN supply

EZX42:9-EZX43:9-EZX44:4-EZ139X1:A17

Vcan

part of: system CAN

backpanel 4512 108 05983 only

V15S

signal bus supply

backpanel 4512 108 05983 only

EZX23:13/25-EZX44:5-EZX45:7-EZ130X1:AC6-EWAX51:7-EWAX52:7-EWA100X2:AC27

+15V Vsgn

part of: signal bus

VO CR IF 0

density voltage correction II format dependent 10"

EWBX13:3-EWB100X1:C22

0V/+26V low active

VO CR IF 1

density voltage correction II format dependent 5" / 6"

EWBX13:9-EWB100X1:A23

X ACT/

X-ray active signal bus

EZ139X1:A5-EZX23:5-EZX45:6-EWAX51:6-EWAX52:6-EWA100X2:C24-EWBX51:6-EWBX52:6

-EWB100X2:C24

driven by CU if X_ACT_S/ was sent from FU-kV or during fluoro, old: EXON signal

measuring point: EZX86

part of: signal bus

0V/+15V

X_ACT_S/

X-ray active signal

kV > 75% nominal value driven by FU-kV or fluoroscopy high tension on driven by CU

EZ119X2:A8-EZ130X2:A8-EZ139X2:A8-EZ150X2:A8-EZX52:5-EZX77

0V/+5V

measuring point EZX77

part of: XS/XRG bus, controls X_ACT/ status

XG_RD_EX_1

X-ray generator ready for exposure request

grid / sync release signal

EWA100X1:C3-EWAX1:2

EWAB100X1:C3-EWBX1:2

0V/+26V low active

partner of RQ_XG_EX for grid sync (20-21)

XG RD EX 2

X-ray generator ready for exposure request

grid / sync release signal

EWA100X1:A6-EWAX2:2

EWB100X1:A6-EWBX2:2

0V/+26V low active

partner of RQ_XG_EX for grid sync (20-21)

XG_RD_EX_3

X-ray generator ready for exposure request

grid / sync release signal

EWA100X1:A8-EWAX3:2

EWB100X1:A8-EWBX3:2

0V/+26V low active

partner of RQ_XG_EX for grid sync (20-21)

XG_RD_EX_4

X-ray generator ready for exposure request

grid / sync release signal

EWA100X1:A10-EWAX4:2

EWB100X1:A10-EWBX4:2


0V/+26V low active

partner of RQ_XG_EX for grid sync (20-21)

3-86 (b/02.0)OPTIMUS R/F OPTIMUS_RF_3_b020

13. Optimus AEC switch-OFF philosophy

For explanation the following figure:



Every APR using the AEC technique as the preferred technique must have mAs, mA-s or mAs-s parameters in the background. These should almost match the typical organ-related dose to the selected film/screen combination.

A film which got at least 40% of the desired density can be used for diagnosis.

FAULT FINDING OPTIMUS R/F

If the AEC exposure starts now two supervisions are active with the aim of not giving unnecessary dose (or simply a proper AEC exposure is obtained):

- 1. The organ-dependent background mAs value is multiplied with 9.5. If the exposure is not finished at 9.5 x mAs_{backup} the generator stops. One must expect that something went wrong if the exposure exceeded 9.5 times the typical mAs value. This exposure has not been cut OFF by the supervision 2.
- 2. With the 9.5 x mAs_{backup} a kV and filament load dependent backup time is calculated by DRC (dose rate control). At 10% of this time value DRC checks if at least 4% of the desired dose has been detected by the measuring chamber.

If the 4% limit does not increase, the exposure switches OFF. The minimum of 40% density cannot be obtained during the remaining backup time.

This 4% dose detection is automatically OFF, if the film/screen combination is too sensitive (> 400 speed systems). The 4% value is too small to be reliable for a measurement.

With overriding the supervision switches OFF.

How to test the limits of 600mAs or 4000ms in AEC technique

One has to bypass the 4% detection and the background mAs value must be high enough to reach 600mAs. The 4% detection can be switched OFF with modifying the value dose of FSC [μGy]:

- Type in a value of 1 (which is equal to a 1000 speed system) in the dose of FSC data field of any of the programmed film/screen combinations.
- Now select any APR and increase the background mAs value to 100mAs.
- Close the collimator or cover the chamber with lead.

The AEC exposure stops at a value which is always below 600mAs, a typical limit is 588mAs.

With the modified parameters the 4000ms test can be carried out:

- Select the modified APR on the control desk and go to <SELECT APR> and <CHANGE APR> with the PC.
- Reduce the le max factor to 5% and transmit the APR screen.
- Select the APR button again, the modified data are active now.
- · Select the small focal spot.
- Switch an AEC exposure. It should last 4000ms.
- Change all modifications back to normal.

The supervision can be switched ON or OFF, programming path:

XRGSCOPE - Optimus (XRG90) - Programming - Dose Rate Control - AMPLIMAT - Fault Exposure Detection AEC or TDC - ON/OFF

(Explanation see documentation.)

Precalculation tables of the exposure which is actually displayed on the control desk can be seen on the PC under:

XRGSCOPE - Optimus (XRG90) - Faultfind - X-Ray Log - Dose Rate Control Logging - etc.

14. AEC fault exposure detection strategy

The major intention of having a fault exposure detection is to prevent unnecessary radiation for the patient in case of a malfunction of the installation or a mistake when handling the X-ray equipment.

AEC fault exposure detection = ON

The factors determining whether the 4% dose value at 10% of the APR backup time are checked are

- the 10% backup time value > 10ms

- the expected 4% density voltage value > 20mV

In case of APR100 the check could be performed because the density voltage values are high enough.

The density voltage at 10% of the backup time would be too small to be measured for APR800, therefore the exposure continues up to the 9.5 x APR mAs value. The exposure finally terminates at 570mAs if the APR mAs value is \geq 60mAs.

With APR100* the exposure terminates at 10% of the max backup time, which is 4000ms for all AEC exposures after overriding of any APR parameters. (The 600mAs limit does not switch OFF the exposure, 1500mA emission current is not available).

With APR800* the exposure terminates either at 600mAs or 4000ms, depending on which of the limits is reached first.


AEC fault exposure detection = OFF

APR100 and APR800 exposures have the same termination point at 9.5 x APR mAs. The exposure finally terminates at 570mAs, if the APR mAs value is \geq 60mAs.

APR100* and APR800* exposures terminate either at 600mAs or at 4000ms, depending on which of the limits is reached first.

For explanations see figure and list of terms on following 2 pages.

FAULT FINDING OPTIMUS R/F



List of terms

AEC =	Automatic Exposure Control
-------	----------------------------

APR100 = APR program with a less sensible film/screen combination of 100 speed, original parameters

as programmed.

APR100* = as same as APR100, but parameter(s) modified on the control desk (overriding).

APR800 = APR program with a very sensible film/screen combination of 800 speed, original parameters

as programmed.

APR800* = as same as APR800, but parameter(s) modified on the control desk (overriding).

600mAs = programmed mAs limit for AEC exposures (can be changed, must comply with the local

regulations).

4000ms = max time limit of AEC exposures (cannot be changed).

1 = point of the mAs dependent APR backup time, which is calculated from the

9.5 x (typical) organ mAs value of the APR

2 = max mAs limit for AEC exposures (can be changed)

3 = max exposure time limit of 4000ms (cannot be changed)

4 = 10% (of the APR backup) time point

5 = 10% backup time point of the max exposure time limit (4000ms) = always 400ms

To explain the difference in switching the fault exposure detection **ON** or **OFF**, a very sensible (800 speed system) and a less sensible (100 speed system) film/screen combination have been chosen.

FAULT FINDING OPTIMUS R/F

15. Explanations on programming the generator

Select menu: OPTIMUS\ PROGRAM\ DOSE RATE CONTROL\ CONT ...

scantime_TV [ms]: [20.000] = default 20ms

20ms must be programmed for all TV chains with a scantime ≤ 20ms.

In case a TV chain has a longer scantime program the actual scantime value:

for scantime of TV chain ≤ 20ms scantime_TV [ms] = [20.000]

scantime TV [ms] = [xx.xxx] for scantime of TV chain = xx.xxx ms > 20ms

scantime_TV valid: YES (default)


With YES, the programmed variable scantime TV [ms] is used for the control in FU Dose Rate Control (DRC).

At the moment only this value is possible.

With NO, the programmed variable scantime TV [ms] is not used for the control in FU Dose Rate Control (DRC). The variable is delivered by the system controller.

This version is a future option. At the moment it is not applicable.

See figure explaining the program settings below:



AIL RIGHTS RESERVÉD

16. Printed-circuit boards

Low-voltage power supply: EZ102


Also see Z1-2.3 "Low-voltage power supply".

LEDs H2 to H5 indicate whether the supply voltages are present.

The low-voltage power supplies of PCB EZ102 are short circuit proof. Therefore it is most likely that in case one of the LEDs grows dark one of the external consumers and not the PCB itself is the cause of the error.


It is recommended that one after the other all consumers be disconnected from the respective power supply until the LED is illuminated again.

The last consumer that was removed has probably caused the short-circuit.





Central rack, service aid

15V 15V 15V - 24V 15V 15V - 24V 2 2 2 - 5\ 0 0 EZX83 EZX76 EZX85 EZX82 EZX84 measuring points EZX8x only backpanel > = 4512 108 05984 and 4512 108 0936x (23) [8] (50)(19) 4 Signals during an exposure released from the desk handswitch with Bucky Controller (signal bus EZX23) (16) -goes low active even for a free cassette RGDV4, has no function in the BuCo for free cassette (15)(4 4 -triggers the grid release (21) in case of a table/wall bucky or Tomo RGDV1/2/3 (13)the EN_X/ on the signal bus must be confirmed by the CAN message from the desk CPU to be valid signal pulled up to + 15V, driven to low in case of NOT ready -all RGDV 9 တ ω in case of all non AEX exposures ဖ 2 Dose Rate Control in case of AEC exposures (1)(2)with backpanel 4512 108 05983 and CU 4512 108 0902x (2) with backpanel 4512 108 05984 and CU 4512 108 0902x (2) with backpanel 4512 108 0936x and CU 4512 108 0902x 4 acceleration, filament heating က when activated by user RD PR X (1) must be HIGH to get X-ray started NR_PR_X/(2) signal driven by BuCo and/or OPTIMUS A4 97-12-18 Sie./Re READY_EXP_P_T(BuCo) Ŋ SYNC2_T(BuCo) optimus_r_f_3z-10 STOP_X_C/ CTRL_X_C/ RQ SN X/ X_ACT_S/ EN_X_C/ CTRL_X/ EXRQ(21) X_ACT/ EXON EN_X EXOFF PRRQ (98.0)3Z-10 OPTIMUS R/F ©1996 Philips Medical Systems ALL RIGHTS RESERVED



ALL RIGHTS RESERVED





Comparison release decades -CP generators< - \(\text{PPTIMUS} \)

Commands given via release decades not desk handswitch

REPLACEMENT

TEXT

	Contents	4-0.1
1.	Traceable items	4-1
1.1.	Generator cabinet	4-1
1.2.	H.V. tank	4-2
1.3.	Operating desk	4-2
2.	Printed-circuit boards	4-3
3.	Exchange of firmware or update to release 3.5:	4-5
3.1.	PC and generator settings to avoid problems during uploading of CU complete files	4-5
3.1.1.	Preparation of the service PC to guarantee a safe loading process	4-6
3.1.2.	Preparation of the generator	4-6
3.2.	Backup / Installation procedure	4-7
3.2.1.	Backup of all configuration data	4-9
3.2.2.	Loading the new firmware into the generator	4-9
4.	Replacement of parts of "Function Unit kV"	4-10
4.1.	General information	4-10
4.2.	Connecting and setting the scope	4-11
4.3.	Deactivating the kV controller	4-12
4.4.	Setting of exposure data	4-12
4.5.	Adjustment of the "factor for duty cycle"	4-13
5.	Tube replacement	4-16
5.1.	Tube conditioning	4-16
5.1.1.	Preconditions / Programmings	4-16
5.1.2.	Procedure	4-17
5.2.	Tube adaptation	4-21
5.2.1.	General information	4-21
5.2.2.	Preconditions / Programmings	4-22
5.2.3.	Procedure	4-23
5.3.	Final tube adjustment work	4-24
5.4.	Symptoms and solutions if problems occur during tube adjustment	4-25


OPTIMUS R/F REPLACEMENT

1. Traceable items

Trace items are:

- 1. Generator cabinet
- 2. H.V. tank
- 3. Operating desk

They are labelled as follows:



With new traceable items for replacement a separate label is delivered.

This must be affixed to the label bracket on the top left corner of the generator cabinet. See drawing 2Z-10 "Labelling".

The new type number, code number and serial number must be entered on the master card for the generator.

Please, send a copy of the corrected master card as FAX to:

Philips Medical Systems
DMC Hamburg, Germany
Department: GEN-OPERATION
FAX No.: +49 40 5078 1247

Generator cabinet

The generator cabinet as a traceable item is labelled by a 6-digit serial number:

Example:

1.1.

s/n **01 1234**

Meaning:

01 = year of manufacture, e.g. 20**01** 1234 = consecutive number REPLACEMENT OPTIMUS R/F

1.2. H.V. tank

H.V. tanks have a 7-digit serial number which has the following meaning:

Example:

s/n 01 04 123

Meaning:

```
01
                  year of manufacture, e.g. 2001
   04
                  power class, e.g. 65/80kW, 2 tubes
                  consecutive number
      123
```

Power classes:

50kW, 1 tube 01 02 50kW, 2 tubes 03 65/80kW, 1 tube = 65/80kW, 2 tubes 04

Caution!

An exchange of a H.V. tank requires a new alignment of "Function Unit kV". For alignment work refer to chapter 4 in this section.

1.3. Operating desk

The operating desk is labelled by a 8-digit serial number:

Example:

s/n 01 02 1234

Meaning:

01 year of manufacture, e.g. 2001 internal number of subcontractor 02 1234 =consecutive number

OPTIMUS R/F REPLACEMENT

Printed-circuit boards 2.

РСВ	HW programming	SW programming via XRGSCOPE	Tube adaptation	Remarks	
EZ backpanel	• see Z2-5.1/.2/.3			To attend to: X4 emergency OFF X10 EN_X/ X42 system CAN termination X44 function programming plug X45 generator CAN termination X52 shall not be present	
EZ102 low voltage supply					
EZ119 mA control	 exchange PROM or insert new PROM set battery jumper to ON see 5Z-1 	load tube data set(s)	all tubes	Bucky TH, Digital Diagnost, Thoravision systems: Set RGDV according to adaptation. See section 2, chapter 8.3.2.	
EZ130 kV control	 exchange PROM or insert new PROM see 5Z-1 			Carry out alignment of "Function Unit kV".	
EZ139 CU	 exchange BOOT PROM or insert new PROM see 5Z-1 	set date and time restore CU complete or start programming from beginning		Carry out alignment of "Function Unit kV" if no CU complete files are present.	
EZ150 basic interface	• see 5Z-1	check AMPLIMAT sensitivity according to jumper W4		 Set jumpers W2 + W3 according required chamber supply. Set jumper W4 according to programmed AMPLIMAT sensitivity. 	
EN100 power ON circuit					
EG100 measuring circuit				Exchange is not allowed. Requires alignment which is not possible in the field. • Exchange the whole tank.	

OPTIMUS R/F REPLACEMENT

РСВ	HW programming	SW programming via XRGSCOPE	Tube adaptation	Remarks
EWA 1WA / 2WA bucky-tomo backpanel EWB R/F adapter	Address W1W3Supply + ground W11W13			
backpanel WA102	• see Z1-15.1			Can be used in:
universal I/O	• see 5Z-2			WA / 1WA / 2WA / WB / 1WB
EY100 rotor control high speed	swap PROM or insert new PROMsee 5Z-2			
EYA100 rotor control low speed				
C300 desk CPU	swap PROM or insert new PROMsee 5Z-2			

OPTIMUS R/F REPLACEMENT

3. Exchange of firmware or update to release 3.5:

NC: 9890 000 0251x Firmware OPTIMUS Rel. 3.5 4512 114 2083x Central Unit

The CU firmware is no longer available as an EPROM (EZ139 D4/D5).

It must be loaded from the PC into the respective flash PROMs.

For loading firmware release 3.5 the following firmware levels must be present in the generator:

- CU-Boot : ≥ 4512 113 2073x

- FU-kV : \geq 4512 113 2013x OPTIMUS RAD

- FU-kV : ≥ 4512 113 2621x OPTIMUS R/F

- FU-mA : \geq 4512 113 2021x - FU-CIE : \geq 4512 113 2031x - FU-HI : \geq 4512 113 2052x - FU-Adap : \geq 4512 113 2061x - FU-RoCo : \geq 4512 113 2232x

Caution!

Before changing the release, save all configuration data of the generator! Refer to chapter 3.2.1 "Backup of all configuration data".

3.1. PC and generator settings to avoid problems during uploading of CU complete files

Optimus R/F release 3.x CMOS data are uploaded in one string without handshake.

Any kind of interruption can cause the loading process to fail.

Problems occur mainly during the download to the PC.

A download file which is not complete cannot be used as a safety backup file.

REPLACEMENT **OPTIMUS R/F**

3.1.1. Preparation of the service PC to guarantee a safe loading process

Start XRGSCOPE always from DOS if possible.

In case of any WINDOWS version:

- · Switch OFF all screensavers.
- · Do not run other programs.
- Do not insert any CD in the drive.

Any kind of power management of the PC hardware (BIOS) as well as the windows power management should be switched OFF.

If the PC is connected to mains power some of these might be automatically OFF.

3.1.2. Preparation of the generator

Preparation of generators without a CAN interface:

· Switch OFF the generator.

Preparation of generators which are connected via a CAN interface:

- BuckyDiagnost TH and TH2
- DigitalDiagnost
- Thoravision
- EasyDiagnost with bucky unit
- Switch OFF the generator.
- Disconnect the following plugs:


System	Connector				
	EZX23 signal bus	EZX42 or EZX42-1 system CAN	EZX43 or EZX43-1 system CAN		
BukyDiagnost TH / TH2	Χ		X		
DigitalDiagnost	Х	Х	Х		
Thoravision	Х	Х	Х		
EasyDiagnost with Bucky unit	Х	Х	Х		

OPTIMUS R/F REPLACEMENT

3.2. Backup / Installation procedure

Provide the service PC with the hardware key and switch it ON.
 The hardware key provides access to special program settings and to menu "Faultfind".
 Standard programming is possible without a hardware key.

 Connect the PC to X5 on EZ139 CENTRAL UNIT CU via a serial data cable: (A 5m long data cable can be ordered via 12NC: 4512 130 56931)



• Insert the floppy disks containing the self-unpacking exe-files of the firmware in the disk drive of the PC:

OMA: 4512 116 022xxOMB: 4512 116 023xx

For unpacking on the harddisk of the PC about 5MB are needed.

- Generate a directory e.g. [C:\OPT_R_35] on the PC by typing <md C:\OPT_R_35> or use WIN commands.
- Copy the firmware from both floppy disks to the PC into the same directory C:\OPT_R_35 by typing <copy A:*.* C:\OPT_R_35> or use WIN commands.
- Start unpacking the programs by typing <OMAxxxxx.exe> and <OMBxxxxx.exe> or doubleclick on <*.exe> -files.

The programs unpack all files needed for the update of the firmware and the newest service tools.

- After unpacking [OMAxxxxx.exe] and [OMBxxxxx.exe] can be deleted on the harddisk by typing
 del OMAxxxxx.exe> and <del OMBxxxxx.exe>.
- For the current contents [of OMA/OMBxxxxx.exe] read [OMA/OMBxxxxx.txt] on the floppy disks.
- For further installation of the firmware in the generator read README.txt.

REPLACEMENT OPTIMUS R/F

General information:

- Button <F1> <help> Call help / cancel help.

- Button **<F2> <transmit>** Store screen contents / data set in the generator ==> transmit to generator.

Button <F3> <save> Store data screen on disk.

For an open data screen the path desired can be selected.

- Button <F4> <load> Load data set from disk. The desired path can be selected.

Button **<ESC>** Commands one step back. Can be used repeatedly.


Fields with ↓ Select the possible range of values by pushing <RETURN>.

The data are specified by the generator as fixed values.

- Fields with [...] Input of data via the keyboard.

Error numbers which appear at the beginning of the programming procedure must be erased from the screen by pushing the **<RETURN>** key.

Current data files, for instance, for online help, tube types, APR programming are available in the PHILIPS-Intranet. Use path *http://technet.best.ms.philips.com*/ and pull down menu as shown below.



If you call the installation program with **<xrgscope** ?> the possible starting parameters for the service program are listed.

OPTIMUS R/F REPLACEMENT

3.2.1. Backup of all configuration data

· Switch ON the generator.

To save the configuration data use the CONFIGURATION BACKUP disk supplied.

Save the complete SW programming of the generator on the floppy disk by using the menu:
 ACCEPT/ BACKUP/ CU complete

A disk space of 700kByte is required.

It takes about 8min to save the data to the disk.

The default backup name:

Backup File Name: CUBACKUP.TDL

The name can be changed into any other filename.

The path (harddisk) is automatically taken into account.

It is also possible to type:

A:\"filename" <RETURN>

to load the backup files directly to the floppy disk.

3.2.2. Loading the new firmware into the generator

- · Switch OFF the generator.
- · When installing the firmware

from COM 1 enter FLASH1.BAT from COM 2 enter FLASH2.BAT "Attempting link" appears on the screen.

· Switch ON the generator.

Depending on the type of PC data transmission takes 15 ... 30min. During this process all red LEDs of the function unit are blinking.

Caution!

When the data transmission to the generator is completed, the scope program is still active. This is unfortunately not displayed on the screen. For several minutes, while the screen is blank, the flash PROMs are loaded into the generator. This process must under no circumstances be disturbed! At the end of this sensible procedure "Flash loaded o.k." appears on the screen. Only now the scope program can be terminated.

· Reset the generator.

REPLACEMENT **OPTIMUS R/F**

Replacement of parts of "Function Unit kV" 4.

In case one of the following assemblies:

- PCB kV-control 3 / 4 (EZ130)
- Converter (EQ)
- H.V. tank (EG)
- PCB Central unit (EZ139)

of "Function Unit kV" has been exchanged, the alignment of the "Function Unit kV" must be repeated.


4.1. **General information**

The actual value of the set kV must be attained at least after 2ms. At kV rise phase there must be neither kV break-in nor a kV overshoot.

The Factor Duty Cycle is based on an adapted tube and determines at local mains voltage and mains resistance conditions:

- the kV rise phase
- the kV behavior during the exposure in falling load technique

The Factor Duty Cycle is stored in the memory of PCB CU EZ139. If the CU has to be replaced the CU complete backup can be reloaded (with the actual factor) to the NVRAM memory or the Factor Duty Cycle must be re-aligned. Refer to figure 1:




During alignment this Factor Duty Cycle must be entered via XRGSCOPE. The influence of this factor as a correction value for the Z-Data Table is monitored as the U_{COMP} signal, since the PI-Controller is deactivated by the grounded GAIN_IN signal.

OPTIMUS R/F REPLACEMENT

4.2. Connecting and setting the scope

For connections see figure 2:



Channel 1 = EZ130 X3 ---> AV HT ---> 20kV/V ---> 1V/div --> Zero-line at bottom of screen

Probe GND = one of the drilling holes at the front cabinet chassis

Channel 2 = EZ130 X26 ---> U_{COMP} ---> 1V/div ---> Zero-line 2 div from bottom of screen

Probe GND = one of the drilling holes at the front cabinet chassis

Trigger = external (preferred) ---> CTRL_X_C/ ---> backpanel EZX74 / negative slope or = internal channel 1 ---> AV HT ---> EZ130 X3 / positive slope at +3V

Probe GND = one of the drilling holes at the front cabinet chassis

Time base = 5 or 10ms/div ---> trigger delay -1div

Note

A digital scope should not have any other ground connection than the ground of the 3 probes at the drilling holes at the front generator chassis.

A mains-driven scope must be isolated from earth potential, otherwise it might display artefacts.

REPLACEMENT **OPTIMUS R/F**

4.3. Deactivating the kV controller

• Connect EZ130 X23 GAIN IN and X6 GNDA with a short link (use a short wire).

Caution!

This alignment requires exposures with high kV. Be sure the tube has been warmed up before.

4.4. Setting of exposure data

a) Set 141kV in case of

- 65/80kW
- the tube limit (of at least one tube) is 150kV perform this adjustment at the tube which has the highest kV limit programmed.

b) Set 125kV in case of

- 50kW

and

- 65/80kW if the programmed application limit of the tube limit is 125kV.

Any tube arcing during this adjustment requires the execution of the tube conditioning next as described in section 2 "INSTALLATION".

Disconnect the short link between X23 and X6.

Start over this adjustment from chapter 2.3 onwards if the tube conditioning was successful.


- Set kV and mA values according to the programmed tube limits:
 - a) 141kV: 200mA at kV_4 (65/80kW)
 - b) 125kV: at kV 3 100mA (50kW)

200mA at kV 4 (65/80kW)

• Set exposure time: 40ms OPTIMUS R/F REPLACEMENT


4.5. Adjustment of the "factor for duty cycle"

- Adjust the Factor Duty Cycle via service software XRGSCOPE by measuring U_{COMP} with the scope.
- Connect the service PC and start XRGSCOPE:
 XRGSCOPE
 SELECT UNIT ---> FU-kV ---> ADJUST ---> IGBT Pulse Width Correction ---> Factor Duty Cycle
- Set the starting value Factor Duty Cycle to 1.00:



- If the *U_{COMP}* value does not match the requirements type in another Factor Duty Cycle value, <TRANSMIT> the screen and push the active RGDV button to get the new value validated.
- Switch an exposure.
 The values are measured in the stationary condition. The transient behavior at the beginning of the exposure is not taken into account.

Result: In standby the U_{COMP} value is at about +11V, during exposure the mean value U_{COMP} must be as given in table 1 or 2, refer to figure 4:



OPTIMUS RF 4 a021

REPLACEMENT OPTIMUS R/F

a) 141kV setting (65/80kW only)

Read the mean value of U_{COMP} for 141kV (see scope figure 4 or 5), correct the Factor Duty Cycle till U_{COMP} meets the required reference of +1V.


kV setpo	oint	mA setpoint	PCB type	U _{COMP}	Tolerance	•	Factor Duty Cycle:	Date:
141	kV	200mA	PCB kV_control 4:	+1V	±0.5V	138kV		

Table 1: Factor Duty Cycle, settings 141kV (150kV limit)

Example how to correct the Factor Duty Cycle:

PCB kV control 4:

- If the mean value of U_{COMP} is: > +1.5V increase the Factor Duty Cycle in steps of 0.01 decrease the Factor Duty Cycle in steps of 0.01
- Check also the kV peak value AV HT (not the overshoot), it must be 138kV for 141kV setpoint. (see scope figure 5)
- Remove short link EZ130 X23 GAIN_IN.
- · Record the findings in table1.



OPTIMUS R/F REPLACEMENT

b) 125kV setting (50/65/80kW)

- Read the mean value of U_{COMP} for 125kV (in principle figure 4 or 5).
- Correct the Factor Duty Cycle till *U_{COMP}* meets the required reference of 0V.

kV setpoint	mA setpoint	PUBIVE HOLD HOLD HOLD I		Factor Duty Cycle:	Date:		
125kV	100mA	PCB kV_control 3:	+0V	+1V / -0,5V	125kV		
125kV	200mA	PCB kV_control 4:	+0V	±0.5V	125kV		

Table 2: Factor Duty Cycle, 125kV limit

Example how to correct the Factor Duty Cycle:

PCB kV_control 3:

• If the mean value of U_{COMP} is: > +1V increase the Factor Duty Cycle in steps of 0.01 decrease the Factor Duty Cycle in steps of 0.01

PCB kV_control 4:

- If the mean value of U_{COMP} is: > +0.5V increase the Factor Duty Cycle in steps of 0.01 decrease the Factor Duty Cycle in steps of 0.01
- Check also the kV peak value AV HT (not the overshoot), it must be 125kV for 125kV setpoint.
- Remove short link EZ130 X23 GAIN_IN.
- · Record the findings in table 2.

REPLACEMENT **OPTIMUS R/F**

5. **Tube replacement**

Any new tube require a new adjustment procedure consisting of:

- 1. Tube conditioning
- 2. Tube adaptation



Warning!

Radiation is released during the adjustment procedure!

The generator must be in the READY state, i.e. the green ring at the desk must be illuminated!

5.1. **Tube conditioning**

5.1.1. Preconditions / Programmings

• Switch OFF the generator.

Preparation of generators which are connected via a CAN interface:

- BuckyDiagnost TH and TH2
- DigitalDiagnost
- Thoravision
- EasyDiagnost with bucky unit
- Disconnect the following plugs:

System	Connector						
	EZX23 signal bus	EZX42 or EZX42-1 system CAN	EZX43 or EZX43-1 system CAN				
BukyDiagnost TH / TH2	Х		X				
DigitalDiagnost	Х	Х	Х				
Thoravision	Х	Х	Х				
EasyDiagnost with Bucky unit	Х	Х	Х				

· Switch ON the generator.

Note

The programming procedure must not be started before relays ENK1 has been energized at least 2 minutes after the generator has been switched ON.

OPTIMUS R/F REPLACEMENT

 Perform the following programmings temporarily for each tube connected to one of the assigned RGDV = Free cassette

Select menu XRGSCOPE:

PROGRAM/ REGISTRATION DEVICES/ RGDV#/ DATA SET A

Programming	Temporarily	Original Tube		
Enable handswitch	YES			
Syncmaster present	NO			
Exposure switch type	Double Step	Verify the customized entries in 2Z-2.x		
Exposure series / Tomo	YES			
Mounted radiographic	NONE			

- · RESET the generator.
- Select appropriate programmed RGDV = "Free cassette" for the tube to be conditioned.

5.1.2. Procedure

· Select large focus only.

Note

The generator must be in the READY state.

- Run conditioning procedure for a new or non-adapted tube, refer to following table "Exposure parameters for conditioning".
- It is recommended that the high tension be monitored during conditioning.

Connect the scope:

Channel1: kV AV HT at EZ130 X3 (1V/div), scale: 20kV/V Trigger external: CTRL X C/ at backpanel EZ X74, negative slope

Time base: 2ms/div

• In case of problems like tube arcing see the following flowchart EXPOSURE SEQUENCE as an example. The flowchart applies for applicable kV range only, e. g.:

109kV is the max. kV value for normal application, perform just up to next higher kV step = 117kV.

Note

Refer to flowchart EXPOSURE SEQUENCE.

If the tube arcs at a certain kV value, switch another 3 exposures with same parameters and 10s pause between subsequent exposures. In case of success (no arcing anymore) continue with next kV step of the following table.


If the last exposure still arcs go one kV step back and follow normal procedure. If this routine has been performed three times without improvement: ==> **Replace the tube!**

REPLACEMENT OPTIMUS R/F

Exposure p	Exposure parameters for conditioning							
kV	mAs	# exposures						
80	0.5	<1>						
80	5	<1>						
80	50	<1>						
10 secon	ds pause							
80	100	< 1 >						
1 minute	e pause							
90	0.5	<1>						
90	5	<1>						
90	50	< 1 >						
10 secon	ds pause							
90	100	< 1 >						
1 minute	e pause							
100	0.5	<1>						
100	5	< 1 >						
100	50	< 1 >						
10 secon	ds pause							
100	100	<1>						
1 minute	e pause							
110	0.5	<1>						
110	5	< 1 >						
110	50	< 1 >						
10 secon	ds pause							
110	100	< 1 >						
1 minute	e pause							
120	0.5	<1>						
120	5	<1>						
120	50	<1>						
10 secon	ds pause							
120	100	<1>						
1 minute	-							
130	0.5	<1>						
130	5	<1>						
130	50	<1>						
10 secon	ds pause							
130	100	<1>						
1 minute	e pause							

Exposure parameters for conditioning						
kV	mAs	# exposures				
140	0.5	< 1 >				
140	5	<1>				
140	50	< 1 >				
10 secon	ds pause					
140*	100	<1>				
1 minute	e pause					
145	0.5	<1>				
145	5	< 1 >				
145	50	< 1 >				
10 secon	10 seconds pause					
145	100	< 1 >				
1 minute	e pause					
148	0.5	<1>				
148	5	< 1 >				
148	50	< 1 >				
10 secon	ds pause					
148	100	< 1 >				
1 minute	e pause					
150	0.5	<1>				
150	5	< 1 >				
150	50	< 1 >				
10 secon	ds pause					
150	100	< 1 >				
1 minute	e pause					

REPLACEMENT OPTIMUS R/F



Note

If a tube arcs at any kV value which is not required for application the max. kV (e.g.117kV) program this new limit value by XRGSCOPE:

PROGRAM/ TUBES/ TUBE LIMITS/ MAX. TUBE VOLTAGE LIMIT [kV]/ [117]

As the limit value decreases for this reason, a following re-adaptation procedure sets the field ADAPTED TO [kV] to this value as well.

- Set RGDV programming to original status if no adaptation procedure has to be executed.
- · RESET the generator.

OPTIMUS R/F REPLACEMENT

5.2. Tube adaptation

5.2.1. General information

Tube adaptation is an automatic process which includes:

- 1. The measurement of the mA offset value that is caused by:
 - the kV measuring circuit
 - the emission current feedback circuit (VCO)
- 2. The measurement of the individual standby filament current (based on $100\mu A$).
- 3. The emission current characteristic as f (kV, filament current).
- 4. The dynamic behavior (positive and negative boost adaptation) where the inertia of the filament with respect to heating up and cooling down is registered.

For more information refer to section 3: FAULT FINDING.

Note

In case of problems check the symptom / solution list at the end of this adjustment chapter. Repeat the adaptation for this particular focus.

REPLACEMENT OPTIMUS R/F

5.2.2. Preconditions / Programmings

· Reset the generator.

Note

The adaption procedure must not be started before relay ENK1 has been energized at least 2 minutes after the generator has been switched ON.

- The tube must be conditioned as described in chapter 5.1.
- · Check the upper kV limit

Select menu XRGSCOPE:

PROGRAM/ TUBES/ TUBE LIMITS/ MAX. TUBE VOLTAGE LIMIT [kV]

The programmed value should match the nominal value of the tube connected or in case of older tubes the upper kV limit should be set to the max. application kV.

Once an adaptation is completed the new limit value indicates as ADAPTED TO [kV].

· Perform the following programmings temporarily for each tube connected to one of the assigned RGDV = Free cassette

Select menu XRGSCOPE:

PROGRAM/ REGISTRATION DEVICES/ RGDV#/ DATA SET A

Programming	Temporarily	Original Tube
Enable handswitch	YES	
Syncmaster present	NO	
Exposure switch type	Double Step	Verify the customized entries in 2Z-2.x
Exposure series / Tomo	YES	
Mounted radiographic	NONE	

OPTIMUS R/F REPLACEMENT

5.2.3. Procedure

· RESET the generator.

• It is recommended that the high tension be monitoredduring adaptation.

Connect the scope:

Channel1: kV AV HT at EZ130 X3 (1V/div), scale: 20kV/V Trigger external: CTRL X C/ at backpanel EZ X74, negative slope

Time base: 2ms/div

• Select the RGDV = Free cassette for the tube to be adapted.

Select menu XRGSCOPE:
 OPTIMUS XRG/ ADJUST/ TUBE ADAPTATION

• Select the tube and focus to be adapted, start with small focus!

Note

To avoid any malfunction make sure that READY is displayed on the desk before transmitting data by pushing <F2>:

READY state disappears, ADAP is displayed on the desk, WAITING is displayed on the PC screen. Wait until the generator turns back to READY state.

Start the adaptation process by pushing the handswitch in PREP and EXP position continuously.

The generator switches about 125 exposures for each focus. The radiation sign at the desk indicates exposures but there is no beep at the end of each exposure.

The actual kV parameters are displayed during adaptation.

The generator carries out the adaptation automatically. The procedure for one focus is completed when the desk indication changes from ADAP to TEST. The WAITING message disappears from the PC screen together with a PC beep, followed by the screen: BEFORE CONTINUING THE GENERATOR MUST BE RESET.

- RESET the generator.
- Run the adaptation for each focus (small and large) and tube.

Note

As there is no tube type with a physical third (middle) focus yet, the third focus cannot be adapted. VARIOFOCUS values are calculated by adapted small and large focus. APR programs using VARIOFOCUS can only be selected until small and large focus are both adapted.

Set RGDV(s) programming to original status according table "RGDV programming" 2Z-2.x in the end of this
chapter.

OPTIMUS R/F **REPLACEMENT**

Final tube adjustment work 5.3.

- 1. BuckyDiagnost TH with CAN interface, DigitalDiagnost, Thoravision:
 - Switch OFF the generator.
 - Reestablish signal bus connector EZX23.
 - Reestablish CAN connectors EZX42-1 and EZX43-1.
 - Switch ON the generator.
- 2. All other systems:
 - · Reset the generator

© 1997 Philips Medizin Systeme AIL RIGHTS RESERVÉD

OPTIMUS R/F REPLACEMENT

5.4. Symptoms and solutions if problems occur during tube adjustment

1. A warning cannot be displayed on the control desk, the WAITING screen on the PC is flickering instead during this event and logged in the error log index.

2. If the tube has already been at a high temperature level (but the tube load indication still indicates green or green-yellow for 100% power) it might happen that the tube load indication changes straight to red and the adaptation is on hold. WAITING is flickering on the PC.


Solution: Keep the handswitch pushed, once the temperature is down adaptation continues automatically.


Note


An increment of one of the temperature levels inhibits the 100% power condition. This event is always logged as warning 00BV in the error log index.

- 3. An error message just flashes for a very short moment and instantly covered by ADAP on the desk afterwards. WAITING is flickering on the PC.
- 4. All buttons at the control desk including the RESET button are inactive during adaptation. The only way to RESET an error is to release the PREP switch which causes an interrupt similar to the RESET command.
- 5. After letting of the PREP switch wait until the desk indicates READY. If READY does not appear at least after 20 seconds run a warmstart of the generator by pushing the RESET button on CU EZ139.
- 6. If adaptation seems to do nothing for more than 30 seconds let go of the PREP switch. Wait until the desk indicates READY. If READY does not appear at least after 20 seconds run a warmstart of the generator by pushing the RESET button on CU EZ139.
- 7. If a constant READY indication appears for more than 2 seconds while PREP and EXP is activated by the handswitch during adaptation let go of the handswitch.
 Wait until the desk indicates READY. If READY does not appear at least after 20 seconds run a warmstart of the generator by pushing the RESET button on CU EZ139.
- 8. If adaptation does not carry on with or without READY indication check whether one of the function units indicates a FATAL error by turning on the red LED. Let go of the handswitch and warmstart the generator by pushing the RESET button on CU EZ139.
- 9. If adaptation has been interrupted by a generator warmstart check the error log index before restarting adaptation:
 - kV errors 02WG and/or 02WH indicate tube arcing. Run conditioning of the tube as described in section 2 "INSTALLATION" and/or reduce the max. kV value to the required application value.












W4	Application of RoCo HS
10 20 open 30	normal
20 2 - 3 closed	MRC solo



D1 = 4512 113 22321 = 4512 113 22322 = 4512 113 22323





D10 = 4512 113 20611

OPTIMUS R/F

ADJUSTMENTS

TEXT

	Contents	6-0.1
1.	Area Dose Calculator (option)	6-1
1.1.	Checking the default adjustment	6-2
1.2.	Correction of the specific yield	6-3
1.3.	Correction of the filter values	6-5
2.	Alignment of "Function Unit kV"	6-8
2.1.	General information	6-8
2.2.	Connecting and setting the scope	6-10
2.3.	Deactivating the kV controller	6-11
2.4.	Setting of exposure data	6-11
2.5.	Adjustment of the "factor for duty cycle"	6-12
3.	Dose rate alignment for continuous fluoroscopy when using the dose-rate signal DR TV (formerly ADC) // dose-rate control via TV-chain	6-15
3.1.	Aim	6-15
3.2.	Principle	6-16
3.3.	Alignment of DR TV with the oscilloscope	6-18
4.	Maximum patient entrance dose rate: Entrance Dose rate Limiter (EDL)	6-23
4.1.	Aim	6-23
4.2.	Principle	6-23
4.3.	Alignment	6-23

1. Area Dose Calculator (option)

Special tools:

- calibrated dosemeter, e.g. DALI with measuring cell 77334 or PMX3
- 1mm lead plate

The following parameters are relevant for calculation:

- SID (Source Image Distance)
- diaphragm aperture
- added filters
- specific yield of tube
- mAs product
- number of exposures

SID, diaphragm aperture and type of filters are supplied by the diagnostic unit, where they are also adjusted. In the generator default values are given for the specific yield of a tube and filter correction.

These default values can be found as reference files on floppy disk in order to recreate the original settings if need be.

Reference files:

ref_yiel.tdl specific yield of tube

ref 2al.tdl filter 2mm Al

ref_01cu.tdl filter 1mm Al + 0.1mm Cu ref_02cu.tdl filter 1mm Al + 0.2mm Cu

The specific yield curve relates to tungsten anodes and 2.5mm primary filters.

Display on the desk is in [cGycm²].

ADJUSTMENTS OPTIMUS R/F

Checking the default adjustment

• Place the lead plate and the measuring cell of the measuring instrument on the table in the central radiation beam. The purpose of the lead plate is to reduce radiation scatter of the table top. Without the plate the test result would be approximately 10% higher using a table top made, for example, of resin bonded paper.

- · Perform the following settings:
 - 1m between the focus and the measuring cell (= SMD Source Measuring Distance)
 - free cassette technique
 - kV-mAs-s technique
 - 10mAs
 - 0.1s
 - collimation 10 x 10cm at the height of the measuring cell
 - no filter
- Determine area dose at the following kV settings and compare it with the respective value displayed on the desk:

	50kV	80kV (81kV)	120kV (117kV)
displayed product	cGycm ²	cGycm ²	cGycm ²
measured dose	сGу	cGy	cGy
measured area (X x Y)	cm ²	cm ²	cm ²
calculated product	cGycm ²	cGycm ²	cGycm ²
difference	%	%	%

Example:

- displayed area dose product: 8.8cGycm²
- measured dose: $890\mu Gy = 0.089cGy$
- calculated area dose product: measured dose × exposed area = 0.089cGy × 100cm² = 8.9cGycm²
- 8.9 8.8 - difference in %: ----- x 100 = 1.12%
- If there are any deviations of over 5% it is recommended that the yield curve be corrected in accordance with the procedure described in 1.2.

1.2. Correction of the specific yield

Prerequisite:

- Test setup and settings in accordance with section 1.1.:
 - 1m between the focus and the measuring cell (= SMD)
 - free cassette technique
 - kV-mAs-s technique
 - 10mAs
 - 0.1s
 - collimation 10 x 10cm at the height of the measuring cell
 - no filter

Principle:

• For each kV specified a dose measurement is taken under the same conditions. If the distance between the focus and the measuring cell deviates from 1m, all the dose values must be corrected with the square of distance (unit of measurement is [m]). Dividing the dose values by the mAs product set gives the respective current yield.

Procedure:

• Measure dose at each kV checkpoint and use it to calculate specific yield.

The values determined must be higher at higher kVs settings and produce a characteristic with a slight curve on the graph. If considerable fluctuations are detected, the measurements must be repeated at the points in question.


Range: 0.00 ... 400.00μGy/mAs

The values can only be stored in the generator if they are within the range specified and rise uniformly with kV.

ADJUSTMENTS OPTIMUS R/F

Specific yield

kV checkpoint	40	50	60	70	80	90	100	110	120	130	140	150
default yield [μGy/10mAs]	11.0	20.8	32.2	45.7	58.0	73.5	88.7	106.1	124.4	143.6	163.1	181.8
measured dose [μGy]												
distance ² factor	If the	If the distance focus – measuring cell (= SMD) differs from 1m correct the dose with the factor: distance ² factor = $(SMD [m] / 1m)^2 = \dots$ (e.g. = 1.44 for a SMD of 1.2m)										
corrected dose [μGy]												
specific yield					specifi	c yield	= correcte	ed dose/	I0mAs			
[μGy/10mAs]												



Correct the default values of the specific yield for all the kV checkpoints.
 Select menu:

ADJUST/ AREA EXPOSURE PRODUCT/ **SPECIFIC YIELD OF TUBE 1 ... 3** with the factor determined and save with <TRANSMIT>.

• Save the specific yield curve with the SAVE function of XRGSCOPE <F3> key on the backup disk. Recommended file name: act_yiel.tdl

1.3. Correction of the filter values

Prerequisite:

- Test setup and settings in accordance with section 1.1.
 - 1m between the focus and the measuring cell (= SMD)
 - free cassette technique
 - kV-mAs-s technique
 - 10mAs
 - 0.1s
 - collimation 10 x 10cm at the height of the measuring cell
 - no filter

Principle:

At otherwise identical settings the dose is determined for the kV values specified with and without filter. The ratio
of dose values with / without filter produces the respective current correction factor.

Procedure:

- Accept measured dose values (not the corrected ones!) for the respective kV checkpoints from yield measurement or measure them again if any changes have been made to the test-setup or settings.
- Move the filter to be checked into the radiation beam.
- Measure dose at each kV checkpoint and enter it in the respective table.

Note

The 40kV range is not used in practice so it does not have to be corrected.

If in the lower kV range the considerably reduced dose can no longer be measured or read perfectly, at that point a higher mAs product must be selected. Then the repeat measurement must be performed without filter.

• Using the ratio between dose with and without filter, determine the respective correction factor.

The values determined must be higher at higher kVs settings and produce a characteristic with a slight curve on the graph. If considerable fluctuations are detected, the measurements must be repeated at the points in question.

Range: 0.000 ... 1.000

The values can only be stored in the generator if they are within the range specified and rise uniformly with kV.

Perform the procedure for each selectable filter type.

ADJUSTMENTS OPTIMUS R/F

Filter correction - 2mm Al


kV-checkpoint [kV]	50	70	100	150
default factor	0.47	0.56	0.66	0.75
measured dose [μGy] with filter				
measured dose [μGy] without filter				
a con factor	new facto	r = dose with 1	filter / dose wi	thout filter
new factor				

Filter correction - 1mm Al + 0.1mm Cu

kV-checkpoint [kV]	50	70	100	150
default factor	0.25	0.37	0.50	0.65
measured dose [μGy] with filter				
measured dose [μGy] without filter				
and facility	new facto	r = dose with 1	filter / dose wi	thout filter
new factor				

Filter correction - 1mm AI + 0.2mm Cu

kV-checkpoint [kV]	50	70	100	150	
default factor	0.123	0.230	0.370	0.53	
measured dose [μGy] with filter					
measured dose [μGy] without filter					
	new factor = dose with filter / dose without filter				
new factor					



· Read out the default values of the filter tables for each kV checkpoint, correct with the factor determined and write back into the generator with <Transmit>.

Select menu:

ADJUST/ AREA EXPOSURE PRODUCT/ ADD FILTER CORRECTION TABLES/ ...

- ... 2mm AL
- ... 1mm AL + 0.1mm CU
- ... 1mm AL + 0.2mm CU
- Save the specific correction tables with the SAVE function of XRGSCOPE <F3> key on the backup disk. Recommended file names:

act 2al.tdl filter 2mm Al

act 01cu.tdl - filter 1mm Al + 0.1mm Cu act_02cu.tdl filter 1mm Al + 0.2mm Cu **ADJUSTMENTS** OPTIMUS R/F

Alignment of "Function Unit kV" 2.

General information 2.1.

The actual value of the set kV must be attained at least after 2ms. At kV rise phase there must be neither kV break-in nor a kV overshoot.

The Factor Duty Cycle is based on an adapted tube and determines at local mains voltage and mains resistance conditions:


- the kV rise phase and
- the kV behavior during the exposure in falling load technique as it takes into account the tolerances of the following FRUs (Field Replaceable Units):

```
1. PCB EZ 130
  kV_control_3 = 50kW
                              1 converter
                                              4512 108 0908x
  kV control 4 = 65/80kW
                              2 converters
                                              4512 108 0910x
                                                                **
2. A complete power converter unit Q
  kV power PCB(s) Q100
                              (part of the power converter unit)
  IGBT transistors
                              (part of the power converter unit)
3. Resonance capacitors
                              (part of the power converter unit)
4. High tension transformer
                                                                **
```

An exchange of one of the ** marked parts requires a realignment of the Factor Duty Cycle.

6-8 OPTIMUS R/F (b/02.0)© 1997 Philips Medical Systems OPTIMUS_RF_6_b020 ALL RIGHTS RESERVED

The Factor Duty Cycle is stored in the memory of PCB CU EZ139. If the CU has to be replaced the CU complete backup can be reloaded (with the actual factor) to the NVRAM memory or the Factor Duty Cycle must be re-aligned. Refer to figure 1:




During alignment this Factor Duty Cycle must be entered via XRGSCOPE. The influence of this factor as a correction value for the Z-Data Table is monitored as the U_{COMP} signal, since the PI-Controller is deactivated by the grounded $GAIN_IN$ signal.

ADJUSTMENTS OPTIMUS R/F

2.2. Connecting and setting the scope

For connections see figure 2:



Channel 1 = EZ130 X3 ---> AV HT ---> 20kV/V ---> 1V/div --> Zero-line at bottom of screen

Probe GND = one of the drilling holes at the front cabinet chassis

Channel 2 = EZ130 X26 ---> U_{COMP} ---> 1V/div ---> Zero-line 2 div from bottom of screen

Probe GND = one of the drilling holes at the front cabinet chassis

Trigger = external (preferred) ---> CTRL_X_C/ ---> backpanel EZX74 / negative slope or = internal channel 1 ---> AV HT ---> EZ130 X3 / positive slope at +3V

Probe GND = one of the drilling holes at the front cabinet chassis

Time base = 5 or 10ms/div ---> trigger delay -1div

Note

A digital scope should not have any other ground connection than the ground of the 3 probes at the drilling holes at the front generator chassis.

A mains-driven scope must be isolated from earth potential, otherwise it might display artefacts.

2.3. Deactivating the kV controller

• Connect EZ130 X23 GAIN IN and X6 GNDA with a short link (use a short wire).

Caution!

This alignment requires exposures with high kV. Be sure the tube has been warmed up before.

2.4. Setting of exposure data

a) Set 141kV in case of

- 65/80kW
- the tube limit (of at least one tube) is 150kV perform this adjustment at the tube which has the highest kV limit programmed.

b) Set 125kV in case of

- 50kW
 - and
- 65/80kW if the programmed application limit of the tube limit is 125kV.

Note

Any tube arcing during this adjustment requires the execution of the tube conditioning next as described in section 2 "INSTALLATION".

Disconnect the short link between X23 and X6.

Start over this adjustment from chapter 2.3 onwards if the tube conditioning was successful.

- Set kV and mA values according to the programmed tube limits:
 - a) 141kV: 200mA at kV_4 (65/80kW)
 - **b) 125kV:** 100mA at kV 3 (50kW)


200mA at kV 4 (65/80kW)

• Set exposure time: 40ms

ADJUSTMENTS OPTIMUS R/F


2.5. Adjustment of the "factor for duty cycle"

- Adjust the Factor Duty Cycle via service software XRGSCOPE by measuring *U_{COMP}* with the scope.
- Connect the service PC and start XRGSCOPE:
 XRGSCOPE
 SELECT UNIT ---> FU-kV ---> ADJUST ---> IGBT Pulse Width Correction ---> Factor Duty Cycle
- Set the starting value Factor Duty Cycle to 1.00:



- If the *U_{COMP}* value does not match the requirements type in another Factor Duty Cycle value, <TRANSMIT> the screen and push the active RGDV button to get the new value validated.
- Switch an exposure.
 The values are measured in the stationary condition. The transient behavior at the beginning of the exposure is not taken into account.

Result: In standby the U_{COMP} value is at about +11V, during exposure the mean value U_{COMP} must be as given in table 1 or 2, refer to figure 4:



a) 141kV setting (65/80kW only)

Read the mean value of U_{COMP} for 141kV (see scope figure 4 or 5), correct the Factor Duty Cycle till U_{COMP} meets the required reference of +1V.


kV setpoint	mA setpoint	PCB type	U _{COMP}	Tolerance	-	Factor Duty Cycle:	Date:
141kV	200mA	PCB kV_control 4:	+1V	±0.5V	138kV		

Table 1: Factor Duty Cycle, settings 141kV (150kV limit)

Example how to correct the Factor Duty Cycle:

PCB kV_control 4:

- If the mean value of U_{COMP} is: > +1.5V increase the Factor Duty Cycle in steps of 0.01 decrease the Factor Duty Cycle in steps of 0.01
- Check also the kV peak value AV HT (not the overshoot), it must be 138kV for 141kV setpoint. (see scope figure 5)
- Remove short link EZ130 X23 GAIN_IN.
- · Record the findings in table1.



ADJUSTMENTS OPTIMUS R/F

b) 125kV setting (50/65/80kW)

- Read the mean value of U_{COMP} for 125kV (in principle figure 4 or 5).
- Correct the Factor Duty Cycle till U_{COMP} meets the required reference of 0V.

kV setpoint	mA setpoint	PCB type	U _{COMP}	Tolerance	•	Factor Duty Cycle:	Date:
125kV	100mA	PCB kV_control 3:	+0V	+1V / -0,5V	125kV		
125kV	200mA	PCB kV_control 4:	+0V	±0.5V	125kV		

Table 2: Factor Duty Cycle, 125kV limit

Example how to correct the Factor Duty Cycle:

PCB kV_control 3:

• If the mean value of U_{COMP} is: > +1V increase the Factor Duty Cycle in steps of 0.01 < -0.5Vdecrease the Factor Duty Cycle in steps of 0.01

PCB kV_control 4:

• If the mean value of U_{COMP} is: > +0.5V increase the Factor Duty Cycle in steps of 0.01 decrease the Factor Duty Cycle in steps of 0.01 < -0.5V

- Check also the kV peak value AV HT (not the overshoot), it must be 125kV for 125kV setpoint.
- Remove short link EZ130 X23 GAIN IN.
- · Record the findings in table 2.


Dose rate alignment for continuous fluoroscopy when using the 3. dose-rate signal DR TV (formerly ADC) // dose-rate control via TV-chain

3.1. Aim

The amplification factor of the control circuit must be adjusted such that

- the TV image becomes visible within about 400ms after fluoroscopy is switched ON with RQXGFL
- there is no flashing of the TV image (overshoot < 20%).

Also see sketch below:




ADJUSTMENTS OPTIMUS R/F

3.2. **Principle**

Adjusting value Transition function of DR TV

Variables TV_pos_limit / TV_neg_limit / Gain factor



• Set the gain factor to default value = 1.00 with FW XRGSCOPE:

OPTIMUS

--> Adjustments

--> Dose rate control

--> Amplification gain

--> CONT

Gain factor 0.5...2 (default = 1): [1.000]

• With fluoroscopy switched OFF, measure the black level of the dose-rate signal of the TV chain. Suitable test points:

```
- WN 25 DRDFP (with a TV chain XTV 11)
```

The kink of the characteristic curve TV pos limit to be adjusted must be:

TV pos limit = DR TV PO (black level measured) - 500mV

DR TV PO = 7.3V → adjusting value TV pos limit = 6.8V **Example:**

Set TV_pos_limit via XRGSCOPE:

OPTIMUS

--> Program

--> Dose rate control

--> CONT:

Program settings for continuous fluoroscopy:

```
+-----+
| scantime TV [ms]: [ 20.000] |
                                   \rightarrow 50Hz = 20ms
                     Yes
 scantime TV valid:
                                  → yes with control via TV chain.
                                     With photo diode = no.
P_max_EDL [W]: [ 250]
                                  → 250 default: Effective only when SID
                                     WB X12:1 DRLM is switched ON.
TV pos limit [V]:
                                  \rightarrow DR TV PO - 500mV.
                   [6.8]
                                    (TV black voltage - about 500mV)
| TV_neg_limit [V]: [-6.8] |
                                   \rightarrow = - TV pos limit [V]
  < OK >
|-----|
```

In case the switch ON procedure is fast enough with this setting

T1 \leq 400ms in quick mode or T1 \leq 500ms in normal mode and the TV chain does not flash when fluoroscopy is switched ON,

overshoot ≤ 20%

the alignment with the oscilloscope described below is **not** required.

ADJUSTMENTS OPTIMUS R/F

3.3. Alignment of DR TV with the oscilloscope

Settings of the oscilloscope:

Y1 to EZ 150 X6 EN X C/ [10V/div]

Y2 e.g. with XTV11

> to WN 25 **DRDFP** [5V/div]

Ground to EZ 150 X5 **GNDA** Time base 200ms/div Y1 / - DC Trigger

The alignment should not be carried out in the quick mode since the differences of the transition functions are not big enough in the quick mode. However, in the normal mode the overshoot is slightly higher than 20%.

Since fluoroscopy is automatically in the quick mode for up to 30s after the end of fluoroscopy, after each oscillogram fluoroscopy must be switched OFF for at least 30s:

```
Pause > 30s → Fluoro ON → Oscillogram → Fluoro OFF
→ Pause > 30s → Fluoro ON → Oscillogram → Fluoro OFF
  → ...
```

Fluoroscopy adjustments for all oscillograms:

- Dark voltage of the TV chain measured: DR TV PO dark = 7.3V

XRGSCOPE:

```
OPTIMUS
```

--> Program

--> Dose rate control

--> **CONT**:


```
scantime TV [ms]
                       : 20.000 (\rightarrow 50Hz)
scantime TV valid
                       : Yes
P max EDL [W]
                        : 250
```

Operating point of fluoroscopy in the automatic mode (example XTV 11):

- 9" I.I. format
- Filter in the path of radiation: 1.5mmCu + 20mmAl ==> 74kV/ 2.2mA

First reference oscillogram:

```
- TV_pos_limit [V]
                   : [ 7.3 ]
                                     ==> T1 ≈ 700ms --> too long --> reduce TV_pos_limit [V]
- TV_neg_limit [V]
                    : [-7.3]
- Gain factor
                     : [ 1.00]
```




ADJUSTMENTS OPTIMUS R/F


Second reference oscillogram:

```
- TV_pos_limit [V] : [ 6.8 ] 


- TV_neg_limit [V] : [-6.8 ] 

- Gain factor : [ 1.00 ] ==> T1 \approx 420ms --> o.k. and gain factor can remain as they are
```





Third reference oscillogram:



ADJUSTMENTS OPTIMUS R/F


Forth reference oscillogram:

```
- TV_pos_limit [V] : [ 6.8] 

- TV_neg_limit [V] : [ -6.8] 

- Gain factor : [0.50] 

=> T1 \approx 600 ms --> too long --> Raise gain factor to 1.00!
```



OPTIMUS R/F ADJUSTMENTS

4. Maximum patient entrance dose rate: Entrance Dose rate Limiter (EDL)

4.1. Aim

In some countries the maximum patient entrance dose rate must be limited to comply with the respective national regulations, for instance $10R/min = 1.4567\mu Gy/s$ in the USA.

4.2. Principle

The maximum patient entrance dose is measured with maximum kV. By reducing the maximum power P_max_EDL provided by the generator, the dose rate is reduced until it complies with the national regulations or the individual wish of the user.

Note

The dose rate reduction programmed via P_{max} _EDL is effective only when the SID contact (**S**ource **I**mage **D**istance) is switched ON: WB X12:1 DRLM \rightarrow closed!

4.3. Alignment

- Cover the I.I. with lead such that with fluoroscopy maximum kV (110kV) are adjusted (protection of the I.I. against excessive radiation).
- Adjust the exam./aux. unit in such a way that the SID contact is closed.
 If necessary, place a shorting plug on decade: WB X12:1 <---> X12:10 (ground).
- Adjust the exam./aux. unit in such a way that the smallest distance between tube and patient is effective.
- Position the dose rate measuring chamber in such a way that it is at the same height where the radiation would enter the patient.
- · Switch ON fluoroscopy.
- Reduce/raise **P_max_EDL [W]** with the service software XRGSCOPE until the dose rate is identical with the maximum permissible value:

```
XRGSCOPE OPTIMUS
```

```
--> Program
```

--> Dose rate control

--> CONT:

Range: 0 ... 9000W (250 default)


 $\,\rightarrow\,$ Effective only when SID contact WB X12:1 DRLM is switched ON.

ADJUSTMENTS OPTIMUS R/F

• Note the measuring values:

P_max_EDL: [W]							
Dose rate: [μGy/s]							

Example of a measuring setup:



OPTIMUS R/F

ACCEPTANCE

	Contents	7-0.1
1.	Preface	7-1
2.	Test equipment	7-1
3.	Setup	7-1
4.	Test	7-2
5.	Exposure counter	7-3

OPTIMUS R/F ACCEPTANCE

1. Preface

The national rules for accepting an X-ray system are very different. Therefore in the following an example is given for checking the generator in the USA.

OPTIMUS generators are factory-calibrated and checked for compliance with the parameter readout tolerances as stated in the relevant Operator's Manuals.

Provided that these generators are installed and set to work in accordance with the Installation Manuals, only the following limited field compliance testing is required.

2. Test equipment

- Keithley voltage divider model No. 35080 with filter packs 32867C, 5C, 9C or equivalent.
- Oscilloscope (storage)
- Digital mA, mAs meter.

Note

Do not start test until generator has been switched ON for at least one hour.

Direct (invasive) kVp measurements on OPTIMUS generators with HV divider tanks normally available to the field service organization are not permitted.

Measurements of kV using instruments other than the Keithley instrument may lead to larger measuring tolerances. The causes are to be found in the specific frequency response and transient response of each test instrument.

3. Setup

- Switch OFF generator and also switch OFF mains.
 Disconnect breaker to system.
- Connect digital mA meter as per instructions in the relevant Service Manual.
- Set up the Keithley voltage divider complete with the appropriate filter as per Keithley Instructions Manual No. 3294 OIM.
- · Connect the oscilloscope to the Keithley divider.

Note

Make sure that the oscilloscope has been calibrated with the aid of the Keithley divider as described in the Keithley Instructions Manual before starting any testing (par. 3.6. Internal Calibration).

• Calculate rejection limits based on the exposure parameter "Specification Limits" shown in the table below.

The "Specification Limits" are based on the actual tolerances as listed in the generator Operator's Manuals. These "Specification Limits" must be restricted to include the actual measuring instrument error. See also section 6, chapter 3.2 of the "COMPREHENSIVE COMPLIANCE TESTING MANUAL", No. 4535 800 2034x. regarding how to calculate rejection limits.

4. Test

- Switch the system ON.
- Measure the mains voltage on ENF1.

Reference voltage: Mains voltage programmed ± 10%

Actual values: L1 - L2: V

L1 – L3: V

L2 - L3: V

- · Select the largest focus.
- Release exposures according to the table below and compare the values measured with the reference values.

Technique	Parameter	Reference range	Measured value	Corrected value
	81kV ±5% ±1kV	76 86kV	kV	_
3-knob technique	250mA ±5% ±0.5mA	237 263mA	mA	mA
	100ms ±5% ±0.5ms	94.5 105.5ms	ms	
2-knob technique	125kV ±5% ±1kV	118 132kV	kV	_
	80mAs ±3% ±0.5mAs	77.1 82.9mAs	mAs	mAs

Owing to an offset current in the measuring circuit of the HV generator the measured values for mA / mAs must be adjusted using the following formulas:

$$I_{corrected} \ [mA] \ = \ I_{measured} \ [mA] \ - \ \frac{U \ [kV]}{R_{calc} \ [M\Omega]} \qquad \qquad Offset \ \approx \ 0.2 \ ... \ 0.75 mA$$

$$Q_{corrected} \text{ [mAs] = } Q_{measured} \text{ [mAs] - } \frac{\text{U [kV]} \times \text{t [s]}}{\text{R}_{calc} \text{ [M}\Omega]} - \frac{\text{4.55 [nF]} \times \text{U [kV]}}{\text{1000}}$$

$$Cable \text{ charge for 20m HV cable}$$

R_{calc} = calculated measuring circuit resistance.

Typical value: $\approx 200 M\Omega$

Read out R_{calc} via service menu: FU_MA/ FAULT FIND/ READ I_e CORRECTIONS

Focus assignment: Focus 1 = tube 1 large focus 2 = tube 1 small focus

3 = tube 2 large focus 4 = tube 2 small focus 5 = tube 3 large focus

6 = tube 3 small focus

t = exposure time according to desk display

OPTIMUS R/F **ACCEPTANCE**

5. **Exposure counter**

Before handing over the generator to the customer, read the exposure counter.

ACCEPT/ INSPECT/ TUBE STATISTIC/ TUBE 1 ... 3 STATISTIC/ SHOW TUBE STATISTIC Record the figure in the table below.

Tube load statistic variable	Unit	Tube1	Tube 2	Tube 3
Reset date	dd.mm.yy			
Last update	dd.mm.yy			
Preparation time large focus	S			
Preparation time small focus	s			
Preparation time vario focus	S			
Preparation counter large focus	1			
Preparation counter small focus	1			
Preparation counter vario focus	1			
Fluoro time	min			
Fluoro counter	1			
Exposure counter large focus	1			
Exposure counter small focus	1			
Exposure counter vario focus	1			
Overload exposures counter large focus	1			
Overload exposures counter small focus	1			
Overload exposures counter vario focus	1			

The tables should be reset whenever the tubes are being replaced.

ACCEPT/ INSPECT/ TUBE STATISTIC/ TUBE 1 ... 3 STATISTIC/ RESET TUBE STATISTIC Record the figure in the table above.

ACCEPTANCE OPTIMUS R/F

Explanation:

Reset date / Last update:

Reset date and date of last update of the tube statistic.

Preparation time:

The sum of all preparation times per focus.

Preparation counter:

Counts the occurrences of transition STANDBY or FLUORO to PREPARATION per focus.

Fluoro time:

The sum of all fluoro times.

Fluoro counter:

Counts the fluoro commands.

Exposure counter:

Counts the exposures per focus (including the overload exposures).

Overload exposures counter:

Counts the exposures at overload conditions of the tube.

Field Change Order Checklist


oment:		Туן	pe No.:
			al No.:
		Implemented	
FCO No.	Date	Signature	Remarks
			<u> </u>


FCO No.		Implemented	Remarks	
	Date	Signature	nemarks	
	<u> </u>			


Schematic drawings

Basis

Block diagram basis	
Power supply	Z1-2.2
kV control	Z1-3.3
mA control	
Central unit	Z1-5.1
Basic interface	Z1-6
Options	
Operating panel C Button and display arrangement	
EY Dual speed rotor control 9890 000 0268x	Z1-13.2
Tube extension overview	
Adapter 4 auxil. units WA/1WA/2WA	








Block diagram Expansions


EZX45


1WA X52 2WA X52 1WB X52

OPTIMUS R/F
© Philips Medical Systems




Power supply







OPTIMUS R/F © Philips Medical Systems Z1-2.3



OPTIMUS R/F (a/01.0)
© Philips Medical Systems





OPTIMUS R/F (a/99.0)© Philips Medical Systems

H.V. generator R/F

EZ139

4512 108 0920x


EZX45

EZX46 Operating panel


4512 983 06551


© Philips Medical Systems


(97.1)



Basic interface









42/A3 98-08-17 C

(b/98.0)






WA / 1WA / 2WA Adapter 4 auxil. units

OPTIMUS R/F
© Philips Medical Systems

(97.0)

Z1-15.1



WB Adapter 4 auxil. units

OPTIMUS R/F

(97.1)

© Philips Medical Systems

OPTIMUS R/F Section Z2

Wiring diagrams

Cabinet E	Z2-1.0
Cabinet wiring E, 50 kW	
Cabinet wiring E, 65/80 kW	Z2-1.2
Earthing diagram E	
Power supply N, 50 kW	Z2-2.1.1/.2
Power supply N, 65/80 kW	Z2-2.2.1/.2
Backpanel EZ / Basis rack-2 Z 4512 108 0936	
Backpanel EZ / Basis rack-2 Z survey of components	Z2-5.4
Dual speed rotor control	
9890 000 0268x	Z2-13
Cabinet wiring: Tube extension WG and air coolers H.V. generator	Z2-14.1
Tube extension WG	
Tube extension 1WG/2WG	
Air coolers H.V. generator	Z2-14.4
Cabinet wiring: Decade adapter 4 auxil. units WA/WB	
Adapter Photomultiplier (SEV) WP	Z2-15.1
Cabinet wiring: 26V DC / 230V AC Adapter WR	Z2-16
Cabinet wiring: Control desk C	

Option rack W Tube extension WG/ 230V/24V Adapter WR Adapter decade Basis rack Z cable WA Power supply N Rotor control Y 2 nd converter 2 Q Converter Q H. V. generator G


Cabinet E

Rear side


OPTIMUS R/F

Front side

front




Cabinet wiring 50kW R/F




Cabinet wiring 65/80 kW R/F


Earthing diagram


(c/98.0)







N Power supply 50kW R/F











N Power supply 65/80 kW R/F






















For survey of the components see Z2-5.4

Backpanel Basis rack-2 Z 4512 108 0936.

SW UN EX

RM DR CT


V226 V223 V222 V212 V221

Room decade

PO 26V


CU CT12

PO 15V




For survey of the components see Z2-5.4

Backpanel EZ Basis rack-2 Z 4512 108 0936.



Backpanel EZ Basis rack-2 Z 4512 108 0936.




Backpanel EZ
Basis rack 2 Z
survey of components

00-09-04

(d/00.0)


© Philips Medical Systems





OPTIMUS R/F © Philips Medical Systems


Schr.

Z2-14.1












2WG

1WG/2WG Tube extension







M Air coolers H.V. generator

WP

X41 □ = Z X19

Cabinet rear side

WA X52

WA X51

or WB X51

₩A X42

X42

2WA

41 Z X19 (2x WA) or WB X42

■ 2WA X41 (2x WA)

or WB X52

WA

Z X45 (2x WA)

■ 2WA X51 (2x WA)

CAN/XS

Cabinet rear side ₩B X42 X42 ■ 2WA X41 **■** WB X52 = 2WA X51 termination resistor W W WP WA 2WA X41 □ = Z X19 X42 ☐ WA X41 WB X51 **■** Z X45 **⇒** WA X51 X19 WB X41 ➡ **WB** 3 Adapters: WA/2WA/WB both WA adapters installed in the option rack WB adapter always installed beside the backpanel

> WA/WB/WP Cabinet wiring: Decade adapter 4 auxil. units Adapter Photomultiplier (SEV)

(97.0)


X1 = Z X14


= Z X21/X22/X31/X32/X41

WP

Cabinet rear side

WA X51 or WB X51






Note:

Adapter 26V-/230V is possible only in connection with decade adapter WA/WB (adapter "Old world").

WR Cabinet wiring 26V DC/230V AC Adapter WR



Cabinet wiring Control desk